# **WALLACH®**

**WA1000B™ General Purpose Cryosurgical System (N2O)  
(LL100™) 900506-1-N2O, 900506-2-N2O**

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**CE0086**

**READ ALL SAFETY INFORMATION AND INSTRUCTIONS BEFORE USING THIS PRODUCT!**

**DESCRIPTION**

The WA1000B™ cryosurgical console system consists of a freezer handle and probe with attached tubing, console and stand plus a cylinder yoke connector. The 900506-1-N<sub>2</sub>O is the console system without a Digital Temperature Indicator and Timer. The 900506-2-N<sub>2</sub>O is the console system with a Digital Temperature Indicator and Timer. A variety of cryosurgical tips are available for use with the unit. The WA1000B N<sub>2</sub>O unit requires Medical Grade N<sub>2</sub>O (Nitrous Oxide) gas.

Not made with natural rubber latex.

**CAUTIONS**

- U.S. Federal law restricts this device to sale by or on the order of a physician.
- Make sure console switch is turned off and pressure is depleted before removing the probe.
- Operating cylinder must be at room temperature in order to stabilize the gas pressure (68 °F to 74 °F, 18 °C to 22 °C). Cold will decrease the system's efficiency and significantly increases defrost time. Heat causes excess gas consumption.
- Allow time for temperature recovery if extremes have been experienced. Do not store gas cylinders near heaters, radiators, or other heat sources, as heat will increase the gas pressure within the cylinder.
- Do not mix or use gases other than N<sub>2</sub>O.
- Do not use wrenches, hammers, or other tools to connect Wallach equipment to gas cylinders.
- Use only non-siphon type gas cylinders.
- To ensure optimal performance, only medical grade gas free of contaminants should be used.

**INDICATIONS FOR USE**

The WA1000B console system is used for a variety of gynecology and dermatology cryosurgical procedures.

**WARNINGS**

- **Do not use overfilled cylinder; gas pressure may be excessive and present a hazard. If gas pressure gauge reads below 40 kg/cm<sup>2</sup>, replace the cylinder. If cylinder pressure is above the 70 kg/cm<sup>2</sup>, bleed excess gas from cylinder. Do not operate freezer if gauge is in the red.**
- **Move cylinder to a well-ventilated area when venting N<sub>2</sub>O. Gas exhaust port will accept an easily attached scavenger hose. Nitrous oxide in excessive concentrations may be hazardous to health.**
- **Visually inspect cryo instrument and tips prior to each use. Discontinue use of any product suspected of being damaged. Such product should be returned to Wallach Surgical Devices for inspection and service.**
- **At no time should chlorine mixtures (bleach) or iodine-based products be used for soaking or cleaning Wallach cryotips. Cryotips will degrade if exposed to such solution.**
- **In case of malfunction, close gas cylinder master valve.**
- **In the event of a defrost failure, turn tank off; leave tip in place. Remove tip when defrost is evident or use warm water to expedite.**
- **Do not depress triggers while changing tips. If triggers are depressed, injury may occur due to tip separation from cryosurgical freezer.**
- **Do not soak cryotips as lines can become clogged.**
- **Moisture in lines and orifices of the cryotips, and cryo console can cause sputtering, poor freeze and other malfunctions. Whenever using fluids or steam to clean or sterilize, all gas lines and cryotip orifices must be plugged.**

**WA1000B™ CRYOSURGICAL SYSTEM COMPONENTS**

1. GAS CYLINDER PRESSURE Gauge	Once cylinder is turned “ON” this gauge should read in the 40 to 70 kg/cm <sup>2</sup> range.
2. CONSOLES a) CONSOLE SWITCH  b) CONSOLE RECEPTACLE	a) When turned to the “ON” position permits gas flow into the probe.  b) Freezer connects to the console by screwing the knurled freezer connector to the receptacle.
3. GAS EXHAUST PORT	A gas scavenger hose may be attached to this port.
4. THE TIMER/TEMPERATURE INDICATOR CONSOLE (WA2000)	Monitors: <ul style="list-style-type: none"> <li>• Time (in seconds) during which the probe tip experiences freezing temperatures.</li> <li>• Temperature of the probe tip, when used with probes which have an integral thermocouple.</li> </ul> <p>You will be able to obtain a digital temperature reading from only one thermocouple at a time.</p> <p>NOTE: 4 “D” cell batteries are needed for this console. (Console will display “Low Batteries” when batteries are in need of replacement.)</p>
5. LL100™ Freezer	Multi-tip or fixed-tip freezer plugs into console.
6. POWER Switch	Used to turn console power “ON” or “OFF.”
7. SELECTOR Switch	Used for determining which of two thermocouple readings will be shown on the visual temperature readout panel.
8. TEMPERATURE INDICATOR Display	Readout panel for tip or tissue temperature.
9. TIMER Display	The timer segment monitors the time (in seconds) during which the probe tip experiences temperatures of –20 °C or colder. Once the probe tip reaches –20 °C, the timer will be activated and it will emit audible beeps at 15 second intervals, thus aiding the practitioner in determining the total freezing time. As the probe tip begins to freeze, the timer automatically activates the digital counter, and the timer continues to run until the probe begins a defrost mode, warming up to 0 °C, at which point the timing stops and the digital readout holds its current value. The practitioner has the option of resetting the timer to “zero” for a new count, or continuing the current count (cumulative) with an additional freeze. To reset to “zero” at any time, depress the “Reset” switch once.
10. Mobile Cart (WA1500)	Compact cabinet secures and conceals gas cylinder and allows transport of unit.

**DIRECTIONS FOR ASSEMBLY AND USE of THE WA1000B™ MOBILE FREEZER SYSTEM**

1. Attach the two console units to the cart by aligning the holes in the bottom of each unit with the locating pins on the top of the cart. Press each console down to firmly seat it in place on the cart. (The Timer/Temperature Indicator Console mounts on the upper rear part of the cart and the General Purpose Probe Console mounts on the upper front part of the cart.) To remove the console from the cart, simply pull straight up and the consoles will detach themselves from the cart.
2. Insert the stainless steel holder into the hole indicated on the General Purpose Probe Console.
3. To install a gas cylinder in the cart, separate the Velcro straps at the rear of the cart, place the cylinder in position (with the hook-up valve facing out) firmly against the cart frame and tightly secure the straps. Attach the gas connector hose of the Probe Console to the gas cylinder.  
**NOTE:** Use only “medical grade” Nitrous Oxide in Non-Siphon Gas cylinders.
4. Connect the LL100™ Freezer to the CONSOLE RECEPTACLE. Push the connector in and secure it in place by tightening the knurled ring of the probe connector. (Make it snug; hand tighten only.)

5. If the LL100 Multi-Tip Freezer is being used, select the probe tip to be used and screw it onto the Freezer shaft. The tip need only be finger tight. Do not use tools.
6. If a freezer with an integral thermocouple is selected, attach the thermocouple plug to the thermocouple receptacle (TC1 or TC2) on the Timer/Temperature Indicator Console. Be sure to match the color code of the thermocouple prongs to that of the receptacle, in order to assure a proper fit. Move the selector switch to the correct position for the thermocouple receptacle selected.
7. To reduce concentration of N<sub>2</sub>O in the immediate cryosurgery area, you may vent the exhaust gas by attaching one end of the exhaust hose to the GAS EXHAUST PORT. The other end of the hose goes to a discharge vent.  
NOTE: Hose must have minimum of 0.635 cm internal diameter with a maximum length of 609.6 cm. Wallach REF 901015 scavenger hose is recommended. Care must be taken to avoid kinking the hose.
8. Slowly open the gas cylinder master valve by turning the hand wheel counter-clockwise. Listen for any gas leaks. (There should be none.) Proceed to open the cylinder valve completely.
9. Once the gas pressure is stable, turn CONSOLE SWITCH to the “ON” position. The Freezer is now ready to operate.  
**PRECAUTION:** Make sure locking/freeze trigger is not engaged before turning console switch to the “ON” position.  
**NOTE:** It is recommended that an agent be used to enhance lubrication and moisture between the freezer and tissue. The surface of the probe should be positioned so as to face the surface of the tissue.
10. Depress the freeze (FRZ) trigger on the left of the LL100 Freezer. Within seconds you will observe frost on the cryotip. The freeze (FRZ) trigger will lock in position until the defrost (DEF) trigger on the right is depressed. When defrosting, depress and hold the defrost (DEF) trigger down for 2-3 seconds. You will see the frost disappear almost immediately. The sound of escaping gas is the normal gas exhaust coming out of the GAS EXHAUST PORT. (As mentioned, this gas may be vented via an optional exhaust hose.)  
**NOTE:** After a prolonged freeze, the white silicone hose exiting from the freezer handle may become cold. A light frosting may occur, but this is normal. Though cold to the touch, the hose will not freeze or damage tissue.
11. Within 3 seconds after releasing the defrost (DEF) trigger, the instrument is ready to refreeze.
12. When the procedure has been completed, shut the gas cylinder master valve “OFF” by turning the hand wheel clockwise.
13. Turn the CONSOLE SWITCH “OFF.”
14. You are now ready to remove the LL100™ Freezer. (If a freezer with an integral thermocouple has been used, remove the probe thermocouple plug from the thermocouple receptacle.) Remove the LL100 Freezer by loosening the knurled ring on the probe connector and pulling the probe connector gently, straight out of the CONSOLE RECEPTACLE. **DO NOT PULL ON THE FLEXIBLE HOSE.** Always use the body of the LL100 Freezer connector when connecting or disconnecting the probe.
15. Remove the gas connector hose of the Probe Console from the gas cylinder.

### TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSE	CORRECTIVE ACTION
1. Gas Cylinder Pressure gauge reads in the RED	The cylinder is over pressurized.	Take immediate action to reduce the gas pressure by closing the cylinder valve, loosening the Cryo connector to the cylinder valve, and carefully turning the cylinder valve “ON” again to bleed off the excess pressure. This procedure is best performed in a wide-open or outdoor area (in order to disperse the gas).
2. Gas Cylinder Pressure Gauge reads in the YELLOW	There is not enough gas pressure to start a procedure.	Either change the cylinders, or, if the cylinder was stored in a cold area, allow the cylinder to warm up to ambient room temperature.
3. Insufficient Freeze to Probe Tip	a. Gas leaking from device b. Non-Medical grade gas used c. Functional failure d. Repair done by unqualified personnel e. Operating in the yellow zone on the pressure gauge	a-d. Refer to the service and repair section located in the instruction manual.  e. Replace tank
4. Pressure Gauge Failure	Pressure lines clogged inside console.	Purge the gas lines to make sure that they are clear of moisture or foreign matter.
5. No Timer/ Temperature Display	a. Probe is not equipped with an integral thermocouple b. Low batteries	a. Use a probe that is equipped with a thermocouple b. Replace batteries

**INSTRUCTIONS FOR CLEANING LL100™**

The Wallach LL100 Cryosurgical Freezer can be disinfected by wiping down with commercial disinfectant. This wipe down procedure can include the black plastic handle and white silicone gas hoses. Other methods involving heat, fluid (cold soaking) or steam are not recommended. The Cryofreezer holder must be sterilized after each procedure.

**CARE of CONSOLE**

Consoles must never be sterilized. Consoles can be wiped down with an appropriate disinfectant. Prevent fluids from accumulating on the surface, or from spilling into the internal mechanism.

**INSTRUCTIONS FOR HIGH LEVEL DISINFECTION of INSULATION PROBE**

The white insulation probe which accepts the cryosurgical tip could have mucous membrane contact and should, therefore, be disinfected “by wrapping with a cloth soaked in a high-level disinfectant to allow the recommended contact time. After disinfection, the probe should be rinsed with tap water and dried before use. As with other high-level disinfection procedures, proper cleaning of probes is necessary to ensure the success of the subsequent disinfection.”<sup>1</sup>

**INSTRUCTIONS FOR CLEANING STERILIZABLE CRYOTIPS**

1. Remove disposable shield from cryotip and discard.
2. Plug cryotip before cleaning and sterilizing. (Plug Orifice does not require sterilization).
3. Scrub tip to remove blood and particles.
4. Rinse with water.
5. Cryotips can be steam sterilized by Gravity or Pre-vacuum using the following procedures:

	Stage / Process	Parameter Set-Points	
<b>Set-up</b>	Cycle Profile	<b>Gravity</b>	<b>Pre-vacuum</b>
<b>Pre-Exposure</b> (Conditioning)	Chamber Purge	1 minute	1 minute
<b>Exposure</b> (Steam Injection)	Sterilization Temperature	132 °C (270 F)	132 °C (270 F)
	Sterilization Time	15 minutes	4 minutes
<b>Post-Exposure</b> (Exhaust / Dry)	Dry Time	15 minutes	20 minutes

6. Attach new disposable plastic shield to the cryotip, ensuring that the tabs on the shields are correctly fitted into the grooves on the metal tip.
7. Cryotip is now ready for use.

**WARNING: DO NOT USE ETO, E-BEAM OR GAMMA STERILIZATION.**

<sup>1</sup>Rutala, William A., Ph.D., M.P.H., David J. Weber, M.D., M.P.H., and the Healthcare Infection Control Practices Advisory Committee (HICPAC), *CDC Guideline for Disinfection and Sterilization in Healthcare Facilities*, pg. 19, 2008.

**WARRANTY**

The WA1500 (Mobile Cart), the WA2000 (Timer/Temperature Indicator Console) and the WA4000 (General Purpose Probe Console) each carry a one-year warranty. **DO NOT ATTEMPT IN-HOUSE REPAIRS AS THEY WILL VOID YOUR WARRANTY.**

**SERVICE AND REPAIR**

In case of malfunction, immediately turn gas cylinder master valve “OFF”. Detach instrument from gas cylinder. Call the factory at 203-799-2000 and ask the Service Department for assistance. If a repair is needed, carefully sanitize and package the unit in a protective carton. Items that are not sanitized will be returned to the customer freight collect. Include a note stating the problem and ask for a repair estimate. All shipments must be made via prepaid parcel post or U.S. Mail. C.O.D. packages will not be accepted. Address package to the attention of the Repair Department at the address below:

**WALLACH®**

95 Corporate Drive  
 Trumbull, CT 06611 USA  
 Phone: 800-243-2463  
 (203) 799-2000  
 Fax: (203) 799-2002  
 www.wallachsurgical.com

**EXPLANATION OF SYMBOLS**

**REF** Reorder Number

**SN** Serial Number

 Consult instructions for use

 Caution

 Air Flow (Gas Exhaust Port)

**F**  
**R** Activate Freeze  
**Z**

**D**  
**E** Activate Defrost  
**F**

**TC1** Thermocouple-1  
**TC2** Thermocouple-2

 Type B Applied Part

 For professional use only

 System OFF

 System ON

 Positive Terminal

 Negative Terminal

 Direct Current

**R<sub>x</sub> Only** U.S. Federal law restricts this device to sale by or on the order of a physician.

**CE0086** Product conforms to the Medical Device Directive 93/42/EEC

**EC REP** Authorized Representative in the European Community.

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