


TOTAL ABSCESSION[®]

Drainage Catheter




A NEW LOCK ON DRAINAGE CATHETERS



ANGIODYNAMICS[®]

INCORPORATED

TOTAL ABSCSSION® Drainage Catheter



TOTAL Security

The VAULT® locking mechanism securely fixes the pigtail and prevents tampering.

TOTAL Versatility

Aspiration can occur in the locked or unlocked position.

TOTAL Performance

"Kink proof" shaft recovers rapidly, even if severely bent, knotted, or twisted.

Ordering Information

Catalog #	Description	Qty/Box
14000801	6.5F x 15cm TOTAL ABSCSSION Drainage Catheter	5
14000802	6.5F x 30cm TOTAL ABSCSSION Drainage Catheter	5
14000803	8F x 30cm TOTAL ABSCSSION Drainage Catheter	5
14000804	10F x 30cm TOTAL ABSCSSION Drainage Catheter	5
14000805	12F x 30cm TOTAL ABSCSSION Drainage Catheter	5
14000806	14F x 30cm TOTAL ABSCSSION Drainage Catheter	5
14000807	8F x 45cm Total Abscession Drainage Catheter	5
14000808	10F x 45cm Total Abscession Drainage Catheter	5
14000809	12F x 45cm Total Abscession Drainage Catheter	5
14000810	14F x 45cm Total Abscession Drainage Catheter	5

ANGIODYNAMICS®
INCORPORATED

603 Queensbury Ave., Queensbury, NY 12804 USA
Phone: 518-798-1215 Fax: 518-798-1360 Toll Free USA: 1-800-772-6446
www.angiodynamics.com

AngioDynamics, Abscession and Vault are registered trademarks of AngioDynamics, Inc. Covered under US Patent 6,000,000

MLC 104 REV F 07/07

