

Catalog Contents

LeMaitre Valvulotomes	2
VascuTape Radiopaque Marking Tape	3
LeMaitre Catheters	4–5
Pruitt-Inahara Carotid Shunt	6
Pruitt Occlusion Catheters	7
Reddick Cholangiogram Catheter	8
Laparoscopic Instruments	9
Information Materials	10
Return Policy	12
Order Information	IBC



	CATALOG #	PRICE
Expandable LeMaitre® Valvulotome (Sterile)		
NEW 1.8 mm Expandable LeMaitre Valvulotome (blade range 1.8-6.0 mm) Includes: one Expandable LeMaitre Valvulotome and two strips Glow 'N Tell Tape	1000-90	\$525
NEW 1.8 mm Expandable LeMaitre Valvulotome 5 Pack (blade range 1.8-6.0 mm) Includes: five Expandable LeMaitre Valvulotomes and ten strips Glow 'N Tell Tape	1000-95	\$2,395
2.0 mm Expandable LeMaitre Valvulotome (blade range 2.0-5.0 mm) Includes: one Expandable LeMaitre Valvulotome and two strips Glow 'N Tell Tape	1000-80	\$495
2.0 mm Expandable LeMaitre Valvulotome 5 Pack (blade range 2.0-5.0 mm) Includes: five Expandable LeMaitre Valvulotomes and ten strips Glow 'N Tell Tape	1000-85	\$2,195
Disposable LeMaitre® Valvulotome (Sterile)		
1.5 mm Disposable LeMaitre Valvulotome	1001-15	\$330
2.0 mm Disposable LeMaitre Valvulotome	1001-20	\$330
2.5 mm Disposable LeMaitre Valvulotome	1001-25	\$330
3.0 mm Disposable LeMaitre Valvulotome	1001-30	\$330
3.5 mm Disposable LeMaitre Valvulotome	1001-35	\$330
4.0 mm Disposable LeMaitre Valvulotome	1001-40	\$330
4.5 mm Disposable LeMaitre Valvulotome	1001-45	\$330
Reusable LeMaitre® Valvulotome (Non-Sterile)		
1.5 mm Reusable LeMaitre Valvulotome	1002-15	\$950
2.0 mm Reusable LeMaitre Valvulotome	1002-20	\$950
2.5 mm Reusable LeMaitre Valvulotome	1002-25	\$950
3.0 mm Reusable LeMaitre Valvulotome	1002-30	\$950
3.5 mm Reusable LeMaitre Valvulotome	1002-35	\$950
4.0 mm Reusable LeMaitre Valvulotome	1002-40	\$950
4.5 mm Reusable LeMaitre Valvulotome	1002-45	\$950
LeMaitre Sterilization/Storage Tray	1002-00	\$320
Reusable LeMaitre® Valvulotome Systems (Non-Sterile)		
Reusable LeMaitre Deluxe System (2.0, 2.5, 3.0, 3.5, 4.0 w/Tray)	1002-80	\$3,700
Reusable LeMaitre Deluxe System (2.0, 2.5, 3.0, 3.5, 4.0)	1002-81	\$3,500
Reusable LeMaitre Starter System (2.0, 3.0, 4.0 w/Tray)	1002-82	\$2,750
Reusable LeMaitre Starter System (2.0, 3.0, 4.0)	1002-83	\$2,500

	CATALOG #	PRICE
Glow 'N Tell® Tape (Sterile)		
Glow 'N Tell Tape (20 Strips)	1100-20	\$360
Glow 'N Tell Tape (50 Strips)	1100-50	\$750
Glow 'N Tell Tape (100 Strips)	1100-00	\$1,295
LeMaitre® Stent Guide (Sterile)		
LeMaitre Stent Guide (20 Strips)	1102-20	\$360
LeMaitre Stent Guide (50 Strips)	1102-50	\$750
LeMaitre Stent Guide (100 Strips)	1102-00	\$1,295
Glow 'N Tell® Tape for Dorsal Placement (Sterile)		
Glow 'N Tell Tape for Dorsal Placement (20 Strips)	1103-20	\$360
Glow 'N Tell Tape for Dorsal Placement (50 Strips)	1103-50	\$750
Glow 'N Tell Tape for Dorsal Placement (100 Strips)	1103-00	\$1,295
Mini Glow 'N Tell® Tape (Sterile)		
Mini Glow 'N Tell Tape (50 Strips)	1105-50	\$215
LeMaitre® Stent Graft Guide (Sterile)		
LeMaitre Stent Graft Guide (10 Strips)	1104-10	\$300

BALLOON Diameter	BALLOON Volume	DIAMETER	LENGTH	CATALOG #	PRICE
Single Lumen Embolectomy Catheter (Sterile)					
4.5 mm	0.05 ml	2F	40 cm	1601-24	\$65
4.5 mm	0.05 ml	2F	60 cm	1601-26	\$65
8.0 mm	0.20 ml	3F	40 cm	1601-34	\$40
8.0 mm	0.20 ml	3F	80 cm	1601-38	\$40
10.5 mm	0.75 ml	4F	40 cm	1601-44	\$40
10.5 mm	0.75 ml	4F	80 cm	1601-48	\$40
13.0 mm	1.50 ml	5F	40 cm	1601-54	\$40
13.0 mm	1.50 ml	5F	80 cm	1601-58	\$40
13.5 mm	1.60 ml	6F	40 cm	1601-64	\$40
13.5 mm	1.60 ml	6F	80 cm	1601-68	\$40
14.0 mm	1.75 ml	7F	40 cm	1601-74	\$40
14.0 mm	1.75 ml	7F	80 cm	1601-78	\$40

Single Lumen Thrombectomy Catheter (Sterile)

14.5 mm	2.00 ml	6F	80 cm	1602-68	\$50
17.5 mm	3.00 ml	8F	80 cm	1602-88	\$50

Single Lumen Biliary Catheter (Sterile)

13.0 mm	1.50 ml	5F	30 cm	1603-53	\$40
13.5 mm	1.60 ml	6F	30 cm	1603-63	\$40

Single Lumen Irrigation Catheter (Sterile)

4F	80 cm	1607-48	\$17
6F	80 cm	1607-68	\$17
10F	40 cm	1607-04	\$17

	BALLOON Diameter	BALLOON Volume	DIAMETER	LENGTH	OLD CATALOG #	NEW CATALOG #	PRICE
Dual Lumen Distal Irrigation Catheter (Sterile) With over-the-guidewire capabilities on 4F, 5F, 6F, 7F							
	6 mm	0.20 ml	3F	40 cm	DL-3F-40D	2302-34	\$78
	6 mm	0.20 ml	3F	80 cm	DL-3F-80D	2302-38	\$78
	9 mm	0.75 ml	4F	40 cm	DL-4F-40D	2302-44	\$78
	9 mm	0.75 ml	4F	80 cm	DL-4F-80D	2302-48	\$78
	11 mm	1.50 ml	5F	40 cm	DL-5F-40D	2302-54	\$78
	11 mm	1.50 ml	5F	80 cm	DL-5F-80D	2302-58	\$78
	13 mm	2.00 ml	6F	40 cm	DL-6F-40D	2302-64	\$78
Stopcock and syringe included	13 mm	2.00 ml	6F	80 cm	DL-6F-80D	2302-68	\$78
Stylet not included on 40 cm length	14 mm	2.50 ml	7F	40 cm	DL-7F-40D	2302-74	\$78
Guidewires available for dual lumen distal models	14 mm	2.50 ml	7F	80 cm	DL-7F-80D	2302-78	\$78

Guidewires (Sterile)

For 4F/5F	dual lumen distal models	.018"	125 cm	GW-80	2302-80	\$10
For 6F/7F	dual lumen distal models	.028"	125 cm	GW-028	2302-81	\$10

Dual Lumen Proximal Irrigation Catheter (Sterile)

	6 mm	0.20 ml	3F	40 cm	3F-IE-40	2303-34	\$78
	6 mm	0.20 ml	3F	80 cm	3F-IE-85	2303-38	\$78
	9 mm	0.75 ml	4F	40 cm	4F-IE-40	2303-44	\$78
	9 mm	0.75 ml	4F	80 cm	4F-IE-85	2303-48	\$78
	11 mm	1.50 ml	5F	40 cm	5F-IE-40	2303-54	\$78
	11 mm	1.50 ml	5F	80 cm	5F-IE-85	2303-58	\$78
	13 mm	2.00 ml	6F	80 cm	6F-IE-85	2303-68	\$78
Stopcock and syringe included	14 mm	2.50 ml	7F	80 cm	7F-IE-85	2303-78	\$78

All single and dual lumen catheters packaged in tube-style packaging.

Pruitt-Inahara® Carotid Shunt

	DIAMETER	LENGTH	OLD CATALOG #	NEW CATALOG #	PRICE
Pruitt-Inahara® Outlying Shunt (Sterile)					
Outlying Carotid Shunt without T-Port	9F	31 cm	100-10-9F	2000-19	\$159
Outlying Carotid Shunt with T-Port	9F	31 cm	400-40-9F	2000-49	\$159
Outlying Carotid Shunt with T-Port	8F	25 cm	400-40-8F	2000-48	\$159
Inahara-Pruitt® Inlying Shunt (Sterile)					
Inlying Carotid Shunt without T-Port	9F	15 cm	200-20-9F	2000-29	\$159
Inlying Carotid Shunt with T-Port	9F	15 cm	500-50-9F	2000-59	\$159
Inlying Carotid Shunt with T-Port	8F	15 cm	500-50-8F	2000-58	\$159
Saf-T-Spand® Safety Balloon Occlusion Catheter (Sterile)					
	5F	27 cm	SOC-5F	2001-05	\$65

Pruitt Occlusion Catheters

	BALLOON Diameter	DIAMETER	LENGTH	OLD CATALOG #	NEW CATALOG #	PRICE
Pruitt Aortic Occlusion Catheter (Sterile)						
Pruitt Aortic Occlusion Catheter (Sterile)		12F	24 cm	12F-10	2100-12	\$119
Pruitt Aortic Occlusion Catheter and Kit (Sterile)		12F	24 cm	12F-RPK	2101-12	\$143
Kit includes: aspiration needle, guidewire, dilator, peel away sheath & syringe						

Pruitt Irrigation Occlusion Catheter (Sterile)

	9F	23 cm	9F-10	2102-09	\$84
--	----	-------	-------	---------	------

Pruitt Occlusion Catheter (Sterile)

7 mm	3F	27 cm	3F-6	2103-36	\$59
9 mm	4F	27 cm	4F-6	2103-46	\$59
11 mm	5F	27 cm	5F-6	2103-56	\$59

Stopcock and syringe included with all models

Hutson Dual Balloon Occlusion Catheter (Sterile)

9 mm	4F	40 cm	H-4F-DB	2104-04	\$50
11 mm	5F	40 cm	H-5F-DB	2104-05	\$50
13 mm	6F	40 cm	H-6F-DB	2104-06	\$50

Vascular Occluders (Sterile)

1.00 mm	75 mm	OCC100	2201-00	\$265
1.25 mm	75 mm	OCC125	2201-25	\$265
1.50 mm	75 mm	OCC150	2201-50	\$265
1.75 mm	75 mm	OCC175	2201-75	\$265
2.00 mm	75 mm	OCC200	2202-00	\$265
2.25 mm	75 mm	OCC225	2202-25	\$265
2.50 mm	75 mm	OCC250	2202-50	\$265
2.75 mm	75 mm	OCC275	2202-75	\$265
3.00 mm	75 mm	OCC300	2203-00	\$265
3.50 mm	75 mm	OCC350	2203-50	\$265
4.00 mm	75 mm	OCC400	2204-00	\$265
4.50 mm	75 mm	OCC450	2204-50	\$265

5 Packs per box

Pack (3 occluders): 2 occluders at stated size and 1 occluder of next larger size.

Reddick Cholangiogram Catheter

	DIAMETER	LENGTH	OLD CATALOG #	NEW CATALOG #	PRICE
<hr/>					
Reddick Cholangiogram Catheter (Sterile) 5 per box					
Reddick Cholangiogram Catheter	4F	50 cm	RCD-50	2400-50	\$349
Reddick Cholangiogram Stiffer Catheter	4F	50 cm	RCD-50/S	2400-51	\$349
Reddick Cholangiogram Catheter with needle introducer	4F	50 cm	RCD-50/N	2400-52	\$399
Reddick Cholangiogram Stiffer Catheter with needle introducer	4F	50 cm	RCD-50/SN	2400-53	\$399

Laparoscopic Instruments

	DIAMETER	LENGTH	OLD CATALOG #	NEW CATALOG #	PRICE
<hr/>					
Reddick-Saye Screw Retractor Kit (Sterile)					
1 Reddick-Saye Screw and 1 13G aspiration needle	5 mm	30 cm	R/S KIT	2500-01	\$49
<hr/>					
Grice® Laparoscopic Suture Needle (Sterile)					
Slotted trocar needle introducer		20 cm	GSN-3	2500-02	\$44

Information Materials

Reversed & In Situ Audiotape (1998 Edition)	A1000-02	FREE
LeMaitre Vascular Resources CD-Rom (2001 Edition)	A1000-03	FREE
In Situ Bypass Grafting Video (2001 Edition)	A1000-06	FREE
Peripheral Vascular Disease Patient Brochure (2001 Edition)	A1000-07	FREE
In Situ Bypass Grafting Booklet (2001 Edition)	A1000-08	FREE
A Crusade for Stroke Prevention Book (2001 Edition)	A1000-09	FREE

From Invoice Date:

1 – 30 Days Any unopened, sterile or non-sterile product will be accepted for any reason for full credit (excluding freight).

Over 30 Days Any unopened, sterile or non-sterile product, subject to approval, 20% re-stocking fee, and the following restrictions.

The following merchandise is not acceptable for return:

- Sterile items that have less than 18 months of shelf life
- Sterile items that are returned without factory seals intact
- Incomplete kits or cases
- Items with hospital stickers or markings
- Items missing manuals or accessories

Authorization is required for all returns. Authorization may be obtained through the customer service department. A Return Goods Authorization (RGA) number will be issued. This number must be referenced on all packages.

Address for Returned Goods:

LeMaitre Vascular, Inc.

26 Ray Avenue

Burlington, MA 01803

Attention: Returned Goods RGA# _____