

**NEW
Latex-
Free**

NovaSil

Latex-Free Embolectomy Catheters


LeMaitre Vascular now offers a new line of premium quality silicone balloon catheters to provide a latex-free option for latex sensitivity.

NovaSil™ catheters offer the following benefits:

- **Latex-Free** — Silicone balloon gives a latex-free option for latex sensitive patients.
- **Concentricity** — Well-centered balloons, reducing intimal damage.
- **Quality** — Performance and reliability to facilitate effective clot removal.

NovaSil: The Solution to Latex Sensitivity


6%



Estimates of latex sensitivity in the general population range up to 6%.

Over 17 million Americans are potentially affected.

The National Institute of Allergy and Infectious Diseases.
National Institute of Health. January 2002.

	Balloon Diameter	Balloon Volume	Diameter	Length	Catalog #
NovaSil™ Latex-Free Single Lumen Embolectomy Catheter (Sterile)					
					
	4 mm	0.05 ml	2F	40 cm	1801-24
	4 mm	0.05 ml	2F	60 cm	1801-26
	6 mm	0.20 ml	3F	40 cm	1801-34
	6 mm	0.20 ml	3F	80 cm	1801-38
	9 mm	0.60 ml	4F	40 cm	1801-44
	9 mm	0.60 ml	4F	80 cm	1801-48
	11 mm	1.00 ml	5F	40 cm	1801-54
	11 mm	1.00 ml	5F	80 cm	1801-58
	13 mm	1.60 ml	6F	80 cm	1801-68
	14 mm	1.75 ml	7F	80 cm	1801-78
NovaSil™ Latex-Free Single Lumen Biliary Catheter (Sterile)					
	11 mm	1.00 ml	5F	30 cm	1803-53
	13 mm	1.60 ml	6F	30 cm	1803-63

All products packaged in tube-style packaging.




For more information, contact your local
LeMaitre Vascular Sales Representative.

For Customer Service:

Tel: 781 221-2266

Fax: 781 221-2223



SURGEONS' COOPERATIVE®
Largely owned by vascular surgeons, LeMaitre Vascular strives to continue the work of its founder, George D. LeMaitre, MD. The company is dedicated to providing vascular surgeons and interventional radiologists with innovative devices and information to improve vascular procedure efficacy and reduce patient recovery time.

LeMaitre®
VASCULAR

LeMaitre Vascular, Inc.
26 Ray Avenue
Burlington, MA 01803
USA
Tel: 781 221-2266
Fax: 781 221-2223

www.lemaitre.com

LeMaitre Vascular GmbH
Wiesbadener Weg 2a
65812 Bad Soden / Ts.
Germany
Tel: 49-(0)6196-659230
Fax: 49-(0)6196-527072