


# VascuTape® Radiopaque Tape

Precision In Vascular Imaging




 **LeMaitre®**  
VASCULAR

More Than The Valvulotome.

**Precise Vascular Imaging  
For Accuracy and Ease:**

# VascuTape® Radiopaque Tape

VascuTape® — the radiopaque marking tape from LeMaitre® Vascular, provides pinpoint precision during vascular procedures. Stent placements, fistula locations during in situ bypass grafting, and balloon catheter insertions are all made more simple and accurate through the use of these innovative imaging tapes. Designed by a practicing vascular surgeon, VascuTape is acclaimed for its accuracy, ease and convenience.



## **Glow 'N Tell® Tape**

Glow 'N Tell Tape (30 cm) is LeMaitre Vascular's original tape product. It has centimeter markings and is ideal for finding fistulas during in situ bypass grafting.

## LeMaitre Stent Guide


LeMaitre Stent Guide (270 mm) has millimeter markings for greater precision and accuracy in sizing and placing stents.



## Choose VascaTape to Facilitate:

- > In Situ Bypass Grafting
- > Stent Placement
- > Balloon Angioplasty
- > Vena Cava Filter Placement
- > Stent Grafting
- > Atherectomy/Plaque Excision

LeMaitre Stent Guide  
240  
270

A roll of VascaTape is shown in the bottom left corner. The tape is white with black millimeter markings and the number 60 printed on it. The roll is partially unrolled, showing the adhesive side.

## How To Use VascaTape

Apply the adhesive sterile tape to the patient's skin before the radiology procedure. The radiopaque markings — visible on the fluoroscope and on the patient's skin — provide precise, stable cross-referencing. LeMaitre Vascular's innovative marking tapes noticeably improve the way interventional procedures are performed.

- > Flexible and Adhesive
- > Easy-To-Use
- > Minimally Invasive
- > Faster Procedure
- > Reduced Exposure to Contrast Medium

	USABLE LENGTH	QUANTITY	MODEL #
<b>GLOW 'N TELL® TAPE</b>	30 cm	20 strips	1100-20
	30 cm	50 strips	1100-50
	30 cm	100 strips	1100-00
<b>LEMAITRE® STENT GUIDE</b>	270 mm	20 strips	1102-20
	270 mm	50 strips	1102-50
	270 mm	100 strips	1102-00

These specifications are not intended as a warranty. In the interest of product improvement, these specifications may be changed from time to time without notice. Please consult your sales representative for details.



More Than The Valvulotome.

LeMaitre Vascular, Inc.  
63 Second Avenue  
Burlington, MA 01803  
USA

Tel: 781 221-2266  
Fax: 781 221-2223

Web: [www.lemaitre.com](http://www.lemaitre.com)

LeMaitre Vascular GmbH  
Otto-Volger-Str. 5a/b  
65843 Sulzbach/Ts.  
Germany

Tel: +49-(0)6196-659230  
Fax: +49-(0)6196-527072

LeMaitre Vascular KK  
1F Kubodera Twin Tower Bldg.  
2-9-4 Kudan-minami, Chiyoda-ku  
Tokyo, 102-0074  
Japan

Tel: +81-(0)3-5215-5681  
Fax: +81-(0)3-5215-5682

For more information,  
please contact your local  
LeMaitre Vascular Sales  
Representative.

For Customer Service:

Tel: 800 628-9470



LeMaitre, Glow 'N Tell, and VasculTape are registered trademarks of LeMaitre Vascular, Inc.  
©2006 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.



PO Box 910  
Saint Just, PR-USA 00978-0910  
Tel. (787) 761-0000  
Fax (787) 761-0404  
[www.pelegrina.com](http://www.pelegrina.com)

**PELEGRINA MEDICAL, INC.**  
**PELEGRINA EQUIPMENT, CORP.**