

CLEARTONE®

Handheld Doppler Ultrasound System

The Clear Choice for Superior Sound Quality

- Advanced audio output
 - Offers superior volume control
 - Provides clear static-free sound quality
 - Includes a 3.5mm audio jack for headset connection
- Flexible and Versatile
 - Interchangeable fetal and vascular probes
 - Optional waterproof probe
 - All probes are connected via a replaceable cable
- Excellent value
 - A durable and elegant doppler at an affordable price
 - Different models to choose from
 - Long lasting battery life
 - 5 year parts warranty
- Optional digital display
 - Easy-to-read Digital Heart Rate display calculates beats per minute to save counting and calculating time
 - Displays heart rate within 2 – 4 seconds
 - Signal strength indicator
 - Monitors battery life and battery charging (with rechargeable version)



CLEARTONE®

Doppler Ultrasound System

ClearTone Specifications

- Fetal Doppler: 2.1 MHz and 3.2 MHz
- Waterproof Fetal Doppler: 2.1 MHz
- Vascular Doppler: 4, 5 and 8 MHz (not for fetal use)

Outputs

One 3.5 mm jack to connect to a headset or audio recording device

Controls

Separate volume control and push ON/OFF button

Power Source

- Battery voltage 3 - AA Alkaline 1.5 volt (non-rechargeable) or 3 - AA NMHi 1.2 volt (rechargeable)
- Battery life: 1,250 one minute exam (Alkaline) or 1,000 one minute exams (NMHi) rechargeable

Weight

11.29 oz. (320 gm)

Dimensions

5.5 x 2.75 x 1.3 in (140 x 70 x 35 mm)

Service

CooperSurgical, Inc. maintains a full-service facility that can promptly repair products returned to the factory

Return/Warranty

- CooperSurgical, Inc. is pleased to present our customers with a 30-day no hassle evaluation return policy in the event you are not satisfied with the product
- Dopplers are covered by a 5-year warranty on parts, 1-year warranty on labor

The ClearTone Doppler comes with one probe and ultrasound gel

Order Information

Frequency	Application	Order Description	Number
2 MHz	For Labor and Delivery and Later Gestational Use	ClearTone Non-Display with 2 MHz probe	CT150-2
		ClearTone Non-Display with 2 MHz probe - Rechargeable	CT150R-2
2 MHz	Waterproof	ClearTone with Display with 2 MHz probe	CT250-2
		ClearTone with Display with 2 MHz probe - Rechargeable	CT250R-2
		ClearTone Non-Display with 2 MHz probe - Waterproof	CT150-2WP
		ClearTone Non-Display with 2 MHz probe - Rechargeable - Waterproof	CT150R-2WP
2 MHz	Waterproof	ClearTone with Display with 2 MHz probe - Waterproof	CT250-2WP
		ClearTone with Display with 2 MHz probe - Rechargeable - Waterproof	CT250R-2WP
3 MHz	Versatile OB/GYN and Early Gestation	ClearTone Non-Display with 3 MHz probe	CT150-3
		ClearTone Non-Display with 3 MHz probe - Rechargeable	CT150R-3
3 MHz	Versatile OB/GYN and Early Gestation	ClearTone with Display with 3 MHz probe	CT250-3
		ClearTone with Display with 3 MHz probe - Rechargeable	CT250R-3
4 MHz	Vascular Use	ClearTone Non-Display with 5 MHz probe	CT150-5
5 MHz	Designed for Obese Patients	ClearTone Non-Display with 5 MHz probe - Rechargeable	CT150R-5
		ClearTone with Display with 5 MHz probe	CT250-5
5 MHz	Designed for Obese Patients	ClearTone with Display with 5 MHz probe - Rechargeable	CT250R-5
8 MHz	Vascular Use	ClearTone Non-Display with 8 MHz probe	CT150-8
		ClearTone Non-Display with 8 MHz probe - Rechargeable	CT150R-8
		ClearTone with Display with 8 MHz probe	CT250-8
		ClearTone with Display with 8 MHz probe - Rechargeable	CT250R-8
Miscellaneous	Miscellaneous	ClearTone 2 MHz probe for L&D and late gestations	CTP-2
		ClearTone 2 MHz probe - Waterproof	CTP-2
		ClearTone 3 MHz for OB/GYN office and early gestations	CTP-3
		ClearTone 5 MHz probe for Vascular use	CTP-5
		ClearTone 8 MHz probe for Vascular use	CTP-8



Optional Rollstand



Non-Display Version



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products which aid clinicians in the management and treatment of commonly seen conditions.

203.601.5200 | 800.243.2974 | coopersurgical.com | 95 Corporate Dr. | Trumbull, CT 06611
Also on YouTube at www.YouTube.com/CooperSurgical

ClearTone® is a registered trademark of CSI.
© 2011 CooperSurgical, Inc.

Form #81134 • 06/11