

# Marlow Balloon Cannula With Atraumatic Surface Disc

*More control for you.  
Less trauma for your Patient.*

- *Elliptical balloon for internal stabilization*
- *Atraumatic surface disc for external stabilization*
- *Low profile for increased working length*
- *Floating trocar tip for reduced tissue trauma*



# Marlow Balloon Cannula

***More Control for You. Less Trauma for Your Patient.***

*Here's a device that provides ultimate performance; enhanced surgeon control and reduced patient trauma. The uniquely designed Marlow Balloon Cannula incorporates many benefits for you and your patient:*

## **Elliptical Balloon**

The exclusive elliptical shape of the balloon assures that no dislocation will occur during instrument exchanges, thereby providing an internal stabilization system. It also permits a tamponade effect over a large surface area minimizing tissue trauma.

## **Atraumatic Surface Disk**

Once the cannula is in place, the disk is moved down the cannula, gently pressing the skin surface to aid in external stabilization. The combined advantages of the elliptical balloon and the atraumatic surface disk configuration create a stable entry port which delivers an excellent tamponade effect to stop bleeding while reducing the loss of pneumoperitoneum at the port site.

## **Floating Trocar Tip**

The tip remains stationary as the surgeon rotates the unit during insertion. This provides the surgeon with more control and reduced tissue trauma.

## **Built in Flap Reducer**

Accommodates 3, 4 and 5mm instruments. The tight seal prevents gas loss and provides for an easy exchange of instruments.

## **Low Profile Cannula Length**

The cannula is shorter than standard products, so there is minimal cannula length inside the abdomen. This improves visibility and increases instrument working length for deep areas.

## **Pre-Packaged Sharp and Blunt Obturators**

This unique approach reduces inventory demands while providing you with complete flexibility in choice of entry technique.

*The more instrument exchanges you do, the more you need the Marlow Balloon Cannula.  
Discover the advantages.*

## **ORDERING INFORMATION**

<b>Product #ML0500D</b>	<b>Balloon Cannula, 5mm, Box of 12, Sharp Trocar Included</b>
<b>Product #ML1012D</b>	<b>Balloon Cannula, 10/12mm, Box of 12, Sharp Trocar and Blunt Obturator Included</b>
<b>Product #ML1012LD</b>	<b>Balloon Cannula, 10/12mm Long, Box of 12, Sharp Trocar and Blunt Obturator Included</b>

**CooperSurgical**

15 FOREST PARKWAY  
SHELTON, CT 06484

phone: 800-992-5581  
203-929-6321  
fax: 203-926-9435



**International Division**  
phone: 203-929-0640  
fax: 203-926-9496

**...HELPING TO KEEP YOU  
AT THE FOREFRONT OF  
WOMEN'S HEALTH CARE™**