

NEO-fit™

For easy, secure endotracheal tube stabilization

Key Features & Benefits

- Unique design eliminates the use of tape on the endotracheal tube and minimizes the use of tape on the patient
- Easy to use; "one size fits all"
- Demarcations are clearly seen on tube
- Hook-and-loop strap allows rapid release and refit, making tube adjustment simple
- Longer term use; cost-efficient



NEO-fit™

For easy, secure endotracheal tube stabilization

The NEO-fit provides quick, secure endotracheal tube stabilization. Designed to firmly grip a neonatal endotracheal tube, NEO-fit utilizes metal clips to mechanically grip the soft material of the tube. The device is positioned on baby's face using a hydrocolloid adhesive that is hypoallergenic; VELCRO® strap allows for rapid release and readjustment of the tube.

Product Ordering information

NEO-fit includes one plastic tube grip with adjustable hook-and-loop strap, two adhesive foam pads and two foam "lollypop" securing tapes (optional use).

Product Number

42-2540

Packaging

Nonsterile, 20 per box, 80 per case (4 boxes)

Complementary Products:

NEO-prep

Provides protection for the tender neonatal septum

NEO-prep provides an extra skin barrier, minimizing septal breakdown, and can be used with CooperSurgical's NEO-fit and INCA products. The precut hypoallergenic skin barrier is available in either small or large size.

Product Ordering information

Code	Size
Neo-prep	
44-8000	Small
44-8001	Large

Packaging

8 per strip, 10 strips per box, 3 boxes (240) per case



CooperSurgical

CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products which aid clinicians in the management and treatment of commonly seen conditions.

203.601.5200 | 800.243.2974 | coopersurgical.com | 95 Corporate Dr. | Trumbull, CT 06611

©2012 CooperSurgical, Inc.

Form #81564 • 11/12