

# EDAN



## iM20 Patient Monitor

Serving as both multi-parameter module of elite V8 and transport monitor, iM20 can provide a comprehensive and seamless monitoring solution both inside and outside hospital.



Waterproof  
Dustproof



Portable  
Design



Ambulance



Emergency

Transport

Care for Health

# iM20

## Patient Monitor

- 5" high resolution color TFT-LCD with full-touch screen
- No fan design
- 12 volt direct power supply with built-in high-capacity rechargeable battery
- iSEAP™ ECG algorithm optimized for arrhythmia detection, pacemaker detection, and HR measurement
- SEMIP 12-lead diagnosis algorithm with 208 kinds of analysis results
- iMAT™ SpO<sub>2</sub> algorithm with outstanding motion resistance and low perfusion resistance performance
- iCUFS™ NIBP algorithm optimized for cardiac patients, hypertensive patients, and neonatal patients
- SunTech Transport Motion Tolerant NIBP to address noise and motion associated with transporting (optional)
- Nellcor SpO<sub>2</sub>, 12-lead ECG, End-tidal CO<sub>2</sub> and Invasive blood pressure are available as options
- Built-in Wi-Fi/3G supports seamless data transmission from ambulatory cares to hospitalization

**1200**

NIBP Measurement  
Review

**120 h**

Trend Review

**200**

Alarm Review



**MFM-CMS**

**Central Monitoring System**



**EDAN** Edan Instruments, Inc.

■ 3/F – B, Nanshan Medical Equipment Park, Nanshan Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China  
 ■ Tel: +86-755-26898326 ■ Fax: +86-755-26898330 ■ www.edan.com.cn ■ Email: info@edan.com.cn

■ All rights reserved. This catalogue is only a temporary version. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.

ENG-PM-iM20-V0.1-20131112