



Non-Contact Thermometer 433

Adtemp 433 Non-Contact Thermometer

ADC's newest non-contact thermometer, the Adtemp 433, features a trigger-style design and fast measurements. It is perfect for rapid screening in medical environments, offices, restaurants, and more.

- Accurate readings in less than two seconds
- Trigger-style design for rapid screenings at the forehead
- No skin contact required, minimizes the risk of cross contamination
- Color-coded LCD display with fever warning
- Audible beep when measurement is complete
- 50-reading memory
- °F/°C switchable display
- Uses 2 AAA batteries (included) with battery life of about 2,000 readings
- Measures 5.9" x 1.7" x 3.7" (150 mm x 95 mm x 44 mm), weighs 4.6 oz (130.5 g) with batteries
- Meets ASTM E1965-98 standards for accuracy
- 2-year warranty

Specifications

Type:	Non-Contact
Technology:	Infrared
Scale:	Dual (°F/°C)
Range °F:	93.2° to 109.2°
Range °C:	34.0° to 422.9°
Accuracy °F:	+/- 0.5
Accuracy °C:	+/- 0.3
Response Time:	2 Seconds
Beeper:	Yes
Auto-off:	Yes
Memory:	9x
Flex Probe:	No
Needs Sheaths:	No
Warranty:	2 Years
Battery Type:	1.5V (AAA)
Battery Life:	2,000 measurements
Latex Content:	Latex-free
Packaging:	Giftbox
Dimensions:	5.9" x 1.7" x 3.7" (150 mm x 95 mm x 44 mm)
Weight:	4.6 oz (130.5 g)
Display:	LCD
Operating Environment:	59.0°F to 95.0°F (15.0°C to 35.0°C), ≤85% RH
Storage Environment:	-4.0°F to 131.0°F (-20.0°C to 55.0°C), ≤93%,RH