CoperSurgical











Premier Products for Women's Health Care

Clinic and Practice-Based Products

Contraception Cervical Cryosurgery

Hysteroscopy Ophthalmic Cryosurgery

Endometrial Ablation Pessaries

Urodynamics Obstetrical Care

Colposcopy Fertility and Family Planning

Gynecological Specialty Instruments General Examination: Diagnostic Disposables

Endometrial Samplers Sono HSG

LEEP

Delivering on the Promise to Help Clinicians Provide Better Health Care for Women

Our mission is to enable physicians to provide increasingly more efficient and effective health care for women. We leverage our expertise in identifying and developing a steady stream of exceptional products that keep Ob/Gyns and others directly involved in caring for women, at the forefront of women's health care.

An Expansive Portfolio of Trusted Brands And Services

Since our founding in 1990, CooperSurgical has researched, developed and manufactured a wide range of trusted brands that have advanced the standard of women's health care. Our products are used by physicians and health care professionals in hospitals and their practices for a wide variety of procedures, from basic annual gynecological exams to complex surgical procedures.

As a full-service R&D and manufacturing company, we are committed to:

- Improved patient outcomes and satisfaction
- Clinically proven technologies that provide ease of use
- Innovative solutions that reduce procedural costs
- Educational and in-service support of our products by educated professional representatives
- The highest quality and safety standards in the industry
- Technical support of our equipment − 20+ years of experience
- Annual support of medical professional educational organizations
- Personalized service and a commitment to your satisfaction

We appreciate the opportunity to partner with you to provide optimum health care for women.

<u>CoperSurgical</u>







Creating Opportunities to Grow Your Practice

Traditionally, the Ob/Gyn only served as a woman's primary care physician during the reproductive years. Today, this trend is changing as clinicians expand their service offerings to include screening for and treatment of osteoporosis, diagnosis and treatment of urinary incontinence, as well as correction of abnormal uterine bleeding. By extending the relationship between clinician and patient, CooperSurgical is creating opportunities for physicians to grow their practice.

On the following pages of this catalog, you will find a wide range of high-quality in-office instruments and devices that will provide reliable, safe and effective results for improved patient outcomes.

This catalog is one of four catalogs that represent CooperSurgical's products. Additional catalogs feature: Surgical Products, Labor and Delivery Products and Critical Care Products.

CooperSurgical instruments are engineered for optimal efficiency and effectiveness, producing consistently successful results for the physician.

Clinic and Practice-Based Brands

At CooperSurgical, we develop and market innovative solutions for Ob/Gyn procedures. Within this catalog, you will find an extensive range of instruments and devices that can be used in your practice to improve patient care.

PIPELLE®







EURO-MED®







H/S Elliptosphere™



LEEP PRECISION™

frigitronics°

SELECT MEDICAL SYSTEMS

SpermSelect[®]

Fischer Cone Biopsy Excisor®

In addition to Ob/Gyn clinic and practice-based products, we provide an array of high-quality devices that can be used in gynecologic surgery, for labor and delivery, and critical care. Contact your CooperSurgical sales representative to learn more about our other offerings.

Surgical Brands

RUMI*II SYSTEM The KOH

The KOH Colpotomizer System $^{\scriptscriptstyle{\text{M}}}$





CARTER-THOMASON
CloseSure System

Advincula Arch™

Koh-Efficient®

OR/Mobius°

Kronner Manipujector®

Zeppelin[®]

Advincula **DELINEATOR**™ with Koh-Efficient® Technology

SeeClear*
Laparoscopic Smoke Evacuation Systems

Labor and Delivery Brands

TITY VACE OF THE NEW SYSTEMS

MityOne[®]

OB/Mobius[®]

TransWarmer®

Mystic[®]II

mityhook.

SAFE

Actim PROM

Critical Care Brands

INCA® NEO-fit™

TransWarmer®

Bronchitrac®

COMFIT

Esophageal Balloon Catheter Set

Irrijet®

TRAKE-fit[™]

WarmGel®

Table of Contents:

Contraception PARAGARD® Intrauterine Copper Contraceptive	6
Office Hysteroscopy Endosee®	8
Urodynamics	
Lumax TS Pro™ System	10
LEISEGANG OptiK®	12
OptiK Model 1OptiK Model 2	
Gynecological Specialty Instruments	
Euro-Med® Biopsy Punches Euro-Med® Stainless Steel Specula	18 ככ
Stainless Steel Lateral Wall Retractors	
Stainless Steel Sounds and Dilators	27
Stainless Steel Curettes	28
Stainless Steel Endocervical Specula	
Stainless Steel Forceps	30
Stainless Steel Needle Holders	
Stainless Steel Scissors	
Endometrial Samplers	
Disposable Sound and Dilators	l39
LEEP TO THE STATE OF THE STATE	40
LEEP PRECISION™ Integrated System LEEP Electrodes	
LEEP Procedure Essentials.	
LEEP System Accessories	46
Euro-Med® LEEP Instruments	48
Cervical Cryosurgery LEISEGANG® LM-900 Cryosurgery System	E 1
Frigitronics® Cryo-Plus™	51 52
Cryo-Surg® Handheld System.	52
Ophthalmic Cryosurgery	
Frigitronics® CE-2000™	54
Milex® Pessary Fitting Kit	57
Milex® Pessaries for Prolapse and Incontinence	
Obstetrical Care	
F9 Fetal Monitor	63
ClearTone® Fetal Dopplers	64
Fertility and Family Planning Intrauterine Insemination	67
Fertility Assessment	70 70
Diaphragms	
In-Office Tubal Patency and Uterine Structure Evaluation ABBI®	
General Examination: Diagnostic Disposables	
Catheters	74
Breast Aspiration Dual-Port Biopsy Needle	77 ככ
Endocervical Curette with Vacu-Lok Syringe	<i>11</i> 77
Milex® Cervical Spatulas	77
Disposa-Derm®	77
Dilators Cervical and Endocervical Samplers	
Index/Customer Service	
Warranty Information	
-warrancy information	0 1



CooperSurgical's World Headquarters Trumbull, Connecticut, USA

CooperSurgical established a new world headquarters in a three-building, 214,000 square foot campus in January of 2010. The move was designed to satisfy the need for additional manufacturing and product design space and increased staff. In addition to housing CooperSurgical's main manufacturing facility, the complex serves as our primary distribution center. The company also has specialized manufacturing capabilities in Pasadena, CA, and Stafford, TX, as well as in Berlin and Tuttlingen, Germany, and Måløv, Denmark.

CA Available for sale in Canada



Ask Her If She Wants a Birth Control That's Hormone Free

Paragard® (intrauterine copper contraceptive) is the only highly effective, reversible birth control that is completely hormone free.¹ Paragard is a copper-releasing device (IUD) that is placed in the uterus to prevent pregnancy for up to 10 years. Paragard is FDA approved for over 30 years and has been used by millions of women.

- 100% hormone free
- 94% patient satisfaction with bleeding and cramping*2
- Removable whenever she decides§
- >99% effective for up to 10 years

The copper in Paragard enhances contraceptive efficacy by interfering with sperm transport and fertilization, and possibly preventing implantation.

Placement of Paragard

The Paragard non-hormonal IUD may be placed anytime during your patient's cycle when you and the patient are reasonably certain she is not pregnant, or after pregnancy. Please see important details, including contraindications.

Removal of Paragard

Because it is immediately reversible upon removal^{§1}, patients can start trying to conceive the same day Paragard is removed by a healthcare professional.

Ordering Paragard


Paragard Access Center offers a comprehensive range of support services and resources for ordering, including direct purchase, patient benefits verification service, fulfillment through our specialty pharmacy, and a solution for patients who choose to self-pay. More information can be found on the Paragard® Access Center webpage or by calling 1-877-Paragard.

Cost to Patients

Paragard may cost \$0 for insured patients under the Affordable Care Act. This may mean:

- No co-pay
- No coinsurance costs
- No deductibles

It is important to check with each patient's health insurance plan to determine coverage.




Indication

Paragard is intended for intrauterine contraception for up to 10 years.

Important Safety Information

- Paragard must not be used by women who have acute pelvic inflammatory disease (PID); have had a
 postpregnancy or postabortion uterine infection in the past 3 months; have cancer of the uterus or cervix;
 have an infection of the cervix; have an allergy to any component; or have Wilson's disease.
- If a woman misses her period, she must be promptly evaluated for pregnancy.
- Possible serious complications that have been associated with intrauterine contraceptives, are PID, embedment, perforation of the uterus, and expulsion.
- Paragard must not be used by women who are pregnant as this can be life threatening and may result in loss of pregnancy or fertility.
- The most common side effects of Paragard are bleeding and spotting; for most women, these typically subside after 2 to 3 months.
- Paragard does not protect against HIV or other sexually transmitted infections (STI).



^{*}Data are from the Contraceptive CHOICE Project. The study evaluated 3- and 6-month self-reported bleeding and cramping patterns in 5011 Long-Acting Reversible Contraceptive (LARC) users (n=826, Paragard), and the association of these symptoms with method satisfaction. Study participants rated satisfaction with their LARC method as "very satisfied," "somewhat satisfied," or "not satisfied." For the data analyses, "satisfied" and "very satisfied were grouped together as "satisfied."

[§] Paragard must be removed by a healthcare professional.

¹ Kaneshiro B, Aeby T. Long-term safety, efficacy, and patient acceptability of the intrauterine Copper T-380A contraceptive device. Int J Womens Health. 2010;2:211-220.

² Diedrich JT, Desai S, Zhao Q, Secura G, Madden T, Peipert JF. Association of short-term bleeding and cramping patterns with long-acting reversible contraceptive method satisfaction. Am J Obstet Gynecol. 2015;212(1):50.e1-50.e8.

Office Hysteroscopy



Introducing Endosee® —

A Truly Innovative Way to Perform Diagnostic Office Hysteroscopy in Any Room, at Any Site

Endosee Features Include:

- All-in-one, handheld, portable, cordless system
- Simple and quick to set up, needs only minimal staff training
- Lightweight reusable handset
- Single-use, flexible, thin (<5mm) sterile cannula
- New camera and light source every time so there's no degradation in visualization
- Bright, 3.5-inch touch-screen color LCD
- Still image and video capture
- Channel for fluid infusion
- Dock for battery charge and data transfer
- Rechargeable battery lasts more than 2 hours
- Procedure is coded and reimbursed using current hysteroscopy codes


Endosee® hysteroscopy simplifies: portability



setup space cost visualization everything

Blind Endometrial Biopsy vs. Hysteroscopy in Diagnosing Abnormalities

Sensitivity for Diagnosing Polyps¹ Blind Biopsy 11% Hysteroscopy



Sens	sitivity	for Dia	gnosing	Hyper	olasia
Blind	Biopsy				
		25%	/ 0		
Hyste	eroscopy			_	
			74%		



Endosee		CA*
Number	Description	
ES8000	Endosee Hysteroscope - Includes handset, docking station, AC adapter, USB cable	5
ESDX5	Endosee Dx Cannula - Disposable, single-use, sterile diagnostic cannula (Box of	5)

Optional Accessories		4
Number	Description	
ESSYR	50mL Syringe, sterile (40/Case)	
ESTUB	Standard Bore Extension Tubing 30 in (Sterile 100 per case)	

Urodynamics

Lumax TS Pro[™] System


Diagnose Urinary Incontinence Accurately and Quickly



CooperSurgical's Lumax TS Pro Fiberoptic Cystometry System is the leading multichannel office urodynamic system for diagnosing female urinary incontinence. This state-of-the-art urodynamic system is uniquely tailored to meet the needs of the busy Ob/Gyn practice. Its intuitive software and touch screen interface facilitate patient studies, while using the patented fiberoptic Lumax catheters ensures accurate results.

Lumax	TS Pro Package CA * C€
Number	Description
TSELITE	Lumax TS Pro Elite System includes: TS Pro Software, Laptop Computer, USB Cable, Uroflow Transducer, Volume Transducer, Remote Handheld Event Marker, Remote Footswitch Marker, EMG Module and Supplies for two Patient Studies
TSADV	Lumax TS Pro Advanced System includes: TS Pro Software, Laptop Computer, USB Cable, Uroflow Transducer, Volume Transducer, Remote Handheld Event Marker, Remote Footswitch Marker and Supplies for two Patient Studies
TSBAS	Lumax TS Pro Basic System includes: Uroflow Transducer, Volume Transducer, Remote Handheld Event Marker, Remote Footswitch Marker and Supplies for two Patient Studies
TSAUP	Upgrade from TSBAS to TSADV. TS Pro Software, Laptop Computer, USB Cable
EMG	Upgrade from TSADV to TSELITE. EMG Module, Mounting Hardware and Software





High-tech, Built-in Transducer Simplifies Pressure Sensing

The Lumax TS Pro System performs a complete differential diagnosis of lower urinary tract dysfunction including uroflowmetry, cystometry, urethral pressure profile, leak point pressures, pressure flow studies and electromyography. Lumax Fiberoptic Catheters simplify and speed up these procedures, providing highly accurate bladder and abdominal pressure measurements.

- 1. Light generated by the Lumax system travels through an optical fiber
- Diaphragm depresses when pressure is exerted, occluding a percentage of light
- 3. Gold reflector returns unblocked light
- 4. Returning light travels through an optical fiber
- Light received by the Lumax system is converted precisely into pressure (cm H₂O)

THE RESULT: Accurate data for appropriate treatment of your patients.

Lumax TS Pro Disposables		CA * (€		
Number	Description	Number	Description	
10307-010	10 Fr Bladder Catheter (10/Box)	53364	EMG Lead Cable — Reusable (1/Pack)	
10307-005	7 Fr Bladder Catheter (10/Box)	10310-000	S-Series Transmission Cable	
10307-009	10 Fr Abdominal Catheter (10/Box)	10206-000	Thermal Printer Paper (5 Rolls/Box)	
53397	EMG Electrodes (4 Sets/Pouch)	53604	Cystometry Tube Sets (10 Sets/Box)	

Accessorie	es		CA * (€
Number	Description	Number	Description
23486	Lumax TS Pro Software with Hasp Key and Documentation	10300-001	Deluxe bar code scanner
35821	Lumax TS Pro User Manual	60890-000	Bar code scanner brackets (upgraded scanner)
35384	Lumax User Manual	60898-000	Bar code scanner screws (upgraded scanner)
IVBT	IV pole with complete accessories: locking wheels, basket,	60819-000	Funnel
	laptop and monitor tray with catheter hooks	60897-000	Interface low pole mount (height 32" from floor)
IVB	IV pole with locking wheels and basket	10536-001	Uroflow Transducer
23467	Wheelbase only for IV pole (a replacement part only)	10537-001	Volume Transducer
23465	Basket with clamp to secure to IV pole	60996-000	IV pole power strip – deluxe model (4 outlets)
23463	Laptop tray with clamp to secure to IV pole	10440-000	Printer paper for monitor 4114UF (Pro) (5 Rolls/Box)
23462	Monitor pole clamp with mounting plate (for center mount of monitor instead of tray mount)	60997-000	Velcro® tie-down strap for power cord
35799	Monitor tray with catheter hooks for IV pole	60450-000	Clamp to attach laptop tray to IV pole
60598-000	Catheter supplies tray with IV pole mount	60688-006	Beaker
23466	USB cable 3'	60958-001	Transmission cable clip
53326	Remote footswitch marker	60958-002	Catheter bag clip
53298	Remote handheld marker	60450-000	Pole mount clip
53309	Uroflow Transducer	10388-4	Lumax Pro CD-ROM (replacement only)
53310	Volume Infused Transducer	23335	USB/serial port adapter
60993-000	Commode, wide with drop arms	10372-000	S-Series Installation Kit

Colposcopy

Examination of the Cervix Under Direct Magnification

NOW AVAILABLE THROUGH DISTRIBUTION



LEISEGANG celebrates more than 60 years as a world leader in colposcope design and production, recognized for German precision, clarity and durability. LEISEGANG OptiK LED colposcopes shed an uncompromising light and bright optical brilliance on all visual assessments. Select the OptiK Model 1 for routine colposcopy or the OptiK Model 2 for digital colposcopy. Both offer the highest-intensity LED on the market matched with the best optics for enhanced clarity and contrast. The optical brilliance achieved in LEISEGANG is unparalleled.

High-Intensity LED Illumination

LEISEGANG colposcopes use one of the brightest LED lights on the market. It transmits 2.5x more luminance than halogen for greater visibility of cellular and vascular changes, it outperforms other LED-lighted colposcopes in terms of stereoscopic 3-D visualization and matching daylight color temperature, and it delivers a stable color spectrum. This provides for highly improved contrast balance, better clarity and less shadowing for unsurpassed anatomy visualization. Plus, our LEDs have a working life of 20,000 hours or about 20 years of service...ensuring long-lasting dedicated service without compromise.

Precision German Optics

Specially selected, hand-polished and hand-set glass optics are used in all LEISEGANG OptiK colposcopes. The use of internal multiple mirrors is minimized to allow for maximum light exposure. Small apertures provide the greatest depth-of-field for unsurpassed clarity. Precision manufacturing ensures a parfocal magnification for dynamic range of detail. This means that once focused on a region of interest, the colposcope will stay in focus throughout the 3-step magnification. Furthermore, the optics provide for a convergent optical pathway, which mimics the human eye and reduces eye fatigue.

LEISEGANG OPTIK LED






100%

Dimmed

HALOGEN





100%

Dimmed

DAYLIGHT



OptiK® Model 1

Better Contrast, Clarity and Visualization



The Perfect Choice for Routine Colposcopy Exams

LEISEGANG's highly acclaimed colposcopes are renowned for unparalleled precision, uncompromising quality and lasting performance. Specially selected, hand-polished and hand-set German optics, parfocal 3-step magnification, and convergent optical beams contribute to LEISEGANG's remarkable depth and clarity. These light-emitting diode (LED) colposcopes offer the highest-intensity LED light on the market, coupled with improved balance contrast and increased stereoscopic 3-D clarity.

The OptiK Model 1 is available in four base configurations.



Four Different Choices:



Upright Rolling Base – Sturdy and mobile stand



Tilt Stand — Solid, durable and easy to maneuver



Swing Stand —
Standard swing arm for ease of use



Balance-O-Matic Stand —
Best combination of support and stability

LEISEGANG OptiK® Model 1		CA* (€
Number	Description	
OPTIK1-00	LEISEGANG OptiK Model 1 Tilt Colposcope (1DL LED): 3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, Tilt Base	
OPTIK1-01	Same as OPTIK1-00 except with Upright Rolling Base	
OPTIK1-02	Same as OPTIK1-00 except with medium-weight Swing-Arm Rolling Base	
OPTIK1-03	Same as OPTIK1-00 except with heavy-weight Balance-O-Matic Rolling Base	

OptiK® Model 2

Exceptional Optics and Optimized Digital Imagery

The 18-megapixel Canon® 600D T3i dSLR provides digital documentation of your colposcopic images. Its unique adjusting display allows for greater flexibility in your office.



From Routine Office Examinations to Advanced Clinical Applications, We Have a Solution to Meet Your Needs:

Infinite Performance

OptiK Model 2 with a Digital SLR Camera

For conducting research, delivering presentations or meeting detailed photo documentation requirements in your patient records, the OptiK Model 2 with a Digital SLR Camera precisely tracks lesion progression or regression.

Live Viewing

OptiK Model 2 with an Integrated Camera

The OptiK Model 2 with an Integrated Camera provides live image viewing and is an excellent choice for teaching. Designed for conducting research, delivering presentations or meeting detailed photo documentation requirements in your patient records.

Four Different Base Choices:

OptiK Model 2 is available in four base configurations including Upright, Tilt, Swing Stand and Balance-O-Matic.
The OptiK Model 2 features a working distance of 300 mm, field of view of 19 mm to 76 mm, and depth of field of 2 mm to 8 mm.



Upright Rolling Base – Sturdy and mobile stand



Tilt Stand — Solid, durable and easy to maneuver



Swing Stand –
Standard swing arm for ease of use



Balance-O-Matic Stand —
Best combination of support and stability

(See OptiK Model 2 Sales Configuration Options on Page 16)

OptiK® Model 2 Accessories

Cameras and Monitors		CA* (€
Number	Description	
53576	Canon® Rebel T6i Camera Kit for OptiK Model 2 1.0" CCD only	
23676FG	Foot Pedal for Canon Rebel T6i Camera	
23924	24" Digital High Definition LCD Monitor Display with HDMI input for use with Canon Digital Cameras	
22140	10X Oculars - Pair	
22114	20X Oculars - Pair	
40G301	Colposcopy Charts (50/Pad)	



LEISEGANG OptiK® Model 2 Sales Configuration Options:

LEISEGANG	OptiK Model 2 and Base Only	CA * (€
Number	Description	
OPTIK2-03	3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, 1.0" CCD optimized, upright rolling T2 adapter/tubus, for use with Canon® camera kit 53436	g stand,
OPTIK2-04	Same as OPTIK2-03 except with medium-weight Swing-Arm Rolling Base	
OPTIK2-05	Same as OPTIK2-03 except with heavy-weight Balance-O-Matic Rolling Base	

LEISEGANG	OptiK Model 2, Base and dSLR Camera	CA* (€
Number	Description	
OPTIK2-04-002	3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, 1.0" CCD optimized, medium-we Swing-Arm Rolling Base, T2 adapter/tubus, Canon 600D Rebel T3i camera kit with 18 MP camera, battery USB wire, A/C adapter, wired remote and memory card	
OPTIK2-05-002	Same as OPTIK2-04-002 except with heavy-weight Balance-O-Matic Rolling Base	

LEISEGANO	G OptiK Model 2 with Integrated Camera	CE
Number	Description	
OPTIK2-10	LEISEGANG OptiK Model 2 Integrated Analog Video Tilt Colposcope (3MVL LED Y/C): 3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, integrated S-Video analog camera (752 H \times 582 V), Tilt Base, for use with LCD Monitor 37464 $-$ 23598	
OPTIK2-11	Same as OPTIK2-10 except with Upright Rolling Base	
OPTIK2-12	Same as OPTIK2-10 except with medium-weight Swing-Arm Rolling Base	

LEISEGANG OptiK Model 2 with Integrated Camera Packages		
Number	Description	
OPTIK2-10-001	LEISEGANG OptiK Model 2 Integrated Analog Video Tilt Colposcope (3MVL LED Y/C) with LCD Monitor: 3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, Integrated S-Video analog camero (752 H x 582 V), Tilt Base, LCD Monitor	3
OPTIK2-11-001	Same as OPTIK2-10-001 except with Upright Rolling Base	
OPTIK2-12-001	Same as OPTIK2-10-001 except with Swing-Arm Rolling Base	



German Engineered Euro-Med® Instruments

Handcrafted to exacting standards by the world's preeminent steel artisans, Euro-Med stainless steel surgical instruments represent a higher level of engineering, control and precision. Essential for the success of women's health procedures, these instruments are ideal for doctors and practices who use only the best to ensure exceptional results.

 Classic Series Tischler-Morgan™ Biopsy Punch



Euro-Med[®] Stainless Steel Biopsy Punches

Beautifully Balanced for Superior Handling, Response and Comfort



Our most popular design includes ergonomic contours for maximum control and comfort. The standard bite incorporates a single tooth on the upper and lower jaws that draws tissue in for maximum sample retention with minimal trauma.

Number	Length	Bite
64-689	21 cm	7 mm x 3 mm
64-689L	25 cm	7 mm x 3 mm




This classic large rectangular bite includes upper and lower teeth for added traction and maximum sampling.

Number	Length	Bite
64-679	22 cm	9 mm x 3 mm



The petite bite minimizes trauma and is perfect for obtaining samples of lesions just inside the lip of the endocervical canal.

Number	Length	Bite	
64-690	21 cm	4 mm x 2.3 mm	



Our popular Euro-Med handle grants optimal comfort, while the small bite incorporates a tooth permitting sampling at the endocervical canal. This variation of the Mini-Townsend provides samples with minimal patient discomfort and bleeding.

Number	Length	Bite
64-691	21 cm	4.1 mm x 2.3 mm
64-691L	26 cm	4.1 mm x 2.3 mm



A slightly smaller bite and additional anchoring teeth on the bottom jaw ensure ample, well-oriented samples with minimal trauma.

Number	Length	Bite
64-675	22 cm	9 mm x 3 mm



The unique jaw shape yields a distinct sample orientation for histologic analysis.

Number	Length	Bite
64-660	22 cm	7.5 mm x 3 mm

Burke™ Biopsy Punch

CA* (€



The design features an open handle for even the largest of hands and an oval-shaped bite for small samples. Sharp tooth design helps to anchor the tissue for sampling and retention.

Number	Length	Bite	
64-677	22 cm	5.5 mm x 3 mm	

Baggish Biopsy Punch

CA* (€



Another variation of the open-handle design, this punch provides a larger, oval-shaped bite. The jaws have two bottom teeth and a serrated bed to provide optimal retention of tissue samples.


Number	Length	Bite	
64-695	22 cm	6.5 mm x 3.5 mm	

Euro-Med Keyes Vulvar Punch Set and Singles CA★ CE




This essential set includes six stainless steel interchangeable tips ranging from 3 mm to 8 mm in diameter, and a universal handle, all set in a brass tray fitted in a wooden box for simple storage. The individual tips can be placed on the threaded handle as required. The brass insert facilitates set sterilization. Also available: individual "singles" permanently attached to their own handle.

Number	Length
68-180	Set includes: five tips, storage tray and wooden box
68-173	3 mm – Singles
68-174	4 mm – Singles
68-175	5 mm — Singles
68-176	6 mm – Singles
68-178	8 mm – Singles



These rotating punch tips fit either the Euro-Med® closed rotating handle or the Townsend open rotating handle. With eight tips to select from, complete customization is possible.



Euro-Med® Standard Biopsy Punches





CA* (€

One of the classic original biopsy designs, this oval cutting jaw obtains an 8 mm x 4 mm cervical specimen. A special basket in the jaw helps retain the biopsy sample.



Number	Length	Bite
64-687	21 cm	8 mm x 4.5 mm

CA₩ (€

The unique angle of the jaw enables extra access to hard-to-reach areas of the cervix. The large teeth and strong bite provide a maximum tissue sample area.



Number	Length	Bite	
64-681	21.5 cm	8 mm x 4 mm	

CA* (€

Offering the popular Kevorkian bite, a strong and narrow jaw and an added 7.5 cm shaft length, this punch facilitates simple insertion and sampling of the cervical canal, especially in patients for whom the extra length is necessary.



Number	Length	Bite
64-685	26.5 cm	9 mm x 3 mm

CA* (€

One of the more popular biopsy punches used in Europe, the Schubert features a large bite and an angled handle to apply maximum pressure.




Number	Length	Bite	
64-670	21.5 cm	8 mm x 4.5 mm	

Engineered to Exacting Standards Stainless Steel Specula for

Exceptional Visibility and Access

Fundamental? Yes. Basic? Far from it. Euro-Med® provides the essentials for any practice, like our complete array of premium-grade stainless steel specula. Each instrument is designed to ensure your patient's comfort and meet your practice's needs. Ask us about our specially coated stainless steel specula and other instruments used during LEEP procedures (see pages 48-50).





Specula Blades and Yoke Dimensions

Use our convenient speculum sizing guide to help you easily select the proper tools for your practice. Specula are available in a variety of sizes to accommodate all your patients' needs.

In addition, we are now offering a range of extra-large to 4X size specula to facilitate examinations and procedures for larger or obese patients.

See pages 24-25 for a full product selection and model numbers.

Also Available with

Side Opening Yoke

Side Opening Yoke

Vu-Max[™] Yoke – Graves Blades

Medium – W: 33 mm L: 114 mm Large – W: 35 mm L: 125 mm

Vu-Max Yoke – Pederson Blades

Medium – W: 22 mm L: 114 mm Large – W: 25 mm L: 120 mm

Vu-More[™] Yoke – Graves Blades

Mini – W: 34 mm L: 93 mm Medium – W: 33 mm L: 113 mm Large – W: 35 mm L: 127 mm

Standard Yoke – Pederson Blades

Virgin – W: 13 mm L: 78 mm Small – W: 13 mm L: 82 mm Extra Narrow – W: 17 mm L: 108 mm Medium – W: 22 mm L: 114 mm Large – W: 25 mm L: 130 mm

Standard Yoke — Graves Blades
Also Available with Small — W: 19 mm 1: 82 mm

Small – W: 19 mm L: 82 mm Medium – W: 33 mm L: 114 mm Large – W: 35 mm L: 127 mm

Vu-More Yoke – Graves or Pederson Blade

Snowman[™] – W: 35 mm L: 152 mm Snowman – W: 25 mm L: 152 mm

Vu-Max and Vu-More Yokes – Graves Blades

X Large – W: 35 mm L: 150 mm 2X Large – W: 35 mm L: 160 mm 3X Large – W: 35 mm L: 170 mm 4X Large – W: 35 mm L: 180 mm

Euro-Med® Stainless Steel Specula



CA* (€

The popular Euro-Med Graves vaginal speculum features a slightly wider traditional blade design. Available with standard yoke and an open-sided configuration. The Graves blade is also available in longer lengths with a Vu-More[™] or Vu-Max[™] yoke.

Number	Description	Size mm
64-100	Small	W19 L82
64-102	Medium	W33 L114
64-104	Large	W35 L127
64-108	Open-sided	W35 L100

Pederson

CA***** (€

One of the most widely used speculum designs, this model, with its slightly thinner traditional blade, is particularly useful with young or nulliparous women. The Pederson blade is also available with a Vu-Max yoke.

Number	Description	Size mm
391-671	Virgin	W13 L78
64-110	Small	W13 L82
391-674	Extra-narrow	W17 L108
64-112	Medium	W22 L114
64-114	Large	W25 L130

Vu-More[™]

CA***** (€

The Vu-More provides a 25 percent larger viewing area at the yoke than a comparable Graves to afford easier instrumentation access and improved visibility. The Mini's short blade facilitates examinations for women with cervical prolapse.

Number	Description	Size mm
64-130	Mini	W34 L93
63-902	Medium	W33 L113
63-900	Large	W35 L127
63-905	X Large	W33 L150 ₹
63-906	2X Large	W33 L160
63-907	3X Large	W33 L150
63-908	4X Large	W33 L180

(W=Distal End) (L=Penetration Length)







Vu-Max[™]

The Vu-Max vaginal speculum's unique design substantially reduces the chance of obstructing the physician's view or instrument access. Offered in both Pederson and Graves configurations.

in both reactson and draves configurations.			
Number	Description	Size mm	
64-113	Pederson Medium	W22 L114	
64-115	Pederson Large	W25 L120	
64-218	Graves Medium	W33 L114	
64-219	Graves Large	W35 L125	
64-225	X Large	W35 L150	201
64-226	2X Large	W35 L150 W35 L160	-

W35 L170

W35 L180

CA# (€

Snowman specula are built using a 40 mm yoke opening for optimized strength and visualization. Designed to naturally fit the obese patient, their longer blades safely expose and access the cervix.

Number	Description	Size mm
64-107	Graves Blade	W35 L152
64-117	Pederson Blade	W25 L152
64-127	One each/Graves	& Pederson

ExpandaView™ Integrated CA CA **Speculum and Lateral Retractor**

The ExpandaView is designed to enhance clinician visualization. Combines a traditional speculum with a lateral wall retractor to enhance the view of the cervix by retracting redundant vaginal tissue. Especially effective for large or multiparous patients.

Number	Description	Size mm
64-161	Medium	W22 L109
64-162	Large	W22 L133
64-163	Stainless Steel Instrument Holder	



3X Large

4X Large

64-227

64-228



Weisman Graves CA# (€ **Side Opening**

Weisman Graves specula features an open side for improved instrument access and easy removal once instruments are placed. Available in two sizes with left or right openings.

Number	Description	Size mm
391-667A	Left side open, Med	W33 L102
391-667	Left side open, Lg	W36 L121
391-668A	Right side open, Med	W35 L102
391-668	Right side open, Lg	W36 L120



Pederson Side Opening CA* (€

The Pederson open-sided specula features the traditional thinner blade design while providing the advantages of either a left or right side opening for efficient instrument access and removal.

Number	Description	Size mm
64-112L	Left side open, Med	W22 L114
64-114L	Left side open, Lg	W25 L120
64-112R	Right side open, Med	W22 L114
64-114R	Right side open, Lg	W25 L120



CA* (€

This popular 2.5-pound weighted speculum is the standard for posterior vaginal wall retraction in vaginal procedures.

Number	Description	Size mm
64-422	Med (Wide/Short Blade)	W43 L55
64-421	Lg (Narrow/Long Blade)	W35 L90

Stainless Steel Lateral Wall Retractors



Cer-View [™]	CACCA
Lateral Wall Retractor	

Contoured and shaped to fit into the yoke of a standard speculum, the Cer-View retractor expands beyond the width of the speculum blades and retracts the side walls of the vaginal vault to achieve an unobstructed view of the cervical anatomy. Cer-View is especially useful in the presence of vaginal prolapse.

Number	Description	Size mm	
64-310	Medium	W22 L150	
64-320	Large	W22 L160	

Tru-View[™] CA***** C€ Lateral Wall Retractor

Designed specifically for use with obese patients, the Tru-View has longer tapered blades that fit a long vaginal canal without pinching in the yoke area. The unique arm configuration provides added strength and support under prolapse conditions. Tru-View is a natural fit when combined with the Snowman™ speculum.

Number	Description	Size mm
64-307	Extra Large	W32 L190

Jackson CA♦ C€ Vaginal Retractor

This traditional retractor can be inserted and pulled downward to provide maximum posterior viewing. Like all Euro-Med specula and retractors, the Jackson features a polished stainless steel blade for maximum patient comfort.

Number	Description	Size mm
64-322	Medium	W40 L98
64-324	Large	W40 L115

Sounds and Dilators

Sims Uterine Sound CA★ C€



Sims Uterine Sounds are fully graduated for accurate fundal depth determination and feature a bulbous tip to reduce the risk of perforation. Part number 64-600 is made of silver-plated copper, which makes it ideal for situations where a more malleable sound is useful. Part number 64-601 is not as malleable but will resist corrosion due to its stainless steel properties. Both can adapt to virtually any anatomical variation.

Number	Description	Size
64-601	Sims Uterine Sound – Stainless Steel	L12.5"/317 mm
64-600	Sims Uterine Sound — Silver-Plated Copper*	L12.5"/317 mm

^{*}Silver-plated copper may tarnish.

Dil-Os[™] Uterine Sound/Dilator







the internal Os, give accurate determination of uterine depths and provide the tactile determination of intrauterine adhesions and septa.

Number Description Size

Number	Description	Size
64-336	Slight Curve	L12.7"/322 mm
64-337	Strong Curve	L12.5"/317 mm
64-338	Set (64-336, 64-337)	

Stainless Steel Sounds and Dilators



Mini Cervical Dilator Set CA € C€

Ideally suited for nulliparous patients or those with a stenotic Os. Includes a universal handle that holds dilator tips from sizes 1 to 3 mm in half-millimeter increments. Rounded dilator tips help prevent uterine wall puncture.

Number	Description
64-374	5 Dilators, Handle and Storage Pouch



Small Cervical Dilator Set CA € C€

This small cervical dilator set is designed for patients with a stenotic cervical Os. Set of six dilators including sizes 1 through 6 mm. Comes in a fitted and numbered storage pouch.

Number	Description
64-372	6 Dilators and Storage Pouch



Hank Cervical Dilators CA ◆ C€

Hank cervical dilators come in a set of six or individually, ranging from 9 French (3 mm) to 20 French (6.66 mm). The set comes in a custom-fitted and numbered pouch for easy storage and selection.

Number	Size mm
64-370	6 Dilators and Storage Pouch
61201	9 Fr (3 mm)/10 Fr (3.33 mm)
61202	11 Fr (3.67 mm)/12 Fr (4 mm)
61203	13 Fr (4.33 mm)/14 Fr (4.67 mm)
61204	15 Fr (5 mm)/16 Fr (5.33 mm)
61205	17 Fr (5.67 mm)/18 Fr (6 mm)
61206	19 Fr (6.33 mm)/20 Fr (6.67 mm)

Hegar Dilator Set

CA₩ (€

Stainless steel Hegar dilators provide the option of double-ended dilators in even millimeter increments from 3 to 18 for a greater diameter range over the entire set. Each double-ended dilator is custom-labeled accordingly. Sold in a set of eight or individually for added convenience.

Number	Size mm
61207	9 Fr (3 mm)/12 Fr (4 mm) 51 Fr (17 mm)/ 54 Fr (18 mm) (8/Set)
61208	9 Fr (3 mm)/12 Fr (4 mm)
61209	15 Fr (5 mm)/18 Fr (6 mm)
61210	21 Fr (7 mm)/24 Fr (8 mm)
61211	27 Fr (9 mm)/30 Fr (10 mm)
61212	33 Fr (11 mm)/36 Fr (12 mm)
61213	39 Fr (13 mm)/42 Fr (14 mm)
61214	45 Fr (15 mm)/28 Fr (16 mm)
61215	51 Fr (17 mm)/54 Fr (18 mm)

Pratt Dilator Set

CA* (€

Stainless steel Pratt dilators offer the largest diameters in our line of cervical dilators, ranging from 13 French (4.3 mm) to 59 French (19.6 mm). Each is double-ended and custom-labeled accordingly in French scale measurements. Sold in a set of eight or individually for added convenience.

391-830
391-834 17 Fr / 19 Fr 391-835 21 Fr / 23 Fr 391-836 25 Fr / 27 Fr 391-837 29 Fr / 31 Fr 391-838 33 Fr / 35 Fr 391-839 37 Fr / 39 Fr 391-840 41 Fr / 43 Fr
391-835 21 Fr / 23 Fr 391-836 25 Fr / 27 Fr 391-837 29 Fr / 31 Fr 391-838 33 Fr / 35 Fr 391-839 37 Fr / 39 Fr 391-840 41 Fr / 43 Fr
391-836 25 Fr / 27 Fr 391-837 29 Fr / 31 Fr 391-838 33 Fr / 35 Fr 391-839 37 Fr / 39 Fr 391-840 41 Fr / 43 Fr
391-837 29 Fr / 31 Fr 391-838 33 Fr / 35 Fr 391-839 37 Fr / 39 Fr 391-840 41 Fr / 43 Fr
391-838 33 Fr / 35 Fr 391-839 37 Fr / 39 Fr 391-840 41 Fr / 43 Fr
391-840 41 Fr / 43 Fr
391-840 41 Fr / 43 Fr
391-841 45 Fr / 47 Fr
391-842 49 Fr / 51 Fr
391-843 53 Fr / 55 Fr
391-844 57 Fr / 59 Fr

Stainless Steel Endocervical Curettes



A narrow tip permits sampling without dilation, while the stainless steel blade resists clogging and dulling. Includes basket for sample retention.

 Number
 Length

 66-442
 12"/305 mm

Like the Kevorkian with basket, this curette's design is traditionally weighted for optimal balance and control.

Number	Length	
66-440	12"/305 mm	

The rounded and tapered tip facilitates entry through a stenotic Os. The design of the cutting edge yields improved endocervical sampling.

Number	Length
66-445	12"/305 mm



Another Euro-Med exclusive with a unique basket that offers greater assurance of sample retention. Rounded tip facilitates cervical entry.

Number	Length	
66-441	12"/305 mm	



CA₩ CE

The flanged basket facilitates sample collection and ensures optimal sample retention while exiting. The design also incorporates a rounded, tapered tip to permit nontraumatic entry into a stenotic Os.

Number	Length	
64-490	11.7"/297 mm	



Our revolutionary design allows optimal sample collection with a simple twirl motion — no pushing or pulling. Closed blade on both ends eliminates the chance of sample loss and reduces patient discomfort.

Number	Length
66-447	12.3"/312 mm

Endometrial Curettes



Our most popular reusable endometrial sampler. Features a suction attachment for improved sampling.

Number	Size	Length
64-624	1 mm	9.8"/248 mm
64-626	2 mm	9.8"/248 mm
64-628	3 mm	9.8"/248 mm
64-630	4 mm	9.8"/248 mm

The less aggressive Randall endometrial sampler also has a suction attachment for improved sample collection.

Number	Size	Length
64-620	3.9 mm	9.9"/252 mm



The premier choice for uterine curettage. Sharp and malleable, they come in many different sizes and are an industry standard 11 inches long.

Number	Size	Length
392-110	6 mm	11"/280 mm
392-111	7 mm	11"/280 mm
392-112	8 mm	11"/280 mm
392-113	9 mm	11"/280 mm
392-114	11.5 mm	11"/280 mm
392-116	17 mm	11"/280 mm

Stainless Steel Cervical Tenacula




This tenaculum is an industry standard offering maximum cervical stability via its two opposing hooks. The long length provides excellent cervical access.

Number	Length
64-420	9.9"/252 mm



This tenaculum has pincers like the Braun tenaculum but provides additional length and a double curve for more challenging situations.

Number	Length
61352	11"/280 mm



The Schroeder tenaculum is essentially the same as the Braun except that it has more rounded pincers.

Number	Length	
64-410	9.9"/252 mm	



The Goldstein Grasp provides rapid atraumatic stabilization of the cervical Os during placement of catheters, endometrial samplers and other intrauterine devices, with minimal tissue trauma.

Number	Length
64-6026	9.4"/239 mm

Stainless Steel Hooks

Duplay Tenaculum Double Curved



These convenient hooks offer a unique angle to allow better access to and visualization of the cervix during use. Available configurations: straight, left and right.

Number	Description	Length
36-630	Straight Hook	9.8"/249 mm
36-631	Left Hook	9.6"/244 mm
36-632	Right Hook	9.6"/244 mm

Iris Hook CA ♥ C€

CA# (€

This petite and delicate hook readily stabilizes the cervix during procedures and allows greater access.

Number	Length
36-625	10.3"/262 mm

Marlow Corkscrew Hook CA◆ C€



Excellent for stabilizing the cervix during invasive procedures. This unique hook can be screwed into the cervical tissue that needs to be removed so the sample will not be lost.

Number	Length
67-451	9.3"/236 mm

Double IUD Hooks CA★ C€



Dual grasping hooks for better IUD retrieval. With a working length of nearly 27 cm and an angled design, the Double IUD Hook has a rounded tip for easy cervical insertion.

Number	Length
392-151	10.6"/269 mm

Emmett Tenaculum Hook CA ◆ C€

Larger hook affords added stability and control of the cervix during procedures.

Number	Length
66-450	9.8"/249 mm



For easy and effective retrieval and removal of an IUD. Made of malleable stainless steel.

Number	Length	
392-148	10.3"/262 mm	

Stainless Steel Endocervical Specula



Ratchet Handle Speculum

CA* (€

By far our most popular handle design. The easy single-handed operation with quick-release lock delivers the versatility and performance needed for the most delicate of procedures.

Number	Tip Size	Length
66-437	Small 3 mm	10.6"/270 mm
66-435	Medium 4 mm	10.6"/270 mm
66-436	Large 5 mm	10.6"/270 mm

Screw / Knob Handle Speculum CA◆ C€

The screw/knob handle speculum easily adjusts and maintains position using the threaded thumb screw.

Number	Tip Size	Length
66-431	Small 3 mm	10.6"/270 mm
66-430	Medium 4 mm	10.6"/270 mm
66-432	Large 5 mm	10.6"/270 mm

No-Lock Handle

CA₩ (€

This innovative stainless steel endocervical speculum offers precise singlehanded action. For those who prefer total control, the No-Lock Handle allows immediate adjustment as necessary.

Number	Tip Size	Length
66-433	Small 3 mm	10.6"/270 mm
67-435	Medium 4 mm	10.6"/270 mm
66-434	Large 5 mm	10.6"/270 mm

Each Endocervical Speculum Can be Ordered with a Small, Medium or Large Tip







Stainless Steel Sponge Forceps

Foerster Sponge Forceps

Foerster sponge forceps are essential for gynecologic surgery. Featuring oval tips with serrations, the Foerster is used to hold sponges during the procedure. The ratchet handle helps to hold the sponges firmly.

Number	Length
16-182	9.7"/246 mm



Stainless Steel Hemostatic Forceps

Halsted Mosquito (Straight or Curved)

CA# (€



With their fully serrated jaw and 5-inch working length, the narrow-pointed Halsted Mosquito forceps are typically used to clamp off smaller vessels that control the flow of blood. Incorporate ratchet closure for added control.

Number	Description	Length
61902	Straight	5"
61903	Curved	5"

Adson CA**♦** (€



Longer than the Halsted Mosquito, the 7¼-inch Adson's narrow-pointed jaw is typically used to clamp off smaller vessels. Ratchet handle provides added control.

Number	Description	Length
61931	Straight	7.25"



CA* (€



Slightly longer than the Halsted Mosquito, the Kelly forceps have a longer jaw to clamp off vessels that control blood flow. The jaw is half-serrated and can also be used for blunt dissection. Incorporate ratchet closure for added control.

Number	Description	Length
61908	Straight	5.5"
61909	Curved	5.5"

Rochester-Pean (Straight or Curved)

CA* (€



Sometimes referred to as the "Big Hemostat," the Rochester-Pean (available in sizes up to 10¼ inches) has both straight and curved styles, with fully serrated jaws. Typically used to occlude or clamp off larger vessels to control bleeding. The built-in ratchet mechanism holds firm so you don't have to.

Straight

Number	Length
61939	5.5"
61940	6.25"
61941	7.25"
61942	8"
61943	9"
61944	10.25"

Curved

Number	Length
61945	5.5"
61946	6.25"
61947	7.25"
61948	8"
61949	9"
61950	10.25"

Stainless Steel Dressing and Tissue Forceps

Allis

Available in three lengths, the Allis forceps feature an interlocking 4×5 or 5×6 tooth pattern on their distal tip for secure tissue retraction. The built-in ratchet mechanism will hold firm.

Number	Description	Length	
61981	4 x 5 Teeth	6"	
61983	5 x 6 Teeth	7.5"	
61882	5 x 6 Teeth	10"	



CA# (€

Potts-Smith Tissue CA ◆ C€



The Potts-Smith Tissue forceps are the tissue-grasping sibling to the Potts-Smith Dressing forceps. The distal end features a 1 x 2 tooth pattern for excellent tissue manipulation.

Number	Description	Length
61871	1 x 2 Teeth	7"
61872	1 x 2 Teeth	8.2"



Dressing forceps are the instrument of choice to help manipulate and pack dressings and sponges. These versatile forceps can also be used to handle and manipulate tissue.

Number	Length	
61804	5.6"	
61805	6.2"	



Longer and more delicate than the standard Dressing forceps, the Potts-Smith are used to pick up and manipulate dressings and sponges.

Number	Length
10-202	8.25"



Standard tissue forceps with a 1 x 2 tooth pattern to hold and manipulate tissue during suturing.

Number	Description	Length
61841	1 x 2 Teeth	5.75"
61843	1 x 2 Teeth	7"
61844	1 x 2 Teeth	8"

Adson Dressing Serrated CA★ C€



The serrated tips help to hold and manipulate tissue during surgery.

Number	Description	Length
61830	Serrated Tips	4.75"

Adson Tissue 1 x 2 Serrated CA* CE



The serrated tips help to hold and manipulate tissue during surgery.

Number	Description	Length
61864	1 x 2 Serrated Tips	4.75"



To more aggressively hold and manipulate tissue during surgery.

Number	Description	Length
61863	1 x 2 Teeth	4.75"

Adson-Brown	Grasping	CA# (€
-------------	----------	---------------



Features two rows of 7×7 interlocking teeth for less traumatic manipulation of tissue during surgery.

Number	Description	Length	
61866	7 x 7 Teeth	4.75"	

Balance and control are hallmark qualities of all Euro-Med instruments — whether it's a simple specula, a precision-built biopsy punch or a more delicate tool such as a dressing or tissue forceps. Trust Euro-Med instruments to help your procedures go smoothly and successfully.

Stainless Steel Uterine Forceps



Double-curved for optimal uterine insertion, the Bozeman Uterine Forceps are a universal tool for various situations such as holding gauze and sponges, manipulating tissue as well as doubling up as a longer hemostat.

Number	Length
392-660	10.2"/257 mm



Used to remove fibroids from the uterus.

Number	Length
417-131	12.5"/318 mm

Foerster Uterine Polyp Forceps





With a smaller distal tip than the Kelly Uterine Forceps and a Euro-Med locking ratchet mechanism, these forceps are used to remove polyps from the uterus.

Number	Description	Length
415-140	Straight	10"/254 mm
415-141	Curved	10"/254 mm

Stainless Steel Needle Holders



Mayo-Hegar (Wider Jaw) **CA*** (€

This essential Gyn instrument comes in various lengths from 5 to 8 inches. The serrated jaws hold suture needles true and tight for efficient tissue repair.

Number	Length
62809	5.5"
62810	6.3"
62815	7.1"
62820	8.1"

Heaney (Curved)

CA* (€

Known for its curved distal tip and 81/4-inch length, the Heaney needle holder features tungsten carbide inserts for outstanding holding power.

Number	Length	
62825	8.3"	

Crile-Wood (Narrow Jaw)

CA* (€

With a thinner distal tip than the Mayo-Hegar, the Crile-Wood is ideal for more delicate suturing applications.

Number	Length
62840	7.1"

Hand-Finished Stainless Steel Scissors

Operating Scissors (Straight)

The workhorse of surgical scissors. Available in three tip configurations: sharp/sharp (S/S), sharp/blunt (S/B) and blunt/blunt (B/B). Ideal for routine tissue cutting and dissection.

Number	Description	Length
62002	B/B	5.5"
62008	S/B	5.5"
62014	S/S	5.5"
62004	B/B	6.5"
62005	B/B	7.5"
62010	S/B	6.5"
62011	S/B	7.23"



Metzenbaum Scissors (Straight or Curved)

CA* (€



Essential scissors for cutting delicate tissue and blunt dissection. Having a longer handle-to-blade distance, the Metzenbaum is available in a variety of lengths; the 7-inch length is most preferred.

Number	Length
62116	7"
62127	8"
62120	9.1"
62122	11.2"

Spencer Suture Scissors

Curveo

Number	Length
62131	5.9"
62118	7"
62119	8.1"
62124	9"
62126	11"

02127	<i>-</i>
62126	11"
	CA# CE



Standard scissors for the removal of sutures. One blade is notched to slide underneath the suture to facilitate removal.

Number	Description	Length	
62164	Straight	3.5"	



CA# (€



A classic pair of scissors for any gynecologist office. While typically used for cutting more dense fascia, they can also be used to cut sutures. Available in different lengths; the 6%-inch is most preferred.

Straight

Number	Length
62063	5.7"
62065	6.8"
62066	9"

Number	Length
62067	5.5"
62069	6.75"
62070	9"

Iris Scissors (Straight or Curved)

CA* (€



"Must have" scissors for very delicate tissue dissection and the cutting of fine sutures. These petite scissors are available in the classic length of $4\frac{1}{2}$ inches.

Number	Description	Length
62140	Straight S/S	4.5"
62143	Curved S/S	4.5"

Endometrial Samplers

Pipelle® Endometrial Suction Curette

CA# CE



Pipelle's intuitive design provides consistency of sample collection every time. The procedure takes less than a minute to perform and produces little if any post-procedure bleeding or patient discomfort. Completely self-contained, the Pipelle requires no additional equipment to complete the procedure. In most cases the Pipelle can be inserted without the need for dilation. Flexibility of the curette has been optimized to facilitate adaptation to normal uterine curvature, thus promoting contact with the wall for safe, improved sample collection.

Number	Description
8200	Pipelle Endometrial Suction Curette (25/Box)

Combo-Pak Endometrial Sampling System



Ideal for patients requiring dilation before insertion due to a stenotic Os or convoluted endocervical canal. The Combo-Pak includes Pipelle endometrial samplers for optimal retrieval and sterile Os Finders for quick and easy dilation.

Number	Description
CSICP	Combo Pak – 25 Pipelle and 6 sterile Os Finders per box

Pipet Curet™





Doctor-preferred endometrial sampler for smoother sampling and patient comfort. Thicker stylet with double "O-Ring" maximizes suction. Centimeter markings are highly pronounced for better visibility. Larger opening facilitates maximum tissue retrieval for more reliable results. Available in three options. The Standard Pipet comes in 3 mm and 4 mm and the soft yellow 3 mm is pliable, allowing the catheter to bend for easier insertion.

Number	Description
MX140	3 mm (50/Box) — Greater integrity for stenotic Os
MX145	3 mm Soft Yellow (50/Box) — Pliable for molding catheter
MX150	4 mm (50/Box) – Multiparous patient

SelectCells Endometrial Samplers

CA*



SelectCells Endometrial Samplers offer proven clinical performance and reliability. The contoured finger grip, easy-to-read graduation marks and piston suction allow for smooth and efficient sampling. Available in Standard or Mini sizes, both samplers feature a narrow profile with a sharp edge for effective tissue collection.

SelectCells Mini offers the smallest tip diameter on the market. This powerful yet small endometrial sampler outperforms others, with its nontraumatic, easy-to-insert design. SelectCells Mini is the ideal choice for patients with stenotic Os or low pain tolerance.

Number	Description	
SM200	SelectCells Standard® Endometrial Sampler (25/Box)	
SM203	SelectCells Mini® Endometrial Sampler (25/Box)	

Uterine Explora™ Model I CA* C€ with Vacu-Lok Syringe



Utilizing a syringe suction technique for exceptional control and sample recovery, the cutting edge ensures optimal sample collection from all areas of the endometrial cavity. This disposable Explora Model incorporates a long nylon curette with Randall-type cutting edge, stylet and Vacu-Lok syringe that adjusts suction from approximately 15- to 20-inch Hg negative pressure (6 cc to 12 cc).

Number	Description
MX120	Uterine Explora with Vacu-Lok Syringe, 3 mm, Model I (25/Box)
MX124	Uterine Explora with Vacu-Lok Syringe, 4 mm, Model I (25/Box)

Uterine Explora™ Model II Dual-Port with Vacu-Lok Syringe

CA* €€



Uterine Explora Model II offers the same features as the Model I plus it has a dual-port opening that reduces suction force and discomfort for patients with low pain threshold or who are very apprehensive. Tissue is obtained via a scraping and peeling action. The Vacu-Lok syringe allows you to adjust the suction from approximately 15- to 20-inch Hg negative pressure (6 cc to 12 cc).

Number	Description
MX130	Uterine Explora II with Vacu-Lok Syringe Model II (25/Box)

Cannula Curette

The Cannula Curette is a sterile, latex-free disposable device consisting of a Karman-type suction curette for use with the 60 cc Handyvak Locking Syringe. This quiet, portable device obtains a histologic sample from the endometrium on an outpatient basis, usually without anesthesia. The Cannula Curette comes in a variety of sizes — smaller sizes are especially useful in the early detection of endometrial carcinoma or its precursors.

STAINLESS		CA*

With Syringe

Number	Description
MX622	2 mm (10/Box)
MX623	3 mm (10/Box)

SEMI-RIGID	CA₩
With Syringe	Without Syringe

With Syringe		Without Syringe	
Number	Description	Number	Description
MX603	3 mm (10/Box)	MX624	4 mm (10/Box)
MX604	4 mm (10/Box)	MX625	5 mm (10/Box)
MX605	5 mm (10/Box)	MX626	6 mm (10/Box)
MX606	6 mm (10/Box)	MX627	7 mm (10/Box)
MX607	7 mm (10/Box)	MX633	3 mm (10/Box)

FLEXIBLE	CA *
With Syringe	Without Syringe

Number	Description	Number	Description
MX503	3 mm (10/Box)	MX513	3 mm (10/Box)
MX504	4 mm (10/Box)	MX514	4 mm (10/Box)
MX505	5 mm (10/Box)	MX515	5 mm (10/Box)
MX506	6 mm (10/Box)	MX516	6 mm (10/Box)
MX507	7 mm (10/Box)	MX517	7 mm (10/Box)



RIGID: CUR\	/ED SHORT (NO SYRING	E)	CA*
Number	Description	Number	Description
MX588	8 mm (10/Box)	MX591	11 mm (10/Box)
MX589	9 mm (10/Box)	MX592	12 mm (10/Box)
MX590	10 mm (10/Box)		

RIGID PLASTIC: CURVED LONG (NO SYRINGE)			
Number	Description	Number	Description
MX508	8 mm (10/Box)	MX511	11 mm (10/Box)
MX509	9 mm (10/Box)	MX512	12 mm (10/Box)
MX510	10 mm (10/Box)		


RIGID PLAST	FIC: STRAIGHT LONG (NO	O SYRINGE)	CA*
Number	Description	Number	Description
MX558	8 mm (10 Box)	MX561	11 mm (10/Box)
MX559	9 mm (10/Box)	MX562	12 mm (10/Box)
MX560	10 mm (10/Box)		

Endometrial Samplers

Tis-U-Trap™ with Cannula Curette

Designed to work with a suction source such as a pump or wall suction unit, endometrial tissue is collected directly into the Tis-U-Trap collection chamber, eliminating the need to transfer the tissue sample to another container. The collection chamber's unique design facilitates sample size determination. (See page 28 for Reusable Stainless Steel versions.)

The smaller sizes are especially useful in early detection of endometrial carcinoma or its precursors and endometrial dating.



Tis-U-Trap: Flat Shape

WITH CANNULA CURETTE – FLEXIBLE		
Number	Description	
MX403	3 mm (10/Box)	
MX404	4 mm (10/Box)	

WITH CANNULA CURETTE — SEMI-RIGID	
Number	Description
MX413	3 mm (10/Box)
MX414	4 mm (10/Box)

WITH CANNU	WITH CANNULA CURETTE — STAINLESS	
Number	Description	
MX422	2 mm (10/Box)	
MX423	3 mm (10/Box)	

Tis-U-Trap: Cone Shape

WITH CANN	WITH CANNULA CURETTE – FLEXIBLE	
Number Description		
MX405	5 mm (10/Box)	
MX406	6 mm (10/Box)	
MX407	7 mm (10/Box)	

WITH CANNULA CURETTE — SEMI-RIGID	
Number	Description
MX415	5 mm (10/Box)
MX416	6 mm (10/Box)
MX417	7 mm (10/Box)

WITH CANNULA CURETTE – RIGID			
Number	Description		
MX448	8 mm (10/Box)		
MX449	9 mm (10/Box)		
MX450	10 mm (10/Box)		
MX451	11 mm (10/Box)		
MX452	12 mm (10/Box)		

Traps Only

Number	Description	
MX400CONE	Traps only — Cone Shape (10/Pack)	
MX400FLAT	Traps only — Flat Shape (10/Pack)	





Disposable Sound and Dilators

Fundus Sound™ XM



The Fundus Sound XM's clear depth markings and centimeter notches present a quick, easy-to-read display for fast and precise measurement of uterine depth. The disposable sound features a small diameter for greater patient comfort. Made from an advanced polymer material, the sound is malleable and offers the correct amount of pliancy to follow the curve of the cervix and uterus for quick and efficient procedures. Simply dispose of the sound after each use.

Number	Description	
1196	Fundus Sound XM Disposable Uterine Sound (25/Box)	





The Os Finder naturally seeks and gently opens the cervical Os. Made of an extremely smooth and malleable Teflon-like material, this instrument is ideal for use in intracervical dilation or uterine sounding, or for intrauterine procedures.

Number	Description
1176	Os Finder (50/Box)

Milex® Dilator / Sound Set	CA *



The Milex Dilator/Sound Set devices are single-use, disposable products that are ergonomically designed and engineered to provide superb balance and heightened tactile sensation. Specifically designed to safely facilitate: rapid location of the cervical Os, conservative endocervical dilation, gradual opening of the internal Os, accurate determination of uterine depth.

Number	Description
MX11	Milex Os Locator (25/Box)
MX12	Milex Canal Dilator (25/Box)
MX13	Milex Fundal Sound (25/Box)
MX21	Milex Dilator/Sound Set, 3 dilators per set: MX11, MX12, MX13 (10/Box)



Dilateria is ideal for gentle, atraumatic cervical dilation. It is a by-product of Laminaria Japonica, a natural plant from Northern Asia. The plant's stipe, or stem, is dried, sized and subsequently hand-processed to ensure proper rounding of the tip and a more uniform diameter. A string is attached for ease of location and removal, while a disc placed at the end minimizes further migration into the cervical canal.

Number	Description	Size
MX210	Dilateria Thin	3 mm
MX220	Dilateria Medium	4 mm
MX230	Dilateria Thick	5 mm
MX240	Dilateria Extra-Thick	6 mm

LEEP

LEEP PRECISION[™] Integrated System

The Preferred Choice for Functionality, Performance, and Precision for Loop Electrosurgical Excision Procedures

Fully integrated and mobile, the LEEP Precision System combines three components — Cart, Generator and Smoke Evacuator — for maximum efficiency.

- Hand or foot activation for operating precision
- Smoke evacuator eliminates odor at 99.999 percent efficiency level
- Fully integrated LEEP station with storage
- Safe and easy to operate
- Heavy-duty casters facilitate mobility
- Elegantly designed for the modern procedure room

NOW AVAILABLE THROUGH DISTRIBUTION













LEEP PRECISION Inte	grated System

Number	Description
LP-10-120	LEEP PRECISION Integrated System, 120V
LP-10-220	LEEP PRECISION Integrated System, 220V

LEEP PRECISION Generator

Number	Description
LP-20-120	LEEP PRECISION Generator, 120V
LP-20-220	LEEP PRECISION Generator, 220V

LEEP PRECISION Smoke Evacuator

Number	Description
LP-30-120	LEEP PRECISION Smoke Evacuator, 120V, Stand-Alone
LP-30-220	LEEP PRECISION Smoke Evacuator, 230V, Stand-Alone

LEEP PRECISION Cart

Number	Description
LP-40-001	LEEP PRECISION Cart

LEEP PRECISION Generator Disposables

Number	Description
LP-50-101	Hand Switch-Operated Push-Button Cautery Pencils (10/Box)
LP-50-102	Foot Pedal-Operated Cautery Pencils (10/Box)
LP-50-103	Foot Pedal Adapter (1/Bag)
LP-50-201	Dual Dispersive Electrodes with Connector (10/Box)
LP-50-202	Dual Dispersive Electrodes without Connector (Connect to Reusable Cable) (10/Box)
LP-50-203	Reusable Cable (1/Bag)

LEEP PRECISION Smoke Evacuator Disposables

Number	Description
52560	Smoke Evacuator Pre-Filter, Single-Patient Use (10/Bag)
6082	Smoke Evacuator ULPA Filter (1/Box)
6083	Smoke Evacuator Reducer, Single-Patient Use (10/Bag)
6084	Smoke Evacuator Tubing — 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box)

LFFP Flectrodes

Fischer Cone Biopsy Excisor®

Improves Control During Cone Biopsy Procedures




An innovative LEEP electrode, the Fischer Cone Biopsy Excisor (FCBE) comes in a variety of sizes and shapes for precise cone biopsy specimens.

- Smooth, 360-degree rotating motion removes a perfect, cone-shaped cervical specimen in seconds
- Linear wire electrode stays rigid and makes a clean, precise cut with little or no bleeding, and less thermal artifact for easier pathology evaluation and simplified tissue analysis
- Stop Arm provides firm support and added stability for the wire electrode, resulting in better control and smoother cutting
- Arm curvature permits precise instrument positioning and allows for unobstructed visualization during procedure
- Shaft and Stop Arm are fully insulated to protect endocervical canal and underlying vaginal structures

Fischer C	Cone Biopsy Excisor	CA* (€
Number	Description	
900-150	Fischer Cone Biopsy Excisor, Small (5/Box)	
900-151	Fischer Cone Biopsy Excisor, Medium (5/Box)	
900-152	Fischer Cone Biopsy Excisor, Large (5/Box)	
900-154	Fischer Cone Biopsy Excisor, Medium Extended	(5/Box)
900-155	Fischer Cone Biopsy Excisor, Large Shallow (5/8	Box)
900-156	Fischer Cone Biopsy Excisor, Mixed Set (S, M, ME, L & LS) 1 each (5/Box)	
900-157	Fischer Cone Biopsy Excisor, Small Wide-Angle	(5/Box)
900-158	Fischer Cone Biopsy Excisor, Small Wide-Angle Extended (5/Box)	


LEEP Electrodes (Shown Actual Size)

These original tungsten wire electrodes provide exceptional quality and performance. Our complete selection of 21 different offerings includes loops, squares, balls and needle tips specifically designed for loop electrosurgical excision procedures. Made of the finest tungsten wire for minimal charring, these universal electrodes have a standard ½-inch diameter brass connector shaft, for optimal conductivity, that can be used on most qualified generators.





2.0 cm W x 0.8 cm D • 12 cm shaft (5/Box)

R2008



LEEP Electrodes (Shown Actual Size)





LEEP Procedure Essentials



LEEP RediKit includes everything you need to prepare and care for your patient before, during and after the procedure. This single-use kit is set up as an operating tray, organized to present procedure items as you need them. For sale in the USA only.

1 Pantyliner

• 1 Underpad

*Sterile

Each LEEP RediKit contains the following:

- 1 Vial Xylocaine (Lidocaine HCI 1%) 10 ml with Epinephrine*
- 1 Potocky Needle® 27 G (for administration of cervical block)* 1 Alcohol Wipe*
- 1 Lugol's (Strong lodine Solution, USP) 8 ml, Single Dose
- 1 Lugor 3 (Strong louine Solution, OSI / O IIII, Single Dose
- 1 AstrinGyn® (Ferric Subsulfate), 8 grams, Single Dose
- 6 Applicator Swabs, 8 inch
- 20 G Needle
- 12 cc Luer Lock Syringe*

Number	Description
6061	LEEP RediKit, Single-use, disposable. Contact Customer Service for ordering information.

Potocky Needle®



The 27-gauge Potocky Needle was specifically designed to solve the problem of applying injections to cervices of varying length and density. Its reinforced and extended shaft makes cervical application an easy process. The Potocky Needle attaches to a "dental style" cartridge syringe, while the Potocky Needle with Luer Lock is used with a standard medical syringe.

Number	Description
6066	Potocky Needles (6/Box)
6066-4	Potocky Needles – 4 Boxes of 6 (24/Case)
6166	Potocky Needles with Luer Lock (6/Box)
6166-4	Potocky Needles with Luer Lock — 4 Boxes of 6 (24/Case)

Needle Extension



Number	Description
433-132	Needle Extension • 10 cm (1/Box)
433-133	Needle Extension • 13 cm (1/Box)
433-134	Needle Extension • 15 cm (1/Box)

AstrinGyn[®]



AstrinGyn — CooperSurgical's thickened Monsel's solution — quickly controls oozing and bleeding for colposcopic and LEEP procedures. The gel consistency makes this styptic product both easy to apply and essential to your practice. **For sale in the USA only.**

Number	Description
6065	AstrinGyn (12/Box)
6065-4	AstrinGyn – 4 Boxes of 12 (48/Case)

Lugol's



Lugol's strong iodine solution can be used during colposcopy to assist in recognition and delineation of cervical lesions. Conveniently packaged in single-use vials; both beginning and experienced colposcopists appreciate the added visualization. For sale in the USA only.

Number	Description
6064	Lugol's. Contact Customer Service for ordering information. (12/Box)
6064-4	Lugol's. Contact Customer Service for ordering information. (48/Case)

LEEP System Accessories



These nonsterile patient return pads are used to complete the electrical circuit between the patient and the electrosurgical unit. The large size provides optimal safety for your patient.

Number	Description
6050P1	LEEP Return Pad, Single-Patient Use (10/Box)

CA* (€

CA***** (€



This reusable adapter connects the single-use handpiece to the LEEP System 1000 and other generators.

Number	Description
22178	Handpiece Adapter for LEEP System 1000

Foot Pedal, LEEP System 1000

The foot pedal for the LEEP System 1000 is designed for exceptional performance and reliability during the LEEP procedure. Featuring pneumatic, air-operated control, the foot pedal does not use electric current, thus offering maximum safety for you and your patient.

Number	Description
6032	Foot Pedal for LEEP System 1000

Handpiece, Single-Patient Use CA◆



The sterile, single-patient use handpiece works with any of our disposable electrodes. The wand is operated with a footswitch for ultimate hand precision.

Number	Description
6040	Handpiece, Single-Patient Use (10/Box)

Patient Return Pad Adapter LEEP System 1000 CA★ C€



This reusable replacement adapter connects the single-use Patient Return Pad to the LEEP System 1000.

Number	Description
23635	Patient Return Pad Adapter for LEEP System 1000 (for use with 6050P1 pads)

Disposable Smoke Evacuator (DSE) Tubing CA★ C€



Single-use, nonsterile DSE tubing is designed to work specifically with the LEEP rose-coated specula. The 9-inch DSE tubing is flexible and noncombustible and contains no harmful PVC. Simply slide the tubing through the DSE tube guide on your compatible instrument and attach to flexible "S" tubing to complete the setup. Discard used tubing after each procedure.

Also available as a set with 6 feet of tubing. See Tubing Set, Single-Patient Use, product item number 6084.

Number	Description		
6150	Disposable Smoke Evacuator Tubing — 9" straw (50/Bag)		
34371	LEEP System 1000 Manual (No charge)		

Pre-Filter CA♦ (€



Replacement filters for the CooperSurgical Smoke Evacuator. Used to eliminate smoke odor.

Number	Description
52560	Smoke Evacuator Pre-Filter (10/Box)

ULPA Filter CA★ C€



Used with the Smoke Evacuator for maximum particle capture. These highly rated filters capture particles down to 0.14 microns. Effective for approximately six months.

Number	Description
6082	Smoke Evacuator ULPA Filter (1/Box)

Tubing Set, Single-Patient Use CA◆



Each set includes 6 feet of sterile disposable tubing and a DSE (Disposable Smoke Evacuation) tube. This set completes the connection between the reducer and the specula.

Number	Description
6084	Smoke Evacuator Tubing — 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box)

Reducer CA***** C€



Disposable, single-patient use reducer connects the Pre-Filter (52560) with the Patient Tubing Set (6084) on the Smoke Evacuator.

Number	Description
6083	Smoke Evacuator Reducer (10/Box)

LEEP Accessories


Mini-Vac - Stackhouse Flexible Tubing CA♥ C€



Stackhouse Flexible Tubing is intended for brief procedures generating minimal amounts of smoke.

Number	Description
395-703	Smoke Evacuator Filters for MiniVac Stackhouse Evacuator (3/Case)
395-704	Smoke Evacuator Tubing Assembly; Corrugated Tubing Reducer. ¾" Tubing (3/Box)

Compatible Disposable Accessories



The following items are compatible with most of the competitive products in the marketplace. See your local sales representative for more information.

Number	Description	
395-210	Electrode Tip Cleanser, Disposable (24/Box)	C€
395-302	Conductive Electrode Pads, Disposable for ESU (2	25/Box)
395-305	Conductive Electrode Pads, Dual Dispersive, Disposable (25/Box) Plus II Gen.	
395-310	Hand-Switch Pencil and Cord, Disposable (25/Bo	x)
72308	Dual Section Grounding Pads (10/Box)	
72102	Filter/Hose Set for Smoke Evacuator (Pack of 3)	C€
72311	Patient Return Cable	CA * (€
22192	LEEP 6000 Return Pad	
22191	Foot Pedal	


Euro-Med[®] LEEP Instruments


These nonconductive, fully autoclavable LEEP instruments are offered in a variety of styles to facilitate cervical tissue excision and removal.

The first choice for proven performance, CooperSurgical provides the broadest selection of fully autoclavable instruments for LEEP procedures. Our "rose-coated" instruments will provide you with durability and performance unmatched in the industry.





Pederson LEEP Specula Designed specifically for electrosurgical applications, Pederson's narrow blades are more comfortable for smaller patients or those with a narrow vaginal canal. Number Description Size mm F223 Small W13 L81 W22 L114 F221 Medium F220 Large W25 L120





Prima Series® Vaginal LEEP Specula Prima Series LEEP Specula are made from a high-density polymer, so there is no coating. That means no chipping, no cracking and no recoating. Ever. (Note: this is the only Euro-Med Speculum not crafted from stainless steel.) Number Description Size mm PSV-1L W33 L113 Medium W35 L125 PSV-2L Large





(W=Distal End) (L=Penetration Length)

Mini Vu-More™ LEEP Specula

This speculum features very short blades to facilitate examination of patients with cervical prolapse or very short vaginal canals.

Number	Size mm
F250	W34 L93



CA*

Auvard Surgical Specula

This 2.5-pound weighted speculum can be used safely in any operative vaginal procedure where the use of electrosurgery is anticipated.

Number	Description	Size mm
F230	Medium (Wide/ Short Blade)	W43 L55
F232	Large (Narrow/ Long Blade)	W35 L90

y CA♥

Cer-View™ LEEP Lateral Wall Retractor

Retracting redundant vaginal tissue, Cer-View provides increased exposure of the cervix, especially important when vaginal prolapse is an issue or with nulliparous patients. Contoured and shaped to fit into the yoke of a standard speculum and expand beyond the width of the speculum blades for an unobstructed view.

Number	Description	Size mm
F410	Medium	W22 L142
F420	Large	W22 L155



CA*

Graves Specula (with Integrated Smoke Tube)

This popular vaginal speculum features the standard Graves blades and incorporates an integrated smoke tube. Available in three sizes.

Number	Description	Size mm
61513	Small	W22 L87
61514	Medium	W33 L113
61515	Large	W35 L125

Open-Sided LEEP Specula

This vaginal speculum features a sturdy one-sided hinge arrangement for improved instrument access and visibility.

Number	Size mm	
F210	W35 L107	

Jackson LEEP Retractor

This coated retractor can be inserted and pulled downward to provide maximum posterior viewing.

Number	Description	Size mm
F235	Medium	W40 L89
F234	Large	W40 L108

Tru-View™ LEEP Lateral Wall Retractor

CA

CA*

Designed especially for the obese patient, the Tru-View features longer blades to safely expose and access the cervix. Provides added strength and support under prolapse conditions. The perfect complement to the Snowman Speculum.

Number	Description	Size mm
F317	Lateral Wall Retractor	W32 L190



Pederson LEEP Specula (with Integrated Smoke Tube) c

This model is specifically modified for electrosurgical applications. The narrow blades are more comfortable for your smaller patients or those with a narrow vaginal canal. Features integrated smoke tube. Available in three sizes.

Number	Description	Size mm
61521	Small	W16 L84
61518	Medium	W22 L114
61519	Large	W26 L121



LEEP Tenacula, Hooks and Rulers

Instruments are checked for appropriate coating thickness and proper insulation to ensure maximum patient safety during electrosurgical procedures.

LEEP Car	npion Forceps	CA*

Angled shaft and tip provide better visualization of the target tissue being grasped or removed.

Number	Instrument Length
F925	220 mm



The angle at which the prongs are set permits easy manipulation of the cervix while avoiding obstruction of the electrode.

Number	Instrument Length
F802	260 mm



Select this small hook when stabilization of tissue is required during excision.

Number	Instrument Length	
F810	Straight	250 mm
F812	Left-Angled	250 mm
F813	Right-Angled	250 mm

LEEP Ring Forceps CA*



Long-handled forceps facilitate manipulation and retrieval of tissue during and after the procedure.

Number	Instrument Length
F920	220 mm



These long tenacula provide maximum cervical stabilization through their two opposing hooks.

Number	Instrument Length
F900	240 mm

LEEP Three-	Prong Hook	CA*


The prong arrangement is designed to allow complete access, maximum tissue manipulation and stabilization of the cone during the procedure.

Number	Instrument Length
F803	260 mm

LEEP D	epui ni	ilei	CA
			 _
and the same of th			

Ensures accuracy of tissue removal during loop excision with easy-to-read markings.

Number	Instrument Length
F600	240 mm




These tissue forceps are perfect for safe tissue manipulation at the surgical site.

Number	Working Length
F700	210 mm



This larger hook affords added stability to the cervix during LEEP procedures.

Number	Instrument Length	
F805	250 mm	



This popular small hook stabilizes the tissue and allows maximum electrode access.

E800 263 mm	Number	Instrument Length	
203 11111	F800	263 mm	



The easy-to-read millimeter markings are ideal for measuring the width of the excised sample.

Number	Instrument Length	
F610	210 mm	

LEEP Endocervical Specula

Based on our most popular ratchet handle design, these "rose-coated" instruments facilitate accurate examination of the endocervical canal for complete patient care. The 5 mm tip is recommended for the multiparous woman.



Number	Tip Size	Length
F500	5 mm	273 mm




CA*

Number	Tip Size	Length
F510	4 mm	273 mm

Cervical Cryosurgery

LEISEGANG® LM-900 Cryosurgery System

High Performance and Low Cost



Number	Description	
50501	LM-900 N₂O Cryosurgery Unit with Case	C€
50502	LM-900 CO ₂ Cryosurgery Unit with Case	C€

ACCESSORIES

Number	Description	
50000	Empty 20 lb. N₂O	C€
50010	Empty 20 lb. CO ₂	C€
50003	20 lb. Cylinder Cart	C€
50007	E-Cylinder Adapter N₂O	C€

50008	E-Cylinder Adapter CO ₂	CE
50009	Carrying Case	CE
50011	O-Rings (10 μ) CA *	CE
50012	Cleaning Plugs (1 µ) CA◆	CE

LEISEGANG LM-900 Probes













Number	Description
50101	Microderm 2 mm D x 25 mm L

Number	Description	





(
-		

	2 coci.ip tioi.i
50101	Microderm 2 mm D x 25 mm L

Number	Description
50102	Soll Lash Epilating
	2 mm D v 27 mm I

Number	Description
50103	Condyloma
	5 mm D x 25 mm l

Num	ber	Description
5010	4	Derm Bevel 5 mm D x 24 mm L

Number	Description
50105	General Purpose
	19 mm D x 8 mm L











50106 Large Derm Bevel

50107	Curved Anorectal
	8 mm D x 38 mm

50108 General Purpose 19 mm D x 8 mm L













Special Exo/Endo 19 mm D x 10 mm L Exo/Endocervical 19 mm D x 20 mm L

Maxiderm

Large Exocervical 25 mm D x 7 mm L

Jumbo Exo/Endo 25 mm D x 14 mm L

Frigitronics® Cryo-Plus™

Engineered to Provide Effective Cryosurgical Treatment



- Excellent thermal conductivity
- Temperature gauge for both cryo probe tip and thermocouple ensures adequate freeze with each treatment
- Rapid freeze and quick defrost
- Accurate and easy-to-read pressure gauge
- Probe shaft connects at the handle for added patient safety
- Wide choice of probes for every lesion and location
- Console with integrated storage area for five tips
- Autoclavable probes and shaft

Ease of Operation

The Frigitronics Cryo-Plus System has everything you need to perform cryosurgery procedures, from excellent probe thermoconductivity to temperature- and pressure-control gauges. By monitoring the actual tip temperature with every freeze, you are assured an adequate performance with each treatment.

Designed for safety and portability, the Frigitronics Cryo-Plus probes are interchangeable and autoclavable for improved equipment turnaround and to eliminate the potential for cross-contamination.

Gas Tar	nks, Carts, and Accessories	
Number	Description	CA* (€
659	"E" Cylinder, Empty	CA * (€
50000	20 lb. "G" Cylinder, Empty	C€
679	Double "E" Cylinder Cart	C€
681	Single "E" Cylinder Cart	CA * (€
694	Single "G" Cylinder Cart	CA * (€
32873	Exhaust Tube, Rear 20'	CA* (€
2475	Procedural Timer	CA * C€
2476	Probe Storage Clip	CA * (€
2477	Thermocouple Needle (blue) 22-gauge x ¾"	CA * (€
P190	Gas Purifier with One Cartridge (Filter)	CA * (€
P191	Replacement Cartridges for P190 (6/Package)	CA* (€

Cryo-Surg® Handheld System

Cryo-Surg Tubing

CA* (€

This lightweight, non-electric system provides a continuous uniform freeze, with specialized valve and anti-clog tips to ensure no freezing, dripping, or clogging. The self-pressurizing unit maintains a constant low pressure for immediate use. The unique handle is made from a thermally non-conductive polymer to eliminate hand freeze, and is ergonomically angled to reduce hand and wrist fatigue. Available with a wide assortment of interchangeable Luer Lock tips.

Number	Description
P5900	Cryo-Surg Liquid Nitrogen Handheld System with 19-gauge Spray Needle, Exhaust Tube and Spare "O" Ring
ACCESSO	RIES
52240	Exhaust Tubing (20" with Velcro)



20188

Various cart configurations allow the system to be used in numerous practice locations (office and hospital). The system has a 5-probe storage capacity.

Frigitronic	s Cryo-Plus Systems
Number	Description
2400	Cryo-Plus with 20 lb. "G" cylinder yoke
2401	Cryo-Plus with single "E" cylinder yoke
2402	Cryo-Plus with double "E" cylinder yoke



All Probe Cannulas Have Scientifically Measured Air Gap Insulation Within the Cannula

Frigitronics Cryo-Plus Gynecology Probes

No.	Description	No.	Description	No.	Description	No.	Description	No.	Description
2440	General Purpose Probe 6 mm D x 6.5 mm L	2441	Micro Probe 2 mm D x 4.5 mm L	2442	Small Cervix Probe 19 mm D x 9 mm L	2443	Large Cervix Probe 25 mm D x 8 mm L	2444	Small Cone Probe 19 mm D x 15 mm L
=		_		_				_	

Cryo-Surg® Spray Cone Tips



Number	Description
5212	10 mm Spray Cone
5213	15 mm Spray Cone
5214	20 mm Spray Cone
5215	25 mm Spray Cone

Cryo-Surg Spray Needle Tips



Number	Description	
515	18-gauge Spray Needle	CA* (€
516	19-gauge Spray Needle	CA* (€
5232	20-gauge Spray Needle	
5233	22-gauge Spray Needle	
5234	24-gauge Spray Needle	

Cryo-Surg Closed Probe Tips



Number	Description
5235	1.5 mm Closed Probe
5236	3 mm Closed Probe
5237	5 mm Closed Probe
5238	10 mm Closed Probe

Ophthalmic Cryosurgery

Frigitronics® CE-2000™

Continuing to Set the Cold Standard in Your Practice

NOW AVAILABLE THROUGH DISTRIBUTION

The Frigitronics CE-2000 is a completely non-electric, low-pressure, fully exhausted cryosurgical system that uses nitrous oxide gas. The system destroys tissue by applying extreme, controlled temperature down to a maximum of -89° C.

The CE-2000 consists of a control console with foot pedal, mobile equipment carrier, nitrous oxide cylinders and cryosurgical probes. Depress the foot pedal to rapidly freeze and release to defrost. The backlight display activates automatically when the foot pedal is depressed.

An automatic shut-off built into the system will stop the flow of nitrous oxide gas and will shut off the backlight and display approximately four minutes after the probe has been set aside. The system can be reactivated by simply depressing the foot pedal. A wide range of probes are available for use with the CE-2000

Number	Description	
CE-2000	CE-2000 Control Console with Foot Pedal	CA * (€

Gas Tanks			
Number	Description		
659	"E" Cylinder, Empty	CA * (€	
50000	20 lb. "G" Cylinder, Empty		
20370-3	High-Pressure Hose (Synflex – black) with Fittings	CA * (€	

Yoke A	ssemblies for "E" and "G" Cylinders	
Number	Description	
5088	"E" Cylinder Single-Yoke Assembly with Filter	CA * (€
5093	"E" Cylinder Double-Yoke Assembly with Filter	CA €
5014	"G" Cylinder Single-Yoke Assembly	CA * (€
5087	"G" Cylinder Single-Yoke Assembly with Filter	CA* (€
20917	Retrofit "E" Cylinder Yoke	CA* (€
51746	Retrofit "G" Cylinder	CA €
52038	"G" Cylinder Double-Yoke Assembly with Filter	CA* (€





Carts and Carriers for Cylinders		
Number	Description	
20566	Exhaust Tube, Front 6"	
32873	Exhaust Tube, Rear 20'	

Exhaust Output				
Number	Description			
503	Office Cart for Single "G" 20 lb. Cylinder with Yoke Assembly	CE		
504	Office Cart for Single "G" 20 lb. Cylinder with Yoke Assembly with Filter	C€		
506	Mobile Equipment Carrier (with built-in "E" Cylinder Double-Yoke Assembly with Filter)	(€		
507-1A	Mobile Equipment Carrier (with built-in "G" 20 lb. Cylinder Double-Yoke Assembly with Filter)	C€		

Frigitronics® CE-2000™ Probes

Supe	Supercool Retinal and Glaucoma Probes CA CA			
No.	Description			
142	Posterior Retinal – 2.7 mm tip x 2.7 mm dia. x 38 mm shaft			
161	Curved Retinal – 2.5 mm tip x 2.5 mm dia. x 38 mm shaft			
162	Straight Retinal — 2.5 mm tip x 2.5 mm dia. x 38 mm shaft			
163	Spherical Retinal – 3.0 mm sphere x 2.1 mm dia. x 38 mm shaft			
164	Straight Glaucoma – 3.5 mm tip x 3.5 mm dia. x 38 mm shaft			
165	Curved Glaucoma – 3.5 mm tip x 3.5 mm dia. x 38 mm shaft			
166	3.2 mm Straight Retinal — 3.2 mm tip x 3.2 mm dia. x 38 mm shaft			
167	3.2 mm Curved Retinal – 3.2 mm tip x 3.2 mm dia. x 38 mm shaft			

Stan	dard Retinal Probes CA* CE
No.	Description
124	Straight Retinal – 2.5 mm tip x 2.5 mm dia. x 38 mm shaft
125	Curved Retinal – 2.5 mm tip x 2.5 mm dia. x 38 mm shaft
139	Spherical Retinal – 3.0 mm sphere x 2.1 mm dia. x 28 mm shaft
143	Hammerhead, Curved -3.5 mm tip x 2.5 x 6.5 mm dia. x 40 mm shaft
144	Hammerhead, Straight – 3.5 mm tip x 2.5 x 6.5 mm dia. x 40 mm shaft

Cataract Probes CA★ C€		
No.	Description	
126	Straight $-$ 1.0 tip x 1.8 mm tip dia. x 27 mm shaft x 1.8 mm shaft dia.	
127	Curved – 1.0 tip x 1.8 mm tip dia. x 27 mm shaft x 1.8 mm shaft dia.	
128	Micro Straight – 1.0 mm tip dia. x 27 mm shaft x 1.0 mm shaft dia.	
147	$\label{eq:micro} \mbox{Micro Curved} - 1.0 \mbox{ mm tip dia.} \ \mbox{x 27 mm shaft x 1.0 mm shaft dia.}$	

Trichiasis Probes CA CA		
No.	Description	
171	18-gauge Cone Probe – 19 mm shaft x 1.3 mm shaft dia.	
172	3 mm Multi-Lash — 19 mm shaft x 2.8 mm shaft dia. x 5 mm tip x 2.7 mm tip	dia.
173	6 mm Multi-Lash — 19 mm shaft x 2.8 mm shaft dia. x 5 mm tip x 2.7 mm tip d	ia.

Pedia	Pediatric Retinal Probes CA★ C€		
No.	Description		
120	1.8 mm Retinal Curved — 1 mm tip x 2.54 mm dia. x	26 mm shaft	
121	2.08 mm Retinal Curved — 1 mm tip x 2.84 mm dia. x	27 mm shaft	
123	3 mm Disc Probe – 3 x 1.5 mm tip x 1.5 mm dia. x 3	5 mm shaft	

Vitreous Intraocular Probes CA◆ C€			
No.	Description		
134	20-gauge	35 mm x .89 mm diameter	

	u' el est co			
insula	Insulating Sleeves CA★ C€			
No.	Description			
P805	For Probes 124/125/142/161/162/166/167 (Retinal) (6/Pkg)			
805-1	For Probes 143/144 (Hammerhead) (6/Pkg)			
32709	For Probes 128/132/147 (Micro) (1/Pkg)			
32710	For Probes 126/127 (Cataract) (6/Pkg)			
32787	For Probes 139/163 (Retinal Spherical) (1/Pkg)			
33098	For Probes 130/164/165 (Glaucoma) (6/Pkg)			
33555	For Probes 120/121 (Retinal Pediatric) (1/Pkg)			
33966	For Probe 122 (Retinal Spherical Pediatric) (1/Pkg)			

Milex[®] Pessaries

Making Life Easier for You and Your Patient

CooperSurgical offers:

- Office inventory management ensure most commonly used pessaries and quantities in stock
- Superior knowledge of your account reduces over-stock issues
- Experts in specific indications for which the different pessaries are utilized

 In-service/re-training on the different indications and patient profiles

• Excellent product quality and consistency

Superior Customer Service

Patient education materials

Milex® Pessaries Fitting Kit........... 57

Milex® Pessaries for

Prolapse and Incontinence...... 58



Milex® Pessaries

Our Pessaries Are Autoclavable

Milex Pessaries are used to help patients with pelvic organ prolapse (POP) and urinary incontinence. They come in 25 different styles and many sizes to fit each patient's anatomy. Made of soft, autoclavable, medical-grade silicone, Milex Pessaries are an excellent addition to any practice that treats women with lower urinary tract dysfunction.

Milex® Pessary Fitting Kit



The Milex Pessary Fitting Kit was developed to help you in your practice to determine accurate sizing prior to prescribing a pessary for your patient. Note: for fitting only, not for long-term patient use.

Number	Description
MXFIT	Pessary Fitting Set (For sale in the USA only)
MXFIT-CO	Pessary Fitting Set Case and Tray only (no pessaries)
MXFIT-CH	Pessary Fitting Set (China) Note: Limited selection of Gellhorns, Rings with support and Incontinence Dish without Trimo-San
MXFIT-CA	Pessary Fitting Set (Outside the USA) Note: As seen in above photo but without Trimo-San

ACCESSORIES	CA* (€

Number	Description	Number	Description
MXFIT0016	Fitting Gellhorn Short Stem Flexible 21/4"	MXFIT0801	Fitting Incont. Dish with support #2
MXFIT0017	Fitting Gellhorn Short Stem Flexible 21/2"	MXFIT0802	Fitting Incont. Dish with support #3
MXFIT0018	Fitting Gellhorn Short Stem Flexible 2¾"	MXFIT0803	Fitting Incont. Dish with support #4
MXFIT0019	Fitting Gellhorn Short Stem Flexible 3"	MXFIT0804	Fitting Incont. Dish with support #5
MXFIT0012	Fitting Gellhorn Flexible 21/4"	MXFIT0805	Fitting Donut #2
MXFIT0013	Fitting Gellhorn Flexible 2½"	MXFIT0806	Fitting Donut #3
MXFIT0014	Fitting Gellhorn Flexible 2¾"	MXFIT0807	Fitting Donut #4
MXFIT0015	Fitting Gellhorn Flexible 3"	MXFIT0808	Fitting Donut #5
MXFIT0020	Fitting Ring with support #2	MXFIT0809	Fitting Cube #2
MXFIT0021	Fitting Ring with support #3	MXFIT0810	Fitting Cube #3
MXFIT0022	Fitting Ring with support #4	MXFIT0811	Fitting Cube #4
MXFIT0023	Fitting Ring with support #5	MXFIT0812	Fitting Cube #5
MXFIT0024	Fitting Ring with support #6	MXFIT0072	Pelvic Teaching Model
MXFIT0025	Fitting Ring with support #7		

Trimo-San[™] – Vaginal Jelly

CA**₩** (€

Trimo-San helps restore and maintain the normal vaginal acidity and coats the walls of the vagina with a lubricating film that helps reduce odorcausing bacteria. Trimo-San is a recommended adjunct for pessary wearers. For sale in the USA only.



Number Description MX5030 Trimo-San Cream with Jel-Jector (MXAPP Jel-Jector replacement at no charge)

Milex® Pessaries for Prolapse and Incontinence

P = Prolapse I = Incontinence

Ring

For mild prolapse.

Product # MXPER (Pessary Only) Product # MXKPER (Pessary, Trimo-San, Jel Jector)



Ring				
Product #				
MXPER00	13/4	44		
MXPER01	2	51		
MXPER02	21/4	57		
MXPER03	21/2	64		
MXPER04	23/4	70		
MXPER05	3	76		
MXPER06	31/4	83		
MXPER07	31/2	89		
MXPER08	3¾	95		
MXPER09	4	102		
MXPER10	41/4	108		
MXPER11	41/2	114		
MXPER12	43/4	121		
MXPER13	5	127		

Ring with Support

For mild prolapse complicated by a mild cystocele.

Product # MXPRS (Pessary Only) Product # MXKPRS (Pessary, Trimo-San, Jel Jector)





Ring with Support			
Product #			
MXPRS00	13/4	44	
MXPRS01	2	51	
MXPRS02	21/4	57	
MXPRS03	21/2	64	
MXPRS04	23/4	70	
MXPRS05	3	76	
MXPRS06	31/4	83	
MXPRS07	31/2	89	
MXPRS08	3¾	95	
MXPRS09	4	102	
MXPRS10	41/4	108	
MXPRS11	41/2	114	
MXPRS12	43/4	121	
MXPRS13	5	127	

Ring with Knob

For mild prolapse complicated by stress urinary incontinence.

Product # MXPERK (Pessary Only) Product # MXKPERK (Pessary, Trimo-San, Jel Jector)



Ring with Knob				
Product #		mm		
MXPERK00	13/4	44		
MXPERK01	2	51		
MXPERK02	21/4	57		
MXPERK03	21/2	64		
MXPERK04	23/4	70		
MXPERK05	3	76		
MXPERK06	31/4	83		
MXPERK07	31/2	89		
MXPERK08	33/4	95		
MXPERK09	4	102		
MXPERK10	41/4	108		
MXPERK11	41/2	114		
MXPERK12	43/4	121		
MXPERK13	5	127		

Incontinence Ring

For stress urinary incontinence.

Product # MXPCON (Pessary Only) Product # MXKPCON (Pessary, Trimo-San, Jel Jector)



Incontinence Ring					
Product #					
MXPCON00	13/4	44			
MXPCON01	2	51			
MXPCON02	21/4	57			
MXPCON03	21/2	64			
MXPCON04	23/4	70			
MXPCON05	3	76			
MXPCON06	31/4	83			
MXPCON07	31/2	89			
MXPCON08	33/4	95			
MXPCON09	4	102			
MXPCON10	41/4	108			
MXPCON11	41/2	114			
MXPCON12	43/4	121			
MXPCON13	5	127			

Inflatoball™

Inflatable for 3rd degree prolapse/procidentia with or without a mild cystocele and/or rectocele.

Product # MXPINF (Pessary Only) Product # MXKPINF (Pessary, Trimo-San, Jel Jector)



Inflatoball–Size				
Product #				
MXPINFS	2	51		
MXPINFM	21/4	57		
MXPINFL	21/2	64		
MXPINFXL	23/4	70		

Ring with Support & Knob

For mild prolapse complicated by a mild cystocele and stress urinary incontinence.

Product # MXPRSK (Pessary Only) Product # MXKPRSK (Pessary, Trimo-San, Jel Jector)



Incontinence Ring W/ Support & Knob Product # in. mm MXPRSK00 1¾ 44 MXPRSK01 2 51 MXPRSK02 2¼ 57 MXPRSK03 2½ 64 MXPRSK04 2¾ 70 MXPRSK04 2¾ 70 MXPRSK05 3 76 MXPRSK06 3¼ 83 MXPRSK07 3½ 89 MXPRSK07 3½ 89 MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114 MXPRSK11 4½ 114							
MXPRSK00 1¾ 44 MXPRSK01 2 51 MXPRSK02 2½ 57 MXPRSK03 2½ 64 MXPRSK04 2¾ 70 MXPRSK05 3 76 MXPRSK06 3¼ 83 MXPRSK07 3½ 89 MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114							
MXPRSK01 2 51 MXPRSK02 2½ 57 MXPRSK03 2½ 64 MXPRSK04 2¾ 70 MXPRSK05 3 76 MXPRSK06 3¼ 83 MXPRSK07 3½ 89 MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114	Product #	in.	mm				
MXPRSK02 2¼ 57 MXPRSK03 2½ 64 MXPRSK04 2¾ 70 MXPRSK05 3 76 MXPRSK06 3¼ 83 MXPRSK07 3½ 89 MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114	MXPRSK00	13/4	44				
MXPRSK03 2½ 64 MXPRSK04 2¾ 70 MXPRSK05 3 76 MXPRSK06 3¼ 83 MXPRSK07 3½ 89 MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114	MXPRSK01	2	51				
MXPRSK04 2¾ 70 MXPRSK05 3 76 MXPRSK06 3¼ 83 MXPRSK07 3½ 89 MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114	MXPRSK02	21/4	57				
MXPRSK05 3 76 MXPRSK06 3¼ 83 MXPRSK07 3½ 89 MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114	MXPRSK03	21/2	64				
MXPRSK06 31/4 83 MXPRSK07 31/2 89 MXPRSK08 33/4 95 MXPRSK09 4 102 MXPRSK10 41/4 108 MXPRSK11 41/2 114	MXPRSK04	23/4	70				
MXPRSK07 3½ 89 MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114	MXPRSK05	3	76				
MXPRSK08 3¾ 95 MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114	MXPRSK06	31/4	83				
MXPRSK09 4 102 MXPRSK10 4¼ 108 MXPRSK11 4½ 114	MXPRSK07	31/2	89				
MXPRSK10 4¼ 108 MXPRSK11 4½ 114	MXPRSK08	3¾	95				
MXPRSK11 4½ 114	MXPRSK09	4	102				
110.11.11.11.11	MXPRSK10	41/4	108				
MYPRSK12 //3/ 121	MXPRSK11	41/2	114				
IVIALITATION 12	MXPRSK12	43/4	121				
MXPRSK13 5 127	MXPRSK13	5	127				

Donut

Support for 3rd degree prolapse/ procidentia. Not to be deflated during insertion or removal. See Inflatoball if a deflatable pessary is necessary.

Product # MXPDO (Pessary Only) Product # MXKPDO (Pessary, Trimo-San, Jel Jector)



Donut					
Product #					
MXPD000	2	51			
MXPD001	21/4	57			
MXPD002	21/2	64			
MXPD003	23/4	70			
MXPD004	3	76			
MXPD005	31/4	83			
MXPD006	31/2	89			
MXPD007	3¾	95			

Shaatz

For mild prolapse complicated by a mild cystocele.

Product # MXPSH (Pessary Only) Product # MXKPSH (Pessary, Trimo-San, Jel Jector)



Shaatz					
MXKPSH1-1/2	1 1/2	38			
MXKPSH1-3/4	13/4	44			
MXKPSH2-	2	51			
MXKPSH2-1/4	21/4	57			
MXKPSH2-1/2	21/2	64			
MXKPSH2-3/4	23/4	70			
MXKPSH3-	3	76			
MXKPSH3-1/4	31/4	83			
MXKPSH3-1/2	31/2	89			
MXKPSH3-3/4	3¾	95			

Regula

For uterine prolapse (1st and 2nd degree).

Product # MXPREG (Pessary Only) Product # MXKPREG (Pessary, Trimo-San, Jel Jector)



Regula Folding						
Product #	in. (A)	mm (A)	in. (B)		in. (C)	mm (C)
MXPREG00	1	25	1 15/16	49	1%	35
MXPREG01	1 1/16	27	2	51	11/2	38
MXPREG02	11/8	29	21/8	54	15/8	41
MXPREG03	11/4	32	23/16	56	13/4	44
MXPREG04	13/8	35	21/4	57	17/8	48
MXPREG05	17/16	37	2%	60	2	51
MXPREG06	1½	38	25/8	67	21/8	54
MXPREG07	1%	41	21/8	73	21/4	57
MXPREG08	13/4	44	31/8	79	2%	60

Gellhorn

For 3rd degree prolapse/procidentia.

Product # MXPGE (Pessary Only) Product # MXKPGE (Pessary, Trimo-San, Jel Jector)

Also available in: Short Stem Product # MXPGSS Rigid Product # MXPGRS (Pessary Only)

Short Stem Product # MXKPGSS Rigid Product # MXKPGRS (Pessary, Trimo-San, Jel Jector)



Gellhorn					
Flexible	Rigid	Short Stem			
Product #	Product #	Product #		mm	
MXKPGE1-1/2	MXKPGRS1-1/2	MXKPGSS1-1/2	1 1/2	38	
MXKPGE1-3/4	MXKPGRS1-3/4	MXKPGSS1-3/4	13/4	44	
MXKPGE2-	MXKPGRS2-	MXKPGSS2-	2	51	
MXKPGE2-1/4	MXKPGRS2-1/4	MXKPGSS2-1/4	21/4	57	
MXKPGE2-1/2	MXKPGRS2-1/2	MXKPGSS2-1/2	21/2	64	
MXKPGE2-3/4	MXKPGRS2-3/4	MXKPGSS2-3/4	23/4	70	
MXKPGE3-	MXKPGRS3-	MXKPGSS3-	3	76	
MXKPGE3-1/4	MXKPGRS3-1/4	MXKPGSS3-1/4	31/4	83	
MXKPGE3-1/2	MXKPGRS3-1/2	MXKPGSS3-1/2	31/2	89	
MXKPGE3-3/4	MXKPGRS3-3/4	MXKPGSS3-3/4	3¾	95	

Milex® Pessaries for Prolapse and Incontinence

P = Prolapse I = Incontinence

Hodge

For stress urinary incontinence for a woman who has a very small vaginal introitus.

Product # MXPEH (Pessary Only) Product # MXKPEH (Pessary, Trimo-San, Jel Jector)



Hodge w/ Support

For stress urinary incontinence with a mild cystocele for a woman who has a very small vaginal introitus.

Product # MXPHS (Pessary Only) Product # MXKPHS (Pessary, Trimo-San, Jel Jector)



Hodge w/ Knob

For stress urinary incontinence for a woman who has a small vaginal introitus.

Product # MXPEHK (Pessary Only) Product # MXKPEHK (Pessary, Trimo-San, Jel Jector)



Hodge w/ Support & Knob

For stress urinary incontinence complicated by a mild cystocele for a woman who has a small vaginal introitus.

Product # MXPHSK (Pessary Only) Product # MXKPHSK (Pessary, Trimo-San, Jel Jector)



Hodge							
Hodge		w/ Knob	w/ Support & Knob				
Product #	Product #	Product #	Product #				
MXPEH00	MXPHS00	MXPEHK00	MXPHSK00	23/4	70	13/4	44
MXPEH01	MXPHS01	MXPEHK01	MXPHSK01	3	76	11//8	48
MXPEH02	MXPHS02	MXPEHK02	MXPHSK02	31/4	83	115/16	49
MXPEH03	MXPHS03	MXPEHK03	MXPHSK03	31/2	89	2	51
MXPEH04	MXPHS04	MXPEHK04	MXPHSK04	35/8	92	21/8	54
MXPEH05	MXPHS05	MXPEHK05	MXPHSK05	33/4	95	21/4	57
MXPEH06	MXPHS06	MXPEHK06	MXPHSK06	31/8	98	2%	60
MXPEH07	MXPHS07	MXPEHK07	MXPHSK07	41/4	108	21/2	64
MXPEH08	MXPHS08	MXPEHK08	MXPHSK08	45/8	117	25/8	67
MXPEH09	MXPHS09	MXPEHK09	MXPHSK09	5	127	23/4	70

Incontinence Dish

For stress urinary incontinence with mild prolapse.

Product # MXPCOND (Pessary Only) Product # MXKPCOND (Pessary, Trimo-San, Jel Jector)



Incontinence Dish with Support

For stress urinary incontinence with mild prolapse and a mild cystocele.

Product # MXPCONDS (Pessary Only) Product # MXKPCONDS (Pessary, Trimo-San, Jel Jector)



Incontinence Dish					
Dish	Dish w/ Support				
Product #	Product #				
MXPCOND00	MXPCONDS00	21/8	55		
MXPCOND01	MXPCONDS01	2%	60		
MXPCOND02	MXPCONDS02	21/2	65		
MXPCOND03	MXPCONDS03	23/4	70		
MXPCOND04	MXPCONDS04	3	75		
MXPCOND05	MXPCONDS05	31/8	80		
MXPCOND06	MXPCONDS06	3%	85		

Risser

Designed for patient with shallow pubic notch. Has larger weight bearing zone to support vaginal aspects of pubis. Helpful in cases of repeated miscarriages due to an incompetent cervix by repositioning the weight of the growing fetus.



Product # MXPRI (Pessary Only) Product # MXKPRI (Pessary, Trimo-San, Jel Jector)

Risser							
Product #	in. (A)	mm (A)	in. (B)				
MXPRI00	31/8	79	2	51			
MXPRI01	31/4	83	21/8	54			
MXPRI02	31/2	89	21/4	57			
MXPRI03	33/4	95	23/8	60			
MXPRI04	41/4	108	21/2	64			
MXPRI05	41/2	114	25/8	67			
MXPRI06	43/4	121	23/4	70			
MXPRI07	5	127	21/8	73			
MXPRI08	5½	140	3	76			
MXPRI09	53/4	146	31/8	79			

Gehrung with Knob/Folding

For cystocele with or without mild rectocele complicated by stress urinary incontinence.

Product # MXPGSK (Pessary Only) Product # MXKPGSK (Pessary, Trimo-San, Jel Jector)




Gehrung/Folding

Support for cystocele and rectocele. It is also more adaptable in cases of procidentia where the uterus tends to herniate when other pessaries are used.

Product # MXPGS (Pessary Only) Product # MXKPGS (Pessary, Trimo-San, Jel Jector)





Smith

Modification of Hodge pessary. For stress urinary incontinence designed for patients with a well-defined pubic-notch. Helpful in cases of repeated miscarriages due to incompetent cervix by repositioning the weight of the growing fetus.





Smith						
Product #	in. (A)	mm (A)	in. (B)	mm (B)		
MXPES00	31/8	79	2	51		
MXPES01	31/4	83	21/8	54		
MXPES02	31/2	89	21/4	57		
MXPES03	3¾	95	2%	60		
MXPES04	41/4	108	21/2	64		
MXPES05	41/2	114	25/8	67		
MXPES06	43/4	121	23/4	70		
MXPES07	5	127	21/8	73		
MXPES08	5½	140	3	76		
MXPES09	5¾	146	31/8	79		

Cube

Support for 3rd degree prolapse/procidentia. Should be removed nightly.

Product # MXPEC (Pessary Only) Product # MXKPEC (Pessary, Trimo-San, Jel Jector)

Also available in: Cube with Drainage Holes Product # MXPECH (Pessary Only) Product # MXKPECH (Pessary, Trimo-San, Jel Jector)



Cube			
Cube	w/ Holes		
Product #	Product #	in.	mm
MXPEC00	MXPECH00	1	25
MXPEC01	MXPECH01	13/16	30
MXPEC02	MXPECH02	1%	35
MXPEC03	MXPECH03	1½	38
MXPEC04	MXPECH04	15/8	41
MXPEC05	MXPECH05	13/4	44
MXPEC06	MXPECH06	2	51
MXPEC07	MXPECH07	21/4	57

Tandem-Cube

Support for 3rd degree prolapse/ procidentia when single cube provides inadequate support. Should be removed nightly.

Product # MXPECT (Pessary Only) Product # MXKPECT (Pessary, Trimo-San, Jel Jector)





Tandem-Cube			
Product #	in.	mm	
MXPECT2-0	1% & 1	35 & 25	
MXPECT3-1	11/2 & 13/16	38 & 30	
MXPECT4-2	15/8 & 13/8	41 & 35	
MXPECT5-3	13/4 & 11/2	44 & 38	
MXPECT6-4	2 & 15/8	51 & 41	
MXPECT7-5	21/4 & 13/4	57 & 44	

Obstetrical Care

Improving Clinical Outcomes Is Our Mutual Goal

CooperSurgical offers a growing line of products for obstetrics. Select from patient monitors, fetal Dopplers and a full range of market-leading vacuum-assisted delivery products from our Labor and Delivery catalog.

NOW AVAILABLE THROUGH DISTRIBUTION

F9 Fetal Monitor	6	3
ClearTone® Fetal	Dopplers6	4



F9 Fetal Monitor with Deluxe Car

F9 Fetal Monitor

Highly Accurate Technology Provides Greater Ultrasound Detection



CooperSurgical's F9 Fetal Monitor combines the latest ultrasound technology with lightweight, ergonomic design. It is capable of monitoring twins with its two 12-crystal, wide-beam 1 MHz waterproof transducers which provide more ultrasound detection capability. A 12-inch color touch screen ensures data are easily delineated on the screen, with clear color tracings that can be shared with the patient. The unit can be placed on a countertop or mounted to a mobile cart when in use, then folded flat for easy storage.

Other features include:

- Signal Overlap Verification to reduce errors caused by overlapping FHR signals when monitoring twins
- Customizable alarms and controls for every clinical environment
- Auto fetal movement detection and alert

For data capture and sharing, the F9 Fetal Monitor can be equipped with CooperSurgical Insight Software for PC-based data storage and review. The unit is EMR-compatible with this software, and contains a built-in printer and USB/ethernet connections. Its 24-hour CTG memory caching facilitates data storage and playback, even when the paper runs out.



F9	Fetal	Monitor	with	Standard	Accessor

F9 Fetal Monitor with Standard Accessories CA ◆ C€		
Number	Description	
MS9-108022-UL	CooperSurgical Fetal Monitor F9 / 110 v	
MS9-314012	CooperSurgical Fetal Monitor F9 / 220 v	
MS3-31528-A1	Transducer 1 (pink)	
MS3-107705-A1	Transducer 2 (yellow)	
MS3-31527-A1	Toco Transducer (blue)	
MS3-31112	Event Marker	
MS1R-02264	Transducer Belt (reusable)	
M25R-75111	Recording Paper (90 mm \times 152 mm \times 150 p US standard) (Pack of 3)	

Optional Accessories		CA * (€
Number	Description	
MS9-108022-SWP	CooperSurgical Insight Software Package (includes software manual and cables)	
M21-64072	Rechargeable Lithium-Ion Battery	
MS9-107634	Deluxe Cart	
MS9-108166	Base Cart with Basket and Laptop Tray	
MS9-108108	Base Cart with Basket	



The ClearTone is your first-class choice when selecting a new Fetal Doppler. With advanced audio output that provides superior sound quality, the ClearTone has interchangeable fetal probes and is available in a waterproof model for all clinical settings. The durable and elegant Doppler has an optional easy-to-read digital display that calculates beats per minute to save counting and calculating time. ClearTone reliably captures fetal heart rate from early through late gestation in a compact, robust unit.

Specifications:

Fetal Doppler: 2.1 MHz and 3.2 MHz Waterproof Fetal Doppler: 2.1 MHz

Weight: 11.3 oz. (320 g)

Outputs: One 3.5 mm jack to connect to a headset or

audio recording device

Controls: Separate volume control and push ON/OFF button

Power Source: Battery voltage – Three AA Alkaline 1.5-volt (non-rechargeable) or three AA NMHi 1.2-volt (rechargeable). Battery life – 1,250 one-minute exams with Alkaline or 1,000 one-minute exams with rechargeable batteries.

Dimensions: 5.5" x 2.75" x 1.3" (140 mm x 70 mm x 35 mm)

ClearTone	Fetal Dopplers	
WITH 3 MHz Pf	ROBES	CA * (€
Number	Description	
CT150-3	ClearTone Doppler, Non-Display, with 3 MHz Probe	
CT150R-3	ClearTone Doppler, Non-Display, Rechargeable with 3 N	MHz Probe
CT250-3	ClearTone Doppler, With Display, with 3 MHz Probe	
CT250R-3	ClearTone Doppler, With Display, Rechargeable with 3 N	1Hz Probe
WITH 2 MHz PI	ROBES	CA* (€
Number	Description	
CT150-2	ClearTone Doppler, Non-Display, with 2 MHz Probe	
CT150R-2	ClearTone Doppler, Non-Display, Rechargeable with 2 N	/IHz Probe
CT250-2	ClearTone Doppler, With Display, with 2 MHz Probe	
CT250R-2	ClearTone Doppler, With Display, Rechargeable with 2 M	ЛНz Probe
WITH 2 MHz PI	ROBES – WATERPROOF	CA* (€
Number	Description	
CT150-2WP	ClearTone Doppler, Non-Display, with 2 MHz Probe	
CT150R-2WP	ClearTone Doppler, Non-Display, Rechargeable with 2 M	Hz Probe
CT250-2WP	ClearTone Doppler, With Display, with 2 MHz Probe	
CT250R-2WP	ClearTone Doppler, With Display, Rechargeable with 2 M	Hz Probe

ClearTone Probes – Fetal and Vascular

CA***** (€

Interchangeable fetal and vascular probes and optional waterproof fetal model provide the versatility you need in all clinical settings. Probes are available in several ultrasonic frequencies to suit your needs.

Number	Description
CTP-2	ClearTone 2 MHz Fetal Probe
CTP-W	ClearTone 2 MHz Waterproof Fetal Probe
CTP-3	ClearTone 3 MHz Fetal Probe
CTP-4	ClearTone 4 MHz Vascular Probe
CTP-5	ClearTone 5 MHz Vascular Probe
CTP-8	ClearTone 8 MHz Vascular Probe



Vascular Dopplers also available:

ClearTone Vascular Dopplers		
WITH 5 MH	z PROBES CA★ C€	
Number	Description	
CT150-5	ClearTone Doppler, Non-Display, with 5 MHz Probe	
CT150R-5	ClearTone Doppler, Non-Display, Rechargeable with 5 MHz Probe	
CT250-5	ClearTone Doppler, With Display, with 5 MHz Probe	
CT250R-5	ClearTone Doppler, With Display, Rechargeable with 5 MHz Probe	
WITH 8 MH	z PROBES CA★ C€	
Number	Description	
CT150-8	ClearTone Doppler, Non-Display, with 8 MHz Probe	
CT150R-8	ClearTone Doppler, Non-Display, Rechargeable with 8 MHz Probe	
CT250-8	ClearTone Doppler, With Display, with 8 MHz Probe	
CT250R-8	ClearTone Doppler, With Display, Rechargeable with 8 MHz Probe	



Meeting All Your Patients' Needs

CooperSurgical offers a wide array of safe, effective family planning methods including contraception, fertility testing and insemination.



Intrauterine Insemination



Simplify Evaluation, Preparation and Insemination

The Sperm Select family of products is designed to simplify the steps required for sample evaluation, preparation and insemination. Each kit contains all the components necessary to optimize sperm preparation in your office. Sperm Select separates samples in 40 minutes or less, with progressively motile sperm at the top and seminal debris at the bottom. No expensive equipment is necessary and the process is easy to learn and use. For organization and consistent preparation, consider adding the cost-effective Blockheater™ desktop incubator for sample warming.

Number	Description CA*
SM600	Sperm Select System, Syringe, Sperm Medium, Syringe & Vial Holder
SM602-F	Sperm Select System Kit – Complete Sperm Select System with (1) Select IUI® Catheter
SM601-F	Sperm Select System Kit — Complete Sperm Select System with (1) Minispace® IUI Catheter
SM603-F	Sperm Select System Kit – Complete Sperm Select System with (1) Select IU® Catheter
SM700	Sperm Select Blockheater (Incubator)

Kit Contents:

- Luer Lock glass syringe of Sperm Select sterile, membrane-filtered, MEA-assay Sodium Hyaluronate
- Vial of medium (0.5 ml) sterile, membrane-filtered, MEA-assay
- Syringe carrier
- Vial holder
- Name tag
- Directions for Use
- (1) Catheter of your choice



Sperm Select DVD Available Upon Request

Insemination Catheters



TheCurve®	CA # (€

Slight Curve Aids Insertion

TheCurve Intrauterine Insemination Catheter with Luer Lock hub has a rounded, beveled tip which helps to ensure smooth insertion, while the slight curve is designed to facilitate entry. For atraumatic sample delivery.

Number	Description
SM507	17.2 cm x 1.6 mm, .02cc Single-opening (25/Box)

Minispace[®] CA***** C€

The Smallest IUI Catheter

Minispace's small diameter helps minimize uterine cramping and maximize sperm use. This ultra-low-volume catheter leaves minimal residual after delivery. Luer Lock hub.

Number	Description
SM610	17.2 cm x 1.6 mm, .07 cc Single-opening (25/Box)

Select IUI[®] CA**◆** C€

An IUI Catheter for Every Patient

The gold standard of insemination, Select IUI catheters have dual openings that mirror the normal orientation of the fallopian tubes. The wings of the Luer Lock hub indicate the position of the catheter openings. Sterile/single-use.

Number	Description
SM400	17 cm x 2.1 mm, .34 cc Dual-opening (25/Box)



The Milex fine or extra-fine insemination cannula will help maximize your ICI or IUI success rate. Other cannulae/catheters marketed for insemination are longer and larger in lumen size and more apt to stimulate the uterus to contract and expel the semen specimen. The shorter length also minimizes the amount of sperm required for optimal delivery.

- Open end (8 French)
- Closed end, 2 side ports (3.5 French)
- Sterile/Single-use

Number	Description
MX730	Milex Insemination Cannula with Syringe, 8 Fr, Open End (24/Box)
MX732	Milex Insemination Cannula with Syringe, 3.5 Fr, Closed End (24/Box)

Uni-Sem[™] Intrauterine Insemination Catheter





For Smooth, Atraumatic Delivery

Uni-Sem's step-down outer diameter design and unique construction aid in smooth cervical passage for atraumatic placement. Consistently lumened throughout for gentle media transfer, the distal portion is slightly smaller at 2.3 mm O.D. for more atraumatic placement while the remainder of the catheter is 2.6 mm O.D. for added clinician control during insertion. Accurate placement is enhanced through use of the orientation wings on the hub that align with the catheter's ejection ports. As a further aid to placement, depth markings at 5.5 cm and 6.5 cm are indicated on the shaft. The clear tip allows visual confirmation of successful delivery. Uni-Sem is gamma-sterilized, single-use.

Number	Description
7033	Uni-Sem Intrauterine Insemination Catheter (25/Box)

Fertility Assessment

SelectMucus® Endocervical Aspirator CA◆




Designed to sample the contents of the endocervical canal for postcoital testing and any situation where the mucus plug needs analysis or removal. Smooth, rounded tip minimizes patient discomfort.

Number	Description
SM301	SelectMucus Endocervical Sampler (25/Box)

Pro-Ception Oligospermia Cups

CA4



Anatomically shaped to fit the cervix, these hypoallergenic cups provide a simple yet efficient means of ensuring the maximum concentration of semen at the cervical Os. Increased pooling of the seminal fluid during homologous or heterogeneous insemination helps maximize the potential for success in these procedures. Cups are available with and without a stem and in eight sizes: 24 mm, 26 mm, 28 mm, 30 mm, 32 mm, 34 mm, 36 mm and 38 mm.

Number	Description
MXFOCS24	24 mm Pro-Ception Oligospermia Cup with Stem (5/Box)
MXFOCS26	26 mm Pro-Ception Oligospermia Cup with Stem (5/Box)
MXFOCS28	28 mm Pro-Ception Oligospermia Cup with Stem (5/Box)
MXFOCS30	30 mm Pro-Ception Oligospermia Cup with Stem (5/Box)
MXFOCS32	32 mm Pro-Ception Oligospermia Cup with Stem (5/Box)
MXFOCS34	34 mm Pro-Ception Oligospermia Cup with Stem (5/Box)
MXFOCS36	36 mm Pro-Ception Oligospermia Cup with Stem (5/Box)
MXFOCS38	38 mm Pro-Ception Oligospermia Cup with Stem (5/Box)
MXFOC24	24 mm Pro-Ception Oligospermia Cup without Stem (5/Box)
MXFOC26	26 mm Pro-Ception Oligospermia Cup without Stem (5/Box)
MXFOC28	28 mm Pro-Ception Oligospermia Cup without Stem (5/Box)
MXFOC30	30 mm Pro-Ception Oligospermia Cup without Stem (5/Box)
MXFOC32	32 mm Pro-Ception Oligospermia Cup without Stem (5/Box)
MXFOC34	34 mm Pro-Ception Oligospermia Cup without Stem (5/Box)
MXFOC36	36 mm Pro-Ception Oligospermia Cup without Stem (5/Box)
MXF0C38	38 mm Pro-Ception Oligospermia Cup without Stem (5/Box)



This cannula with Vacu-Lok Control Syringe is designed for the following procedures:

- Postcoital testing (Huhner)
- Removal of mucus plug
- Testing interaction of sperm and cervical mucus viscosity
- Fractional sampling of cervical mucus
- Determination of ovulation
- Identification of antisperm antibodies
- Microscopic determination of colonizing microorganisms

Products are sterile/single-use and individually packaged.

Number	Description
MX710	Pro-Ception Fertility Cannula, Individual, No Syringes (36/Box)
MX720	Pro-Ception Fertility Cannula, Individual, with Vacu-Lok Control Syringe (12/Box)

Pro-Ception Fertility Pak

CA*



Patients can easily insert these disposable vaginal tampons themselves to maintain a pool of seminal fluid near the cervical Os and help improve fertilization success rates. They can also be used following IUI or ICI procedures in your office.

Single-use/Disposable.

Number	Description
MXFEPAK	Pro-Ception Fertility Pak — Small (12/Box)
MXFEPAKL	Pro-Ception Fertility Pak — Large (12/Box)

Fertility Assessment

Aspirette® Endocervical Aspirator





Cervical Sampling Made Easy

The Aspirette Endocervical Aspirator is completely self-contained, sterile and ready for use. With a slender (3.1 mm O.D.) body and rounded distal tip, cervical insertion is easy with minimal risk, trauma or patient discomfort. Clear depth markings aid catheter placement. Once collected, the sample is safe from cross-contamination and ready for clinical assessment.

Number	Description
8700	Aspirette Endocervical Aspirator (25/Box)

Diaphragms



The Milex Wide-Seal Arcing Style Diaphragm has a tension-adjusted spring that curves in one place, providing increased suction for added protection. The patented Milex Wide-Seal Diaphragms are made of silicone and completely latex-free. Must be used with spermicidal gel.

MILEX WIDE-SEAL ARCING STYLE DIAPHRAGM		CA* (€
Number	Description	
MXWS60	Arcing Diaphragm 60 mm	
MXWS65	Arcing Diaphragm 65 mm	
MXWS70	Arcing Diaphragm 70 mm	
MXWS75	Arcing Diaphragm 75 mm	
MXWS80	Arcing Diaphragm 80 mm	
MXWS85	Arcing Diaphragm 85 mm	
MXWS90	Arcing Diaphragm 90 mm	
MXWS95	Arcing Diaphragm 95 mm	

ACCESSORIES		CA* (€
Number	Description	
MXWSFS	Arcing Diaphragm Fitting Set	
MXWFFS	Omniflex Diaphragm Fitting Set	
MXFIT0072	Pelvic Teaching Model	

The Milex Wide-Seal Omniflex Style Diaphragm has a distortion-free spring that provides arc no matter where the rim is compressed, providing increased suction for added protection. The patented Milex Wide-Seal Diaphragms are made of silicone and completely latex-free. Must be used with spermicidal gel.

MILEX WIDE-SEAL OMNIFLEX STYLE DIAPHRAGM		CA* (€
Number	Description	
MXWF60	Omniflex Diaphragm 60 mm	
MXWF65	Omniflex Diaphragm 65 mm	
MXWF70	Omniflex Diaphragm 70 mm	
MXWF75	Omniflex Diaphragm 75 mm	
MXWF80	Omniflex Diaphragm 80 mm	
MXWF85	Omniflex Diaphragm 85 mm	
MXWF90	Omniflex Diaphragm 90 mm	
MXWF95	Omniflex Diaphragm 95 mm	

In-Office Tubal Patency and Uterine Structure Evaluation

ABBI®

Air Bubble Based Infuser

- In-office tubal patency evaluation with Sono HSG
- Contrast delivery with a continuous stream of bubbles
- Uses existing ultrasound and a standard intrauterine catheter
- Avoids patient exposure to radiation and X-ray dyes
- Capability of performing an SIS during the same patient visit


Safely Evaluate Tubal Patency Using Ultrasound Versus X-Ray


ABBI (Air Bubble Based Infuser) utilizes ultrasound with a saline and air bubble contrast media to safely evaluate tubal patency. The echogenic air bubbles produced by ABBI will flow continuously through the cornua, fallopian tubes and fimbriae to confirm patency, or pool in the event of tubal occlusion. ABBI provides a minimally invasive procedure by reducing the risk of allergic reaction to radiopaque dyes, and eliminating exposure to fluoroscopy radiation used in traditional HSG.

Dual-Function ABBI also Enables Saline-Only SIS Exam

Additionally, ABBI can be used to perform a saline-only SIS examination to identify abnormalities in the uterine structure including polyps and fibroids that may be interfering with pregnancy.

ABBI® Air Bubble Based Infuser CA	
Number	Description
ABBI	ABBI (Air Bubble Based Infuser) — Single-Use Disposable (1/Box)
ABBI-Kit	ABBI (Air Bubble Based Infuser) with H/S Elliptosphere Catheter (1/Box)





Clinical Support for the Use of Saline and Air For Tubal Patency (HyCoSy/Sono HSG) and Uterine Structure Evaluation


Hysterosalpingo contrast sonography has been introduced in clinical practice as an effective tool for tubal patency and uterine cavity evaluation. HyCoSy shows high overall accuracy in the evaluation of both tubal patency and uterine cavity morphology. Furthermore, HyCoSy versus fluoroscopy avoids both exposure to ionizing radiation and injection of iodinated contrast medium.¹

HyCoSy can be performed as an office procedure. Infusing saline directly into the uterus as a homogeneous, echo-free contrast medium enables good visualization of the uterine architecture. Tubal patency can be assessed by using an air saline mixture creating contrast visible along the fallopian tubes and near the ovarian fossas. In fact, compared with RX-HSG, it is equally accurate, but less painful and with no risk of allergic reactions to the contrast medium and with no patient exposure to radiation.²

This procedure can be performed in the gynecologist's office and enables accurate and complete evaluation of the uterus, uterine cavity and fallopian tubes. The most simple and inexpensive contrast medium used is saline solution mixed with air. Tubal patency is observed by visualizing the hyperechoic air bubbles traversing the tubal course.³



² Graziano A, Lo Monte G, Soave I, Caserta D, Moscarini M, Marci R, et al: Sonohysterosalpingography: a suitable choice in infertility workup. Journal of Medical Ultrasonics July 2013, Volume 40, Issue 3, pp 225-2292




³ Danielle E. Luciano, MD, Caterina Exacoustos, MD, and Anthony Luciano, MD Contrast Ultrasonography for Tubal Patency, Journal of Minimally Invasive Gynecology November/December 2014

General Examination: Diagnostic Disposables

H/S Elliptosphere[™] Catheter Set and Tray

The Safe and Simple Answer for Both Sonohysterography and Hysterosalpingography






The Elliptosphere balloon contour and variable size permit secure placement in the cervical canal or uterine cavity for optimal visualization. The balloon's unique design prevents retrograde leakage of distension media, which facilitates peak uterine distention. A small catheter size enables placement next to ultrasound probe for minimal patient discomfort. The H/S Elliptosphere Catheter can be used with oil- or saline-based medium. Available in two versions — with and without an integrated stylet — the H/S Elliptosphere features latex-free construction.


Number	Description	
61-4005	H/S Elliptosphere Catheter, 5 Fr (10/Box)	CA* (€
61-4005S	H/S Elliptosphere Catheter Set with Integrated Stylet (10/Box)	CA* (€
61-4205	H/S Elliptosphere Tray with 5 Fr Catheter (10/Case)	


H/S Catheters


For Enhanced Visualization of Uterine Pathology











Pre-Saline Infusion

Post-Saline Infusion

H/S Catheter Set CA♦ C€

These easy-to-use devices provide enhanced visualization. Designed for use in sonohysterography and hysterosalpingography, the catheter can be used with oil- or saline-based media depending on the procedure requirements.

Features and Benefits:

- Available with integrated stylet for easy insertion
- Unique placement sheath eliminates the need for a tenaculum
- Comes in both 5 and 7 French
- Ideal for office use
- Latex-free construction

Number	Description
61-5005	H/S Catheter, Kraton, 5 Fr (10/Box)
61-5007	H/S Catheter, Kraton, 7 Fr (10/Box)
61-0020	H/S Stylet Catheter Stiffener (10/Box)

H/S Procedure Tray





- Fenestrated drape
- Disposable speculum
- Disposable cervical dilator
- Three swab sticks
- Sterile wrap
- One pack of povidone-iodine solution
- 36-in. extension tube
- 20 cc syringe with 18-gauge needle
- Lubricating jelly packet
- Sanitary napkin
- Inner tray for povidone-iodine solution

Number	Description
61-5205	H/S Procedure Tray — Round 5 Fr
61-5207	H/S Procedure Tray — Round 7 Fr
61-3205	H/S Procedure Tray — Round 5 Fr with Stylet
61-3207	H/S Procedure Tray — Round 7 Fr with Stylet


Catheters

Tampa Infusion Catheter

For sonohysterography – the simple, cost-effective choice for ultrasound examination of the uterus.

- Unique introducer sheath (33 cm long) facilitates insertion into cervical Os with minimal trauma
- 7 cm marker provides optimal placement of catheter tip in uterus
- Balloonless configuration is ideal for postmenopausal patients with stenotic cervical canals, as well as for many nulliparous patients
- Small catheter size allows easy placement next to ultrasonography probe

Number	Description
61-2005	Tampa Catheter (20/Box)



HSG Procedure Tray

The HSG Procedure Tray is a sterile, single-use product ideal for use in hysterosalpingography and sonohysterography procedures. These convenient trays reduce time spent gathering products because they contain the essential items needed during a procedure.

Number	Description
6050T	Procedure Tray



Breast Aspiration Dual-Port Biopsy Needle



Breast aspiration with a dual-port needle is considered highly accurate for the cytological evaluation of breast cysts and nodules. The aspiration is less traumatic and is well tolerated. Reduces the need for and cost of hospitalization. Procedure can be performed in your office in 15 minutes or less.

Number	Description
MX40	Dual-Port Needle with 3-finger control with 12 cc Syringe, 3" (12/Box)
MX41	Dual-Port Needle with 3-finger control with 12 cc Syringe, 1½" (12/Box)
MX42	Dual-Port Needle, No Syringe, 1½" (5/Box)

Endocervical Curette with Vacu-Lok Syringe



This nylon curette has a sharp Randall-type cutting edge. The outside diameter tapers to ensure adequate endocervical sampling, while the Vacu-Lok syringe allows adjustment of suction from approximately 15- to 20-inch Hg negative pressure. This sterile, single-use disposable is latex-free.

Number	Description
MX111	Endocervical Curette with Vacu-Lok Syringe (25/Box)

Word Catheter



The Word Catheter is an effective office treatment following an incision and drainage of a Bartholin cyst or abscess, providing your patient immediate relief. Bartholin cysts/abscesses have a high rate of recurrence unless an accessory duct is provided to permanently drain the obstructed acini.

Number	Description
53404	Word Bartholin Gland Catheter. Includes syringe and scalpel. Latex-free. (6/Box)

Milex® Cervical Spatulas



Anatomically designed for maximum cell pickup of squamocolumnar junction cells, endocervical cells and ectocervical cells. Blade is mildly serrated for slight cervical scraping and textured on one side for better cell pickup. Sterile and disposable.

Number	Description
MX600	Milex Cervical Spatulas, Small, Yellow (250/Box)
MX610	Milex Cervical Spatulas, Medium, White (250/Box)

Disposa-Derm®

For Fast and Simple Skin Biopsy

Disposa-Derm offers a simple and easy way to biopsy the skin on any exposed area. Deep ribbing and a contoured handle provide additional control during biopsy. These cost-effective punches come presterilized in individual pouches. The punches are color-coded for easy size identification.



GREEN	
Number	Description
DD-2 mm	Disposa-Derm 2 mm (5/Box)



PINK	
DD-5 mm	Disposa-Derm 5 mm (5/Box)



CA*

BLUE	
Number	Description
DD-3 mm	Disposa-Derm 3 mm (5/Box)



ORANGE	
DD-6 mm	Disposa-Derm 6 mm (5/Box)



YFLLOW	
	B
Number	Description

5-PC Sample Pack (1 of each size)	DISPOSA-DERM COMBO PACK		
	5-PC	Sample Pack (1 of each size)	

Dilators

CA*

Made of a malleable Teflon-like material for frictionless natural lubrication, the Comfort-Flex Dilator Set is excellent for use in intracervical dilation, uterine sounding or intrauterine procedures. The gradual dilation, combined with the flexible tip, ensure comfortable introduction and gentle dilation. Set Includes:

- Os Locator (1 ring) naturally seeks and gently opens the external Os
- Canal Dilator (2 rings) gently traces the natural path of the cervical canal while steadily dilating
- Fundus Sound™ (3 rings) opens the internal Os and sounds the



Milex® Vaginal-Hymenal Silicone Dilators

CA*

Vaginal dilators can be used to overcome the pain of sexual intercourse experienced by some patients through muscle-relaxation techniques. The step-by-step instructions describe how to use vaginal dilators to gradually increase patient comfort level. Vaginal dilators are also recommended following reconstructive surgery and radiation therapy to help keep the vaginal canal supple.

Number	Description
MX20	Milex Vaginal Dilator – Set of 4: Small, Medium, Large, Extra-Large with Trimo-San (USA sales only)
MX20-CA	Milex Vaginal Dilator – Set of 4: Small, Medium, Large, Extra-Large without Trimo-San
MX20I	Milex Vaginal Dilator Single — Small 51/8" x 15/16"
MX20IM	Milex Vaginal Dilator Single – Medium 5½" x 1½6"
MX20IL	Milex Vaginal Dilator Single – Large 5¾" x 1¼"
MX20IXL	Milex Vaginal Dilator Single — Extra-Large 6%" x 11/16"



Milex® Vaginal-Rectal Silicone Dilators

CA*

Vaginal-Rectal Dilators are recommended following reconstructive surgery to keep the canal soft and supple. The step-by-step instructions describe how to gradually increase acceptance.

Number	Description
MX25	Milex Rectal Dilator — Set of 4: Small Rectal, Medium Rectal, Large Rectal, Extra-Large Vaginal-Rectal with Trimo-San (USA sales only)
MX25-CA	Milex Rectal Dilator – Set of 4: Small Rectal, Medium Rectal, Large Rectal, Extra-Large Vaginal-Rectal without Trimo-San
MX25IS	Milex Vaginal-Rectal Dilator Single — Small 3" x 11/16"
MX25IM	Milex Rectal Dilator Single – Medium 3½" x ½6"
MX25IL	Milex Rectal Dilator Single – Large 3%6" x 11/6"
MX25IXL	Milex Vaginal-Rectal Dilator Single — Extra-Large 31% " x 15%"



Milex® Vaginal Silicone Dilators

CA*

Vaginal dilators provide gentle, atraumatic dilation of vaginal muscles. A step-by-step guide is included to aid the patient.

Number	Description
MX23	Milex Dilator – Set of 4: Small, Medium, Large, Extra-Large with Trimo-San (USA sales only)
MX23-CA	Milex Dilator – Set of 4: Small, Medium, Large, Extra-Large without Trimo-San
MX23I	Milex Dilator Single – Small 21/16" x ½" x ½"
MX23IM	Milex Dilator Single – Medium 2¾ε" x ¾ε" x ¾ε"
MX23IL	Milex Dilator Single − Large 2¹⅓₅" x ½"
MX23IXL	Milex Dilator Single — Extra-Large 3" x 11/6" x 1/6"



CA*

Cervical and Endocervical Samplers

Medscand® Cytobrush Plus® CA◆ C€



The Classic Cytobrush Plus Cell Collector for endocervical sampling is excellent for nulliparous or menopausal women, or for collecting endocervical cells at conization. The gentle brush collects large amounts of cells that are easily transferred to slides for analysis.

Number	Description
C0021	Cytobrush Plus – 2 Bags, 500 Brushes each (1,000/Box)
C0004	Cytobrush Plus — 10 Bags, 10 Brushes each (100/Box)
C0005	Cytobrush Plus – 25 Bags, 100 Brushes each (2,500/Box)
C0012	Cytobrush Plus Sterile, 1 Brush per Pouch (40/Box)

Medscand® Cytobrush Plus® GT CA★ C€

The Cytobrush Plus GT gentle-touch tip provides all the benefits of the Classic brush plus a protective bulb designed to minimize trauma to the canal, thus yielding a high-quality tissue sample with minimal contamination.

Number	Description
C0121	Cytobrush Plus GT – 2 Bags, 500 Brushes each (1,000/Box)
C0104	Cytobrush Plus GT – 10 Bags, 10 Brushes each (100/Box)
C0105	Cytobrush Plus GT – 25 Bags, 100 Brushes each (2,500/Box)
C0112	Cytobrush Plus GT Sterile, 1 Brush per Pouch (40/Box)
C0305	Cytobrush Plus GT Scored – 25 Bags, 100 Brushes each (2,500/Box)

Medscand® Uterobrush®

Single-use brush to collect all samples from the endometrium.

Number	Description
E0001	Uterobrush, Sterile (500/Box)



Medscand® Pap-Perfect® Spatula

Use the Medscand Pap-Perfect Plastic Spatula in combination with Cytobrush Plus or Cytobrush Plus GT. The spatula gently rubs the mucosa, ensuring a rich, complete sample from the ectocervix. Made of inert polypropylene plastic, the spatula will release collected cells more easily than a wooden spatula.

Number	Description	
1108	Pap-Perfect Endocervical Scrapers (Wooden) (500/Box)	
11080	Pap-Perfect Plastic Spatulas (500/Box)	CE



Medscand® Cytobrush Plus® Kits

An economical solution for one of the most frequently performed gynecology procedures. If you have a high-volume office or clinic, this may be your most effective solution for controlling material costs.



Number	Description	
10115-10/50	Cytobrush Plus Kits – Cell Collector (sets include brush, scraper, fixative, slide and slide transporter); 10 Trays with 50 Pap Smear test sets (500 total Tests/Case)	CA*
60115-10/50	Cytobrush Plus Cell Collector – Gentle Touch, 10 Trays with 50 Pap Smear test sets (500 total Tests/Case)	CA * (€
60115-20/50	Cytobrush Plus Cell Collector – Gentle Touch, 20 Trays with 50 Pap Smear test sets (1,000 total Tests/Case)	CA* (€

Index

Biopsy Punches Disposa-Derm Skin Biopsy Punches77 Euro-Med Classic Series
Cervical Samplers79
Colposcopy OptiK Model 1
Cryosurgery Cryo-Surg Handheld System 52 Frigitronics Cryo-Plus 53 Frigitronics CE-2000 54 LEISEGANG LM-900 Cryosurgery System 51 Ophthalmic 54-55
Cytobrush 79
Curettes-EndocervicalDisposable
Curettes-EndometrialDisposable36Stainless Steel28
Diaphragms Milex Wide-Seal Silicone Diaphragms 71
Dilateria39
Dilators Euro-Med Dilator Sets
DopplersDoppler Probes and Accessories
Electrodes: See LEEP
Endocervical Samplers Cytobrush79 Stainless Steel28
Endocervical Specula LEEP
Endometrial Samplers

Stainless Steel28
Endosee 8-9
Euro-Med Instruments (see also LEEP) Biopsy Punches
Needle Holders 34 Sounds and Dilators 27 Stainless Steel Curettes 28 Stainless Steel Endocervical Specula 30-34 Stainless Steel Forceps 30-34 Stainless Steel Hooks 29 Stainless Steel Lateral Wall Retractors 26 Stainless Steel Scissors 35 Stainless Steel Specula 24-25 Stainless Steel Tenacula 29
Fertility and Family Planning Diaphragms71 Fertility Assessment70-71 Intrauterine Insemination67-69
Fertility Assessment 71 Aspirette Endocervical Aspirator 71 H/S Catheter 75 H/S Elliptosphere Catheter 74 Iul Catheters 68 Pro-Ception Fertility Cannula 70 Pro-Ception Fertility Pak 70 Pro-Ception Oligospermia Cups 70 Sperm Select System 67
Fetal Dopplers – ClearTone 64-65
Fetal Monitors – F9 63
Forceps 50 LEEP 50 Stainless Steel 30-34
Frigitronics Cryo-Plus52-53
Frigitronics CE 200054-55
General Examination: Diagnostic Disposables Breast Aspiration Dual-Port Biopsy Needle
Catheters .74-76 Cervical and Endocervical Samplers .79 Dilators .78 Disposa-Derm .77
Endocervical Curette with Vacu-Lok Syringe77
Miley Conviced Spatules 77

Word Catheter......77

Gynecological Specialty Instruments

Dilators	27 38
Euro-Med Stainless Steel	
Endocervical Curettes Euro-Med Stainless Steel	28
Endocervical Specula	30
Euro-Med Stainless Steel Hooks	
Euro-Med Stainless Steel Specula 24- Stainless Steel Forceps 30-	
Lateral Wall Retractors	
Uterine Sounds	
H/S Elliptosphere Catheter	74
Hooks LEEP	
Stainless Steel	29
Hysteroscopy Endosee	3-9
Intrauterine Insemination	
IUI Catheters	
	0/
Lateral Wall Retractors	49
Stainless Steel	
LEEP	
Accessories46-	
Electrodes42-	44
	44 50
Electrodes 42- Endocervical Specula Fischer Cone Biopsy Excisor Forceps	44 50 42 50
Electrodes	44 50 42 50
Electrodes 42- Endocervical Specula Fischer Cone Biopsy Excisor Forceps Hooks Lateral Wall Retractors	44 50 42 50 50 49
Electrodes	44 50 42 50 50 49 41 45
Electrodes 42- Endocervical Specula	44 50 42 50 50 49 41 45
Electrodes 42- Endocervical Specula	44 50 42 50 50 49 41 45 49
Electrodes 42- Endocervical Specula	44 50 42 50 50 49 41 45 50 50
Electrodes 42- Endocervical Specula	44 50 42 50 50 49 41 45 49 50 47
Electrodes 42- Endocervical Specula	44 50 42 50 50 49 41 45 49 50 47 49
Electrodes 42- Endocervical Specula	44 50 42 50 50 49 41 45 49 50 47 49
Electrodes 42- Endocervical Specula	44 50 42 50 50 49 41 45 49 50 47 49
Electrodes	44 50 42 50 50 49 41 45 49 50 47 49
Electrodes	44 50 42 50 50 49 41 45 49 50 47 49 16 51 71
Electrodes	44 50 42 50 50 49 41 45 49 50 47 49 16 51 71 78
Electrodes	444 50 42 50 50 49 41 45 49 50 47 49 16 71 78 69 61

Needles Potocky Needle45
Obstetrical Care Fetal Dopplers
Ophthalmic Cryosurgery 54-55
Pessaries: See Urinary Incontinence
Pipelle
Pipet Curet36
Retractors 49 Stainless Steel 26
Scissors – Stainless35
Sound 39
Specula 23 Dimensions 23 LEEP 48-49 Stainless Steel 24-25
Sperm Select System67
Tenacula 50 Stainless Steel 29 Tubal Patency
ABBI
Tubing Smoke Evacuation Tubing47 Tubing Set47
Urinary Incontinence Lumax
Urodynamics Catheters
Uterine Infusion Catheters H/S Catheters
Uterine Sounds26
Uterine Structure Evaluation

ABBI72-73

Committed to Our Customers

At CooperSurgical, we strive for the highest level of expertise and knowledge of our products and to develop an ongoing relationship with each of our customers that is built on trust, quality and consistency. We look to provide service that exceeds not only our customers' expectations, but industry standards as well. CooperSurgical's Customer Service Department provides full-service order entry with extended hours to facilitate order placement. As an added convenience we also now offer online order entry. Call now to have your Customer Service Representative set up your account today. Once you are registered for online ordering, you will be able to place orders, track shipments, review prior purchase history and schedule orders for future delivery.

5 Ways to Order



Local Sales Representative

www.coopersurgical.com/pages/salesrepsearch.aspx



Email:

US: Orders@coopersurgical.com Int'l: InternationalCustomerservice@coopersurgical.com



E-Ordering (US Only):

Call 800.243.2974 to set up your account today.



Phone Orders:

800.243.2974 (8:00 AM - 8:30 PM EST) Int'l: 203.601.9818



Fax Orders:

800.262.0105 (24 hours a day) Int'l: 203.601.4747

Warranty Information

CooperSurgical, Inc. ("CooperSurgical") warrants that its products offered in this catalog, **other than** single-use products, consumables, Euro-Med handheld instruments and ClearTone Fetal Dopplers, will be free from defects in materials and workmanship for one year from the date of purchase.

CooperSurgical warrants that its single-use products and consumables offered in this catalog will, through the earlier of 90 days from the date of purchase of each such product and the expiration date of that product stated on that product's package, be free from defects in materials and workmanship that are discovered prior to or during the single use of that product.

CooperSurgical warrants that its Euro-Med handheld instruments offered in this catalog will be free from defects in materials and workmanship for two years from the date of purchase. In addition, CooperSurgical offers complementary sharpening for five years from the date of purchase only on Euro-Med Classic Series biopsy punches and only if CooperSurgical determines that enough material is available on the tip of the punch to accept sharpening.

CooperSurgical warrants that its ClearTone Fetal Dopplers offered in this catalog will be free from defects in materials and workmanship for one year from the date of purchase, and that the parts included in those products will be free from defects in materials and workmanship for five years from the date of purchase. After the expiration of the one-year product warranty applicable to a ClearTone Fetal Doppler, if CooperSurgical determines that a part included in that product fails to conform to the five-year warranty on those parts, CooperSurgical will repair or replace only that part, not the entire product, free of charge, as the sole remedy for that failure to conform. CooperSurgical may charge for time and materials required to remove that part, install a repaired or replacement part and perform related testing.

If CooperSurgical determines that a product offered in this catalog fails to conform to the CooperSurgical warranty applicable to that product, other than the five-year warranty for ClearTone Fetal Doppler parts described above, CooperSurgical will, as the sole remedy for that failure to conform, repair or replace that product, at CooperSurgical's discretion, free of charge.

To return a product to CooperSurgical, a customer must comply with CooperSurgical's Returned Goods Policy described in this catalog. A customer will not have any remedy if a CooperSurgical product does not conform to the warranty for that product unless that product is returned to CooperSurgical in accordance with that Returned Goods Policy. CooperSurgical will ship returned products that it repaired (not including sharpening) or replaced under warranty to the customer who returned those products, at CooperSurgical's expense F.O.B. the customer's facility. Under all other circumstances, including if CooperSurgical sharpens a Euro-Med Classic Series biopsy punch without otherwise repairing the punch under warranty, CooperSurgical will ship returned products to the customer who returned those products at the customer's expense F.O.B. CooperSurgical's facility.

CooperSurgical's warranties, including those in this catalog, do not cover damage caused by misuse, improper care, improper use of chemicals or cleaning methods, loss, theft, servicing or sharpening by non-authorized personnel or negligent or intentional conduct on the part of the owner or user of the product. Except for the sharpening services described above for Euro-Med Classic Series biopsy punches, CooperSurgical's warranties do not cover normal wear and tear or general maintenance, including cleaning and sharpening services required as a result of normal use. Any modifications or changes to a product will void that product's warranty.

THE ABOVE WARRANTIES ARE IN LIEU OF, AND COOPERSURGICAL HEREBY DISCLAIMS, ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, WRITTEN OR ORAL, WITH RESPECT TO COOPERSURGICAL'S PRODUCTS, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO TERMS, CONDITIONS, UNDERSTANDINGS OR AGREEMENTS THAT PURPORT TO MODIFY THE ABOVE WARRANTIES OR THAT MAKE ANY ADDITIONAL WARRANTIES FOR ANY COOPERSURGICAL PRODUCT SHALL HAVE ANY LEGAL EFFECT UNLESS MADE IN WRITING AND SIGNED BY AN AUTHORIZED COOPERSURGICAL CORPORATE OFFICER.

COOPERSURGICAL SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR LOST PROFITS, DAMAGES FROM LOSS OF USE OR LOST DATA, OR INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ITS WARRANTIES OR OTHERWISE FOR ANY CLAIM RELATED TO COOPERSURGICAL'S PRODUCTS, EVEN IF COOPERSURGICAL HAS BEEN ADVISED, KNEW OR SHOULD HAVE KNOWN OF THE POSSIBILITY OF SUCH DAMAGES. COOPERSURGICAL'S LIABILITY WITH RESPECT TO A PRODUCT COVERED BY A WARRANTY OR OTHERWISE SHALL BE LIMITED IN ALL CIRCUMSTANCES TO THE PURCHASE PRICE OF THAT PRODUCT.

Latex Information

CooperSurgical understands the concerns of the medical field for latex sensitivities. We are happy to inform you that our products are LATEX-FREE except for:

• Milex® Latex Inflatoball Pessary

Returned Goods Policy

- 1. Goods will be accepted for return for the following reasons:
 - If shipment was made without Customer's authorization or order
 - If incorrect goods were shipped
 - If defective goods are covered by the standard warranty
- 2. To return goods, you must contact a Customer Service Representative by telephone (1.800.243.2974 or 1.203.601.5200) for a Returned Merchandise Authorization (RMA) number. Items will not be accepted without an RMA number.

Please have the following information:

- Reason you wish to return the goods
- Quantity, description, part number, serial number of the goods
- Date of Receipt of Order
- Customer's Purchase Order and the CooperSurgical Invoice Number
 All used instruments must be cleaned and sterilized prior to being shipped for return.
- 3. Shipment must be sent prepaid by Customer. Freight collect shipments will not be accepted, and goods will be returned to sender.
- 4. If Customer intends to return equipment ordered in error, the following restocking charges and terms will apply:
 - There is no restocking fee if goods are returned within 30 days of shipment.
 - There is a restocking fee of 25% if goods are returned 31-60 days of shipment.
 - Goods must be returned unused, in the original carton, and in marketable condition
 - Refurbishing and replacement charges will be added to the restocking charges for damaged or missing items
 - No return after 60 days
 - No refund on sterile, single-use disposable products

Send all Returns to: CooperSurgical

95 Corporate Drive Trumbull, CT 06611

Trademark Notices

The marks bearing the symbols "TM" and "®" are trademarks and registered trademarks of CooperSurgical, Inc.

Cervex-Brush® is a registered trademark of Futura Nova.

Canon® is a registered trademark of Canon, Inc.

Velcro® is a registered trademark of Velcro Industries B.V.

Valleylab™ is a trademark of Covidien AG.

Cryomedics® is a registered trademark of Cabot Technologies Corporation

Stackhouse[™] is a trademark of Microtek Medical. Inc.

Polocaine™ is a trademark of AstraZeneca, LP

C € Available for sale where CE mark is required

CA Available for sale in Canada

CE and CA labeling is subject to change. Contact Customer Service or your local sales representative for the most-up-to-date information.

Premier Products for Women's Health Care













<u>CoperSurgical</u>

95 Corporate Drive Trumbull, CT 06611