

e-sphyg3™

Digital Blood Pressure Monitor

- **Accurate**
- **Affordable**
- **Clinical Grade**



ADC's all-new e-sphyg 3 captures blood pressure measurements with the touch of a button. Designed for busy professionals who value both accuracy and affordability, the e-sphyg 3 is the feature-packed, clinical-grade solution you've been waiting for.

Three clinically validated modes:

1x

Single measurement mode

with pulse rate and mean arterial pressure (MAP) displayed in as little as 30 seconds.

3x

Customizable, three-reading average mode

with individual measurement display.



Manual mode

for auscultation with a stethoscope.



Last-reading memory




Rechargeable high-capacity battery
STANDARD



The e-sphyg³ is perfect for:



Physician practices



Walk-in clinics



Nursing homes

Desktop design with built-in handle
Or choose one of these optional platforms
to fit your practice workflow.

Oscillometric measurements capture:

[SYS] Systolic blood pressure
(60 to 255 mmHg)

[DIA] Diastolic blood pressure
(30 to 200 mmHg)

[HR] Pulse rate
(40 to 200 BPM)

[MAP] Mean arterial pressure
(35 to 240 mmHg)

[IB] Irregular
Heartbeat

For adult and pediatric patients

INCLUDES: Three-cuff kit
(Small Adult, Adult, and
Large Adult) for accurate
readings on virtually
all patients.

Choose our traditional
Adcuff design or our new
Adcuff+ one-piece cuffs.

Wall Mount
Custom bracket for
easy mounting to
the wall or on an
ADC Adstation system.


*Includes cuff
storage basket.*



Mobile Stand

Our rugged, easy-to-assemble
mobile stand features five
rolling casters (two locking) for
room-to-room portability.

Includes cuff storage basket.



Desk Caddy

A modern, compact
design that helps keep
cuffs and cords
organized.



Dentist offices



Birth centers



Long-term care facilities

e-sphyg³TM

Digital Blood Pressure Monitor

The clinical-grade solution you've been waiting for.

Technical Comparisons

FEATURE	e-sphyg3	Vital Signs Monitors	Home BP Devices	Manual BP Devices
Automated Measurement	Yes	Yes	Yes	No
Clinical Grade NIBP	Yes	Yes	No	N/A
3 Reading Average Mode	Yes	Some	Some	N/A
Automatic and Manual Modes	Yes	Some	No	No
High-Capacity Battery - STANDARD	Yes	No	No	N/A
Includes 3 cuffs	Yes	No	No	No
Multiple Platforms	Yes	Yes	No	Yes
Irregular Heartbeat	Yes	Some	Some	No



Ordering Information

9003K-MCC e-sphyg 3 with two-piece cuffs

9003K-MCC1 e-sphyg 3 with one-piece cuffs

Both configurations include: e-sphyg 3 monitor; built-in, rechargeable, high-capacity battery; 6-foot power cord with four adapters; three Adcuff blood pressure cuffs (small adult, adult, large adult); and 4.5-foot tubing with connector.

Optional Platforms

9003D e-sphyg 3 desk caddy

9003W e-sphyg 3 wall mount bracket (compatible with ADC Adstation systems) and cuff basket

9003M e-sphyg 3 mobile stand and cuff basket

Replacement Parts

9003ADAPT AC adapter set for esphyg 3

9003BAT High-capacity battery for esphyg 3

9003CK Adcuff kit (small adult, adult, large adult) with e-sphyg 3 connectors

9003CT Extension tubing/connectors for esphyg 3



55 Commerce Drive
Hauppauge, NY 11788

Toll Free: 1-800-ADC-2670
(1-800 232-2670)

Email: tvalle@adctoday.com

Web: www.adctoday.com