



SonoTrax
Ultrasonic Pocket Doppler

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122, P.R. China
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn

© 2018 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-OB & GYN-SONOTRAX-V3.3-20180111

SonoTrax

Ultrasonic Pocket Doppler



SonoTrax

Ultrasonic Pocket Doppler



SonoTrax II Pro



Compact Battery Charger



Backlit LCD



Built-in Recorder




Interchangeable Waterproof Probe

Specifications

- Display: 45 x 25 mm LCD display
- Heart Rate Range: 50~210 bpm
- Resolution: 1 bpm
- Recording Time: 240 s
- Audio Output Power: 1W
- Auto Power-off: 1 minute with no signal or operation
- Size: 34 (D) x 89 (W) x 141 (H) mm
- Weight: about 300 g (including the battery)



Carry Case



Rolling Stand

- Easy-to-use
- Compact design
- Support 1.5 V AA battery
- Built-in audio recorder
- Large backlit LCD display
- Accurate fetal heart detection with high-fidelity sound
- Support 2 MHz, 3 MHz, 4 MHz, 5 MHz, 8 MHz waterproof probes
- Automatic and manual counting of FHR
- Interchangeable waterproof probe
- Automatic power-off



SonoTrax Vascular



SonoTrax Lite



SonoTrax Basic



SonoTrax Basic A



SonoTrax Pro



SonoTrax II



SonoTrax II Pro