


M3A

Vital Signs Monitor



Portable Design



Infrared Ear Thermometer



Covidien Quick Temp



M3A Vital Signs Monitor

EDAN M3A Vital Signs Monitor provides essential measurement of patient's SpO₂, NIBP and Temperature. It offers an effective solution and brings high performance vital signs monitoring at your finger tips.

12,000

NIBP Measurement Review

100 h

Trend Review

- 3.5" high resolution color TFT-LCD screen & LED display
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 17 hours of continuous work

- Compatible with PatientCare Viewer PC data management software
- Nurse call
- Suitable for adult, pediatric and neonatal patients
- Barcode scanner support (optional)
- Ethernet connection for bi-directional communication with MFM-CMS central monitoring system (optional)
- Wi-Fi support with optional AP



Web Viewer



MFM-CMS



Mounting Solutions



Wall Mount



Rolling Stand



Bed Rail Clamp

Specification

Configuration	SpO ₂	NIBP	Covidien Quick Temp	Infrared Ear Thermometer
M3A(SpO ₂ only)	✓	—	—	—
M3A(NIBP only)	—	✓	—	—
M3A(SpO ₂ +NIBP)	✓	✓	—	—
M3A(NIBP+Covidien Quick Temp)	—	✓	✓	—
M3A(NIBP+Infrared Ear Thermometer)	—	✓	—	✓
M3A(SpO ₂ +NIBP+ Covidien Quick Temp)	✓	✓	✓	—
M3A(SpO ₂ +NIBP+ Infrared Ear Thermometer)	✓	✓	—	✓