

APPLICATION	EMS PRO	LED+	LEDX	PEDI
Adult IV Access	✓	✓	✓	
Pediatric IV Access	✓	✓		
Emergency	✓	✓	✓	
Neonatology				✓
Sclerotherapy & Phlebectomy			✓	
Oncology & Infusion	✓	✓	✓	
Radiology & Nuclear Medicine	✓	✓	✓	
Home Health	✓	✓	✓	



FEATURES	EMS PRO	LED+	LEDX	PEDI
Number of LEDs	24	28	32	12
Color of LEDs	Orange (12) + Red (4) + White (8)	Orange (22) + Red (6)	Orange (24) + Red (8)	Orange (8) + Red (4)
Color Selection	Yes	Yes	Yes	Yes
Ring Diameter (mm)	21	27	31	14
Battery	2 x AA Alkaline	Rechargeable Lithium Ion	Rechargeable Lithium Ion	1 CR2 Lithium
Battery Run Time (mins)	550	350	130 / 160 (32 / 24)	200 (12 / 8)
Dimensions (mm)	115 x 55 x 21	99 x 57 x 21	102 x 65 x 21	67 x 42 x 21
Weight, inc. battery (g)	108	77	83	41