

Task Seating, Treatment Tables, and Privacy Screens



Task Seating, Treatment Tables & Privess<sup>™</sup> Privacy Screens from Winco provide utility, durability and comfort. Practical Tools from a trusted provider!

Healthcare facilities know Winco for their broad range of medical recliners. That same reliability is found in our task, privacy, and treatment solutions. Our heavy duty steel, lightweight aluminum and upholstered components provide stability and strength so clinicians can confidently focus on their patients. Available in an assortment of colors to coordinate with your decor.



Task Chairs (4350 / 4358 / 4380 / 4388 / 4450 / 4458



## **KEY FEATURES**

- Comfort: Adjustable, upholstered seat-back
- Durability: High-grade vinyls (25+ colors)
- Functionality: Smooth operating gas lift

| Specifications/Features | 4350 / 4358  | 4380 / 4388                            | 4450 / 4458                            |
|-------------------------|--|--|--|
| Overall Height          | 29" - 34"<br>(73 7 cm - 86.4 cm)                             | 36 5" - 46 75"<br>(92 7 cm - 118 7 cm) | 30" - 35"<br>(76 2 cm - 88 9 cm)       |
| Seat Height Range       | 17.75" - 22.75"<br>(45.1 cm - 57.8 cm)                       | 23 5" - 33.75"<br>(59.7 cm - 85.7 cm)  | 18 25" - 23 25"<br>(46.4 cm - 59 1 cm) |
| Back Height Above Seat  | 13" (33 cm)  | 13" (33 cm)                            | 13" (33 cm)                            |
| Seat Width (Diameter)   | 16 5" (41.9 cm)  | 16 5" (41.9 cm)                        | 16 5" (41.9 cm)                        |
| Adjustable Seat Depth   | 15 25" - 17.5"<br>(38 7 cm - 44 5 cm)                        | 15 25" - 17.5"<br>(38 7 cm - 44 5 cm)  | 15 25" - 17.5"<br>(38.7 cm - 44 5 cm)  |
| Overall Weight          | 20 lb (9 1 kg)   | 22 lb (10 kg)                          | 21 5 lb (9 8 kg)                       |
| Weight Capacity         | 250 (113 kg)   | 250 (113 kg)                           | 250 (113 kg)                           |
| Footprint               | 22"x 22" 22" x 22"<br>(55 9 cm x 55 9 cm) (55 9 cm x 55 9 cm |  | 22" x 22"<br>(55 9 cm x 55 9 cm)       |
| Base                    | High Density<br>Plastic                                      | High Density<br>Plastic                | Polished Aluminum                      |

### **OPTIONS**

|             | TB133       |                |                               |
|-------------|-------------|----------------|-------------------------------|
| MODEL       | TB-133 (TB) | FOOT RING (FR) | *HIGH RESISTANCE CASTERS (HR) |
| 4350 / 4358 | Optional    | Optional       | 4358                          |
| 4380 / 4388 | Optional    | Standard       | 4388                          |
| 4450 / 4458 | Optional    | Optional       | 4458                          |

\*High resistance casters allow free chair movement when weighted or occupied, but resist unnecessary chair movement when unweighted or exiting the chair. NOTE: Casters may not release and roll freely if you weigh less than 80 pounds.

WARRANTY Call, email or visit our website for complete warranty details.

**COLORS** Available in 25+ vinyl colors. Visit our website to view all colors or contact our customer care team to request swatches. **NOTE** Winco reserves the right to change specifications without notice.



Task Stools (4300 / 4308 / 4400 / 4408)



4400 / 4408

4300 / 4308

# **KEY FEATURES**

- Comfort: Padded seat
- Durability: High-grade vinyls (25+ colors)
- Functionality: "One-Touch" height adjustment

| Specifications        | 4300 / 4308                          | 4400 / 4408                         |
|-----------------------|--------------------------------------|-------------------------------------|
| Seat Height Range     | 16 5" - 21 5"<br>(41.9 cm - 54 6 cm) | 17.75" - 23"<br>(45.1 cm - 58.4 cm) |
| Seat Width (Diameter) | 16 5" (41.9 cm)                      | 16.5" (41.9 cm)                     |
| Overall Weight        | 14 lb (6.4 kg)*                      | 15 5 lb (7 kg)*                     |
| Weight Capacity       | 250 (113 kg)                         | 250 (113 kg)                        |
| Footprint             | 22"x 22"(55 9 cm x 55 9 cm)          | 22"x 22"(55 9 cm x 55 9 cm)         |
| Base                  | High Density Plastic                 | Polished Aluminum                   |

### **OPTIONS**

|             | TB133<br>Riduire MTM* |                |                               |
|-------------|-----------------------|----------------|-------------------------------|
| MODEL       | тв-133 (тв)           | FOOT RING (FR) | *HIGH RESISTANCE CASTERS (HR) |
| 4300 / 4308 | Optional              | Optional       | 4308                          |
| 4400 / 4408 | Optional              | Optional       | 4408                          |

\*High resistance casters allow free chair movement when weighted or occupied, but resist unnecessary chair movement when unweighted or exiting the chair. NOTE: Casters may not release and roll freely if you weigh less than 80 pounds.

WARRANTY Call, email or visit our website for complete warranty details.

**COLORS** Available in 25+ vinyl colors. Visit our website to view all colors or contact our customer care team to request swatches. **NOTE** Winco reserves the right to change specifications without notice.



Treatment Tables (8500 / 8510 / 8550 / 8570)



| Specifications/Features             | 8500              | 8510              | 8550              | 8570              |
|-------------------------------------|-------------------|-------------------|-------------------|-------------------|
| Overall Height (Standard)           | 30" (76 cm)       | 30" (76 cm)       | 19" (48 cm)       | 30" (76 cm)       |
| Overall Width                       | 28" (71 cm)       | 28" (71 cm)       | 28" (71 cm)       | 28" (71 cm)       |
| Overall Length                      | 72" (183 cm)      | 72" (183 cm)      | 72" (183 cm)      | 72" (183 cm)      |
| Overall Weight<br>(Standard Height) | 74 lb (33 6 kg)   | 74 lb (33 6 kg)   | 73 lb (33 1 kg)   | 91 lb (41 3 kg)   |
| Weight Capacity                     | 400 lb (181.4 kg) |
| Angle of Backrest                   | N/A               | N/A               | N/A               | 6° - 57°          |
| Face Cutout                         | N/A               | Yes               | N/A               | N/A               |

### **OPTIONS**

|       | MEEZS<br>TB133<br>Roduine ments |                            |               |   |   |
|-------|---------------------------------|----------------------------|---------------|---|---|
| MODEL | TB-133 (TB)                     | * SPECIAL HEIGHT LEGS      | ** SHELF (SH) | Additional Specifications                       |   |
| 8500  | Optional                        | Optional (36" / 25" / 19") | 8500-SH       | * Choosiel (Jaight Logo                         | 36″ (91 cm) - Add 2 lb to Overall Weight      |
| 8510  | Optional                        | Optional (36" / 25" / 19") | 8510-SH       | * Special Height Legs<br>(See OPTIONS for model | 25" (63 cm) - Deduct 2 lb from Overall Weight |
| 8550  | Optional                        | N/A                        | 8550-SH       | availability)                                   | 19" (48 cm) - Deduct 4 lb from Overall Weight |
| 8570  | Optional                        | Optional (36" / 25")       | 8570-SH       | ** Optional Shelf (SH)                          | Add 25 lb (11 3 kg) to Overall Weight         |

## **KEY FEATURES**

- Durability: Powder-coated, steel frames
- Functionality: Paper holder standard on all models
- Stability: Adjustable leveling feet



Privess™ Privacy Screens

Winco's line of Privess<sup>™</sup> Privacy Screens provides a broad range of solutions that meet a wide array of requirements.



# PRIVESS<sup>™</sup> UNIPANEL Model 3420 Shown

Perfect for confined spaces!

Available in standard White vinyl (model 3420) or your choice of 5 optional Sure-Chek vinyl colors (model 3429)



### PRIVESS<sup>™</sup> BASIC Model 3130 Shown

Portable privacy at a cost-effective price-point!

Available in standard White vinyl (model 3130) or your choice of 5 optional Sure-Chek vinyl colors (model 3139)



### PRIVESS™ ELITE Model 3170 Shown

Style, size, and stability suitable for multiple settings!

Available in 40+ colors. \*Grade B/C vinyl & Silvertex optional, additional fees apply



#### PRIVESS<sup>™</sup> SWING-AWAY Model 3409 Shown

A fixed placement solution that's convenient to use and store!

Available in standard White vinyl (model 3400) or your choice of 5 optional Sure-Chek vinyl colors (model 3409)



### PRIVESS™ MODULAR Model 3730 Shown

Features a unique snaptogether framework that provides 360° swivel & infinite connectivity.

Two Widths: 3 Panel (3730/3739) or 4 Panel (3740 / 3749),

Available in standard White vinyl (models 3730 & 3740) or your choice of 5 optional Sure-Chek vinyl colors (models 3739 & 3749)



Shown: 3739-91+3730+3739-94

# INFINITE CONNECTIVITY

Multiple modular units are easily connected & can be maneuvered into different shapes to meet your needs

### **MODELS & SPECIFICATIONS**

| SPECIFICATIONS      | BASIC                     | ELITE                     | SWING-AWAY   | UNI-PANEL                  | MOD                        | ULAR                       |
|---------------------|---------------------------|---------------------------|--|----------------------------|----------------------------|----------------------------|
|                     | 3130 / 3139               | 3170                      | 3400 / 3409  | 3420 / 3429                | 3730 / 3739                | 3740 / 3749                |
| Base Size           | 20"x27" (50_8cm x 68_9cm) | 23"x31" (58.4cm x 78.7cm) | N/A  | 20"x42" (50.8cm x 106.7cm) | 15"x40" (38_1cm x 101.6cm) | 15"x70" (38.1cm x 177.8cm) |
| Overall Height      | 70" (177.8cm)             | 70" (177.8cm)             | 72" (182.9cm) - recommended                                  | 69" (175.3cm)              | 69" (175.3cm)              | 69" (175.3cm)              |
| Overall Width       | 79" (200.7cm)             | 84" (213.4cm)             | 40"-92.5" (101.6-235cm)                                      | 42" (106.7cm)              | 60" (152.4cm)              | 80" (203_2cm)              |
| Single Panel Length | 49" (124_5cm)             | 55" (139.7cm)             | 69" (175.3cm)  | 49.5" (125.7cm)            | 46" (116.8cm)              | 46" (116.8cm)              |
| Single Panel Width  | 27" (68.6cm)              | 28" (71.1cm)              | 92.5" (235cm)  | 42" (106.7cm)              | 20" (50.8cm)               | 20" (50.8cm)               |
| Weight              | 25 lb (11.3 kg)           | 45.5 lb (20.6 kg)         | 8 lb (3 6 kg)  | 15 lb (6.8 kg)             | 9 lb (4 <sub>e</sub> 1 kg) | 12 lb (5.4 kg)             |
| FrameType           | Powder-Coated Steel       | Powder-Coated Steel       | Powder-Coated Steel (wall mount)<br>Aluminum (extension rod) | Powder-Coated Steel        | Aluminum                   | Aluminum                   |
| # of Panels         | 3                         | 3                         | 1  | 1                          | 3                          | 4                          |

### **OPTIONS**

| MODELS     | FOOT OPTIONS |          | VINYL OPTIONS         |                       |                 |                 |                 |
|------------|--------------|----------|-----------------------|-----------------------|-----------------|-----------------|-----------------|
|            | CRUTCH TIP   | CASTERS  | WHITE                 | SURE-CHEK             | GRADE A VINYL   | GRADE B/C VINYL | SILVERTEX       |
| BASIC      | Standard     | Standard | Standard (3130)       | Optional (3139)       | N/A             | N/A             | N/A             |
| ELITE      | Standard     | Standard | Standard (3170)       | Standard (3170)       | Standard (3170) | Optional (3170) | Optional (3170) |
| SWING-AWAY | N/A          | N/A      | Standard (3400)       | Optional (3409)       | N/A             | N/A             | N/A             |
| UNI-PANEL  | Standard     | Standard | Standard (3420)       | Optional (3429)       | N/A             | N/A             | N/A             |
| MODULAR    | Standard     | N/A      | Standard (3730, 3740) | Optional (3739, 3749) | N/A             | N/A             | N/A             |

### SURE-CHEK® COLORS

| FAWN (91) | BLUSH (92) | MINT (93) | LIGHT BLUE (94) | DUNE (95) |
|-----------|------------|-----------|-----------------|-----------|
|           |            |           |                 |           |

Winco's Privess<sup>TM</sup> Privacy Screens are not intended for wet environments and should not be used as shower curtains or in outdoor settings.

**WARRANTY** Call, email or visit our website for complete warranty details. **COLORS** Visit our website to view colors or contact customer care team to request swatches **NOTE** Winco reserves the right to change specifications without notice.

004517 Rev A

Additional fees apply to all Optional colors



To arrange a demonstration or to speak with one of our associates, please feel free to contact us at: 800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com